iPM 12/iPM 10/iPM 8

Patient Monitor

Operator's Manual



© Copyright 2011 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. For this Operator's Manual, the issue date is December 2011.



WARNING

Federal Law (USA) restricts this device to sale by or on the order of a physician.

Intellectual Property Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this Mindray product and this manual. This manual may refer to information protected by copyrights or patents and does not convey any license under the patent rights of Mindray, nor the rights of others.

Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

Release, amendment, reproduction, distribution, rental, adaption and translation of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

mindray MINDRAY are the registered trademarks or trademarks owned by Mindray in China and other countries. All other trademarks that appear in this manual are used only for editorial purposes without the intention of improperly using them. They are the property of their respective owners.

Responsibility on the Manufacturer Party

Contents of this manual are subject to changes without prior notice.

All information contained in this manual is believed to be correct. Mindray shall not be liable for errors contained herein nor for incidental or consequential damages in connection with the furnishing, performance, or use of this manual.

Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- all installation operations, expansions, changes, modifications and repairs of this product are conducted by Mindray authorized personnel;
- the electrical installation of the relevant room complies with the applicable national and local requirements;
- the product is used in accordance with the instructions for use.

\bigwedge

WARNING

- This equipment must be operated by skilled/trained clinical professionals.
- It is important for the hospital or organization that employs this equipment to carry out a reasonable service/maintenance plan. Neglect of this may result in machine breakdown or personal injury.

Warranty

THIS WARRANTY IS EXCLUSIVE AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.

Exemptions

Mindray's obligation or liability under this warranty does not include any transportation or other charges or liability for direct, indirect or consequential damages or delay resulting from the improper use or application of the product or the use of parts or accessories not approved by Mindray or repairs by people other than Mindray authorized personnel.

This warranty shall not extend to

- Malfunction or damage caused by improper use or man-made failure.
- Malfunction or damage caused by unstable or out-of-range power input.
- Malfunction or damage caused by force majeure such as fire and earthquake.
- Malfunction or damage caused by improper operation or repair by unqualified or unauthorized service people.
- Malfunction of the instrument or part whose serial number is not legible enough.
- Others not caused by instrument or part itself.

Company Contact

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Mindray Building, Keji 12th Road South, High-tech industrial park,

Address: Nanshan,Shenzhen 518057,P.R.China

 Website:
 www.mindray.com

 E-mail Address:
 service@mindray.com

 Tel:
 +86 755 81888998

 Fax:
 +86 755 26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraβe 80, 20537 Hamburg, Germany

Tel: 0049-40-2513175 Fax: 0049-40-255726

Preface

Manual Purpose

This manual contains the instructions necessary to operate the product safely and in accordance with its function and intended use. Observance of this manual is a prerequisite for proper product performance and correct operation and ensures patient and operator safety.

This manual is based on the maximum configuration and therefore some contents may not apply to your product. If you have any question, please contact us.

This manual is an integral part of the product. It should always be kept close to the equipment so that it can be obtained conveniently when needed.

Intended Audience

This manual is geared for clinical professionals who are expected to have a working knowledge of medical procedures, practices and terminology as required for monitoring of critically ill patients.

Illustrations

All illustrations in this manual serve as examples only. They may not necessarily reflect the setup or data displayed on your patient monitor.

Conventions

- *Italic text* is used in this manual to quote the referenced chapters or sections.
- [] is used to enclose screen texts.
- → is used to indicate operational procedures.

Contents

1 Safety	1-1
1.1 Safety Information	1-1
1.1.1 Dangers	
1.1.2 Warnings	1-2
1.1.3 Cautions	
1.1.4 Notes	
1.2 Equipment Symbols	1-4
2 The Basics ·····	2-1
2.1 Monitor Description	2-1
2.1.1 Intended Use	2-1
2.2 Main Unit	2-2
2.2.1 Front View	2-2
2.2.2 Side View	
2.2.3 Rear View	2-6
2.3 Modules	2-8
2.4 Display Screen	
2.5 QuickKeys	2-11
3 Basic Operations	3-1
3.1 Installation	3-1
3.1.1 Unpacking and Checking	3-1
3.1.2 Environmental Requirements	3-2
3.2 Getting Started	3-2
3.2.1 Turning Power On	3-2
3.2.2 Starting Monitoring	
3.3 Disconnecting from Power	3-3
3.4 Using Keys	3-3
3.5 Using the Touchscreen	3-3
3.6 Setting the Screen	3-4
3.7 Using the Main Menu	3-5
3.8 Removing a Label Conflict	3-5
3.9 Changing General Settings	3-5
3.9.1 Setting up a Monitor	
3.9.2 Changing Language	3-6
3.9.3 Adjusting the Screen Brightness	3-6
3.9.4 Showing/Hiding the Help	
3.9.5 Setting the Date and Time	
3.9.6 Adjusting Volume	
3.10 Operating Modes	
3.10.1 Monitoring Mode	
3.10.2 Night Mode	
3.10.4 Para Made	
3.10.4 Demo Mode	
3.10.5 Standby Mode	
4 Managing Patients	
4.1 Admitting a Patient	
4.2 Quick Admitting a Patient	
4.3 Editing Patient Information	
4.4 Discharging a Patient	
4.5 Transferring a Patient	
4.6 Connecting to a Central Monitoring System	4-4

5 Managing Configurations	5-1
5.1 Introduction	5-1
5.2 Entering the [Manage Configuration] Menu	5-2
5.3 Changing Department	5-2
5.4 Setting Default Configuration	5-3
5.5 Saving Current Settings	5-3
5.6 Editing Configuration	5-4
5.7 Deleting a Configuration	5-5
5.8 Transferring a Configuration	
5.9 Loading a Configuration	5-5
5.10 Restoring the Latest Configuration Automatically	5-6
5.11 Modifying Password	5-6
6 User Screens ·····	
6.1 Tailoring Your Screens	
6.1.2 Changing Measurement Colors	
6.1.3 Changing Screen Layout	
6.2 Viewing Minitrends	
6.2.1 Having a Split-Screen View of Minitrends	
6.2.2 Setting Minitrends	
6.3 Viewing oxyCRG	6-3
6.4 Viewing Other Patients	6-4
6.4.1 Care Group	6-4
6.4.2 Viewing the Care Group Overview Bar	6-4
6.4.3 Understanding the View Other Patient Window	
6.5 Understanding the Big Numerics Screen	6-6
7 Alarms	7.1
7.1 Alarm Categories	
7.1 Alarm Categories	
7.3 Alarm Indicators	
7.3.1 Alarm Message	
7.3.3 Flashing Numeric	
7.3.4 Audible Alarm Tones	
7.3.5 Alarm Status Symbols	
7.4 Alarm Tone Configuration	7-4
7.4.1 Setting the Minimum Alarm Volume	7-4
7.4.2 Changing the Alarm Volume	
7.4.3 Setting the Interval between Alarm Sounds	
7.4.4 Changing the Alarm Tone Pattern	
7.4.5 Setting the Reminder Tones	
7.5 Understanding the Alarm Setup Menu	
7.5.1 Setting Alarm Properties for All Parameters	
7.5.3 Setting Alarm Delay Time	
7.5.4 Setting Recording Length	
7.5.5 Entering CPB Mode	
7.6 Pausing Alarms	
7.7 Switching Off All Alarms	
7.8 Silencing the Alarm Sound	
7.9 Latching Alarms	
7.10 Silencing Technical Alarms	
7.11 Testing Alarms	
7.12 When an Alarm Occurs	
7.13 Using Care Group Alarms	
7.13.1 Care Group Auto Alarms	
7.13.2 Silencing Care Group Alarms	
8 Monitoring ECG	8-1

8.1 Introduction	
8.2 Safety	8-1
8.3 Preparing to Monitor ECG	8-2
8.3.1 Preparing the Patient and Placing the Electrodes	
8.3.2 Choosing AHA or IEC Lead Placement	
8.3.3 ECG Lead Placements	
8.3.4 Checking Paced Status	
8.4 Understanding the ECG Display	
8.5 Changing ECG Settings	
8.5.1 Accessing ECG Menus	
8.5.2 Setting Pacemaker Rate (For Mortara only)	
8.5.4 Setting the ECG Lead Set	
8.5.5 Choosing an ECG Display Screen	
8.5.6 Changing the ECG Filter Settings	
8.5.7 Switching the Notch Filter On or Off	
8.5.8 Changing the Pacer Reject Settings	
8.5.9 About the Defibrillator Synchronization	
8.5.10 Changing ECG Wave Settings	
8.5.11 Enabling Smart Lead Off	8-7
8.5.12 Setting the Alarm Level for ECG Lead Off Alarms	
8.5.13 Adjusting QRS Volume	8-8
8.6 About ST Monitoring	
8.6.1 Switching ST On and Off	
8.6.2 Changing ST Filter Settings	
8.6.3 Understanding the ST Display	
8.6.4 Saving the Current ST Segment as Reference	
8.6.5 Changing the Reference Segment	
8.6.7 Recording the ST Segment	
8.6.8 Changing the ST Alarm Limits	
8.6.9 Setting the ST Alarm Delay Time	
8.6.10 Adjusting ST Measurement Points	
8.7 About Arrhythmia Monitoring	
8.7.1 Understanding the Arrhythmia Events	
8.7.2 Changing Arrhythmia Alarm Settings	
8.7.3 Changing Arrhythmia Threshold Settings	8-14
8.7.4 Setting the Extended Arrhythmia(For Mindray Algorithm only)	8-15
8.7.5 Reviewing Arrhythmia Events	8-15
8.8 ECG Relearning	8-16
8.8.1 Initiating an ECG Relearning Manually	8-16
8.8.2 Automatic ECG Relearning	8-16
8.9 12-Lead ECG Monitoring (for iPM 12 and iPM 10 patient monitors only)	8-17
9 Monitoring Respiration (Resp)	9-1
9.1 Introduction	
9.2 Safety Information	
9.3 Understanding the Resp Display	
9.4 Placing Resp Electrodes	
9.4.1 Optimizing Lead Placement for Resp	
9.4.2 Cardiac Overlay	
9.4.3 Abdominal Breathing	
9.4.4 Lateral Chest Expansion	9-3
9.5 Choosing the Respiration Lead	9-3
9.6 Changing the Apnea Alarm Delay	
9.7 Changing Resp Detection Mode	
9.8 Changing Resp Wave Settings	
9.9 Setting RR Source	
9.10 Setting alarm properties	
9.11 Switching Resp Measurement On/Off	
10 Monitoring PR	10-1

10.1 Introduction	10-1
10.2 Setting the PR Source	10-1
10.3 Selecting the Active Alarm Source	10-2
10.4 QRS Tone	10-2
11 Monitoring SpO ₂	
11.1 Introduction	
11.2 Safety	
11.3 Identifying SpO ₂ Modules	
11.4 Applying the Sensor	
11.5 Changing SpO ₂ Settings	
11.5.1 Accessing SpO ₂ Menus	
11.5.2 Adjusting the Desat Alarm	
11.5.3 Setting SpO ₂ Sensitivity	
11.5.4 Changing Averaging Time	
11.5.5 Monitoring SpO ₂ and NIBP Simultaneously	
11.5.7 Changing the Speed of the Pleth Wave	
11.5.8 Setting the Alarm Level for SpO ₂ Sensor Off Alarm	
11.5.9 Setting the SpO ₂ Tone Mode	
11.5.9 Setting the 5pO ₂ Tone Mode	
11.7 Masimo Information	
11.8 Nellcor Information	
11.8 Neilcor information	11-0
12 Monitoring NIBP ·····	12-1
12.1 Introduction	12-1
12.2 Safety	12-2
12.3 Measurement Limitations	
12.4 Measurement Methods	
12.5 Setting Up the NIBP Measurement	
12.5.1 Preparing to Measure NIBP	
12.5.2 Starting and Stopping Measurements	
12.5.3 Correcting the Measurement if Limb is not at Heart Level	
12.5.4 Enabling NIBP Auto Cycling and Setting the Interval	
12.5.5 Starting a STAT Measurement	
12.6 Understanding the NIBP Numerics	12-4
12.7 Changing NIBP Settings	
12.7.1 Setting the Initial Cuff Inflation Pressure	
12.7.2 Setting NIBP Alarm Properties	12-5
12.7.3 Displaying NIBP List	
12.7.4 Setting the Pressure Unit	12-5
12.8 Assisting Venous Puncture	
12.9 NIBP Leakage Test	12-6
12.10 NIBP Accuracy Test	
12.11 Calibrating NIBP	
13 Monitoring Temp	
13.1 Introduction	13-1
13.2 Safety	13-1
13.3 Making a Temp Measurement	13-1
13.4 Understanding the Temp Display	13-1
13.5 Setting the Temperature Unit	
14 Monitoring IBP	
14.1 Introduction	
14.2 Safety	
14.3 Zeroing the Transducer	
14.4 Setting Up the Pressure Measurement	
14.5 Understanding the IBP Display	
14.6 Changing IBP Settings	14-4

14.6.1 Changing a Pressure for Monitoring	14-4
14.6.2 Setting Alarm Properties	
14.6.3 Changing Averaging Time	
14.6.4 Setting the Pressure Unit	
14.6.5 Setting Up the IBP Wave	
14.7 Measuring PAWP	
14.7.1 Preparing to Measure PAWP	
14.7.2 Setting Up the PAWP Measurement	
14.7.3 Understanding the PAWP Setup Menu	14-7
15 Monitoring Cardiac Output ·····	15-1
15.1 Introduction	
15.2 Understanding the C.O. Display	
15.3 Influencing Factors	
15.4 Setting Up the C.O. Measurement	
15.5 Measuring the Blood Temperature	
15.6 Changing C.O. Settings	
15.6.1 Setting the Temperature Unit	
15.6.2 Setting Alarm Properties	
16 Monitoring Carbon Dioxide	16-1
16.1 Introduction	16-1
16.2 Measuring CO ₂	16-2
16.2.1 Making a Sidestream CO ₂ Measurement	
16.2.2 Making a Microstream CO ₂ Measurement	
16.2.3 Making a Mainstream CO₂ Measurement	
16.3 Changing CO ₂ Settings	
16.3.1 Accessing CO ₂ Menus	
16.3.2 Entering the Standby Mode	
16.3.3 Setting the CO_2 Unit	
16.3.5 Setting up Humidity Compensation	
16.3.6 Setting the Apnea Alarm Delay	
16.3.7 Choosing a Time Interval for Peak-Picking	
16.3.8 Setting the Flow Rate	
16.3.9 Setting up the CO ₂ Wave	16-6
16.4 Setting RR Source	16-7
16.5 Setting Barometric Pressure Compensation	16-7
16.6 Measurement Limitations	16-7
16.7 Leakage test	16-7
16.8 Troubleshooting the Sidestream CO ₂ Sampling System	16-8
16.9 Removing Exhaust Gases from the System	
16.10 Zeroing the Sensor	
16.10.1 For Sidestream and Microstream CO ₂ Modules	
16.10.2 For Mainstream CO ₂ Modules	16-9
16.11 Calibrating the Sensor	16-9
16.12 Oridion Information	16-10
17 Monitoring AG	
17.1 Introduction	
17.2 Understanding the AG Display	
17.3 MAC Values	
17.4 Preparing to Measure AG	
17.5 Changing AG Settings	
17.5.1 Setting Gas Unit	
17.5.2 Setting the Apnea Alarm Delay17.5.3 Changing the Sample Flow Rate	
17.5.4 Setting up the O ₂ Compensation	
17.5.5 Entering the Standby Mode	
17.5.6 Setting up the AG Wave	
17.5.7 Setting RR Source	

17.6 Changing the Anesthetic Agent	17-6
17.7 Measurement Limitations	17-6
17.8 Troubleshooting	17-6
17.8.1 When the Gas Inlet is Blocked	17-6
17.8.2 When an Internal Occlusion Occurs	17-7
17.9 Removing Exhaust Gases from the System	17-7
18 Freezing Waveforms ·····	18-1
18.1 Freezing Waveforms	18-1
18.2 Viewing Frozen Waveforms	18-1
18.3 Unfreezing Waveforms	18-1
18.4 Recording Frozen Waveforms	
19 Review	19-1
19.1 Accessing Respective Review Windows	
19.2 Reviewing Graphic Trends	
19.3 Reviewing Tabular Trends	
19.4 Events	
19.4.1 Marking Events	
19.4.2 Reviewing Events	
19.5 Reviewing Waveforms	
•	
20 Calculations 20.1 Introduction.	
20.2 Dose Calculations	
20.2.1 Performing Calculations	
20.2.3 Titration Table	
20.3 Oxygenation Calculations	
20.3.1 Performing Calculations	
20.3.2 Entered Parameters	
20.3.3 Calculated Parameters	
20.4 Ventilation Calculations	
20.4.1 Performing Calculations	
20.4.2 Entered Parameters	
20.4.3 Calculated Parameters	
20.5 Hemodynamic Calculations	
20.5.1 Performing Calculations	
20.5.2 Entered Parameters	
20.5.3 Calculated Parameters	
20.6 Renal Calculations	
20.6.1 Performing Calculations	
20.6.2 Entered Parameters	
20.6.3 Calculated Parameters	
20.7 Understanding the Review Window	20-8
21 Recording	21-1
21.1 Using a Recorder	
21.2 Overview of Recording Types	
21.3 Starting and Stopping Recordings	
21.4 Setting up the Recorder	
21.4.1 Accessing the Record Setup Menu	
21.4.2 Selecting Waveforms for Recording	
21.4.3 Setting the Realtime Recording Length	
21.4.4 Setting the Interval between Timed Recordings	
21.4.5 Changing the Recording Speed	
21.4.6 Clearing Recording Tasks	
21.5 Loading Paper	
21.6 Removing Paper Jam	
21.7 Cleaning the Recorder Printhead	
21.7 Clearning the necotaer i fillithead	Z1 ⁻⁴

22 Printing	22-1
22.1 Printer	22-1
22.2 Connecting a printer	22-1
22.3 Setting Up the Printer	22-2
22.4 Starting Report Printouts	
22.5 Stopping Reports Printouts	
22.6 Setting Up Reports	
22.6.1 Setting Up ECG Reports	
22.6.2 Setting Up Tabular Trends Reports	
22.6.3 Setting Up Graphic Trends Reports	
22.6.4 Setting Up Realtime Reports	22-4
22.7 End Case Reports	22-4
22.8 Printer Statuses	22-4
22.8.1 Printer Out of Paper	22-4
22.8.2 Printer Status Messages	22-4
23 Other Functions	23-1
23.1 Analog Output	
23.2 Transferring Data	
23.2.1 Data Export System	
23.2.2 Transferring Data by Different Means	
23.3 Nurse Call	
23.4 Network Connection	
23.4.1 Setting the Network Type	
23.4.2 Wireless Network	
23.4.3 Setting the IP Address, Subnet Mask and Gateway	
24 Batteries ·····	
24.1 Overview	
Replacing a Battery	
24.2	24-2
24.3 Battery Guidelines	24-3
24.4 Battery Maintenance	24-3
24.5 Battery Recycling	
25 Care and Cleaning	25-1
25.1 General Points	
25.2 Cleaning	
25.3 Disinfecting	
26 Maintenance	
26.1 Safety Checks	
26.2 Checking Monitor and Module Information	
26.3 Calibrating ECG	
26.4 Calibrating the Touchscreen	26-2
26.5 Calibrating CO ₂	26-2
26.6 Calibrating AG	26-3
27 Accessories ·····	27-1
27.1 ECG Accessories	
27.2 SpO ₂ Accessories	
27.3 NIBP Accessories	
27.4 Temp Accessories	
27.5 IBP/ICP Accessories	
27.6 C.O. Accessories	
27.7 CO ₂ Accessories	
27.8 AG Accessories	
27.9 Others	27-9
A Product Specifications	Δ_1
· · · · · · · · · · · · · · · · · · ·	·········· \\\\\

A.1 Monitor Safety Specifications	A-1
A.2 Power Supply Specifications	A-2
A.3 Physical Specifications	A-3
A.4 Hardware Specifications	
A.5 Data Storage	A-5
A.6 Wireless Network	
A.7 Measurement Specifications	A-5
B EMC And Radio Regulatory Compliance	В-1
B.1 EMC	B-1
B.2 Radio Regulatory Compliance	
C Default Configurations	C-1
C.1 Parameters Configuration	C-1
C.2 Routine Configuration	
C.3 User Maintenance Items	
D Alarm Messages ·····	D-1
D.1 Physiological Alarm Messages	D-1
D.2 Technical Alarm Messages	
E Symbols and Abbreviations	E-1
E.1 Symbols	E-1
E.2 Abbreviations	E-2

1 Safety

1.1 Safety Information



DANGER

• Indicates an imminent hazard that, if not avoided, will result in death or serious injury.



WARNING

Indicates a potential hazard or unsafe practice that, if not avoided, could result in death or serious injury.



CAUTION

 Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor personal injury or product/property damage.

NOTE

• Provides application tips or other useful information to ensure that you get the most from your product.

1.1.1 Dangers

There are no dangers that refer to the product in general. Specific "Danger" statements may be given in the respective sections of this manual.

1.1.2 Warnings



WARNING

- This equipment is used for single patient at a time.
- Before putting the system into operation, the operator must verify that the equipment, connecting cables and accessories are in correct working order and operating condition.
- The equipment must be connected to a properly installed power outlet with protective earth contacts
 only. If the installation does not provide for a protective earth conductor, disconnect it from the power line
 and operate it on battery power, if possible.
- To avoid explosion hazard, do not use the equipment in the presence of flammable anesthetics, vapors or liquids.
- Do not open the equipment housings. All servicing and future upgrades must be carried out by the personnel trained and authorized by our company only.
- Do not come into contact with patients during defibrillation. Otherwise serious injury or death could result.
- Do not touch the equipment's metal parts or connectors when in contact with the patient; otherwise patient injury may result.
- Verify that the DC power supply meets the requirements specified in A.2 Power Source Specifications.
- Do not rely exclusively on the audible alarm system for patient monitoring. Adjustment of alarm volume to a low level or off may result in a hazard to the patient. Remember that alarm settings should be customized according to different patient situations and always keeping the patient under close surveillance is the most reliable way for safe patient monitoring.
- The physiological data and alarm messages displayed on the equipment are for reference only and cannot be directly used for diagnostic interpretation.
- To avoid inadvertent disconnection, route all cables in a way to prevent a stumbling hazard. Wrap and secure excess cabling to reduce risk of entanglement or strangulation by patients or personnel.
- Dispose of the package material, observing the applicable waste control regulations and keeping it out of children's reach.
- Ensure that the patient monitor is supplied with continuous electric power during work. Sudden power failure leads to the loss of patient data.

1.1.3 Cautions



CAUTION

- To ensure patient safety, use only parts and accessories specified in this manual.
- At the end of its service life, the equipment, as well as its accessories, must be disposed of in compliance with the guidelines regulating the disposal of such products. If you have any questions concerning disposal of the equipment, please contact us.
- Magnetic and electrical fields are capable of interfering with the proper performance of the equipment.
 For this reason make sure that all external devices operated in the vicinity of the equipment comply with the relevant EMC requirements. Mobile phone, X-ray equipment or MRI devices are a possible source of interference as they may emit higher levels of electromagnetic radiation.
- Before connecting the equipment to the power line, check that the voltage and frequency ratings of the power line are the same as those indicated on the equipment's label or in this manual.
- Always install or carry the equipment properly to avoid damage caused by drop, impact, strong vibration
 or other mechanical force.
- Dry the equipment immediately in case of rain or water spray.

1.1.4 Notes

NOTE

- Put the equipment in a location where you can easily see the screen and access the operating controls.
- Keep this manual in the vicinity of the equipment so that it can be obtained conveniently when needed.
- The software was developed in compliance with IEC60601-1-4. The possibility of hazards arising from software errors is minimized.
- This manual describes all features and options. Your equipment may not have all of them.

1.2 Equipment Symbols

\wedge	Attention: Consult accompanying documents (this manual).		
0/0	Power ON/OFF (for a part of the equipment)	-+	Battery indicator
~	Alternating current (AC)	X	Alarms paused
*	Alarm silenced	8	Record
M	Freeze/unfreeze waveforms		Main menu
€	NIBP start/stop key	•	Inserted direction
~/==	Alternating/Direct current	===	Direct current (DC)
4	Equipotential grounding	\rightarrow	VGA output
•	USB connector	몶	Network connector
\longrightarrow	Gas outlet	→	Multifunctional output
M	Manufacture date	SN	Serial number
(€ ₀₁₂₃	The product bears CE mark indicating its conformity with the provisions of the Council Directive 93/42/EEC concerning medical devices and fulfils the essential requirements of Annex I of this directive.		
EC REP	European community representative		
	ESD warning symbol for electrostatic sensitive devices.		
4 9 F	Type CF applied part. Defibrillator-proof protection against electric shock.		
4 <u>*</u>	Type BF applied part. Defibrillator-proof protection against electric shock.		
	Dispose of in accordance to your country's requirements.		

NOTE

• Some symbols may not appear on your equipment.

2 The Basics

2.1 Monitor Description

2.1.1 Intended Use


This patient monitor is intended to be used for monitoring, displaying, reviewing, storing and transferring of multiple physiological parameters including ECG, heart rate (HR), respiration (Resp), temperature (Temp), SpO₂, pulse rate (PR), non-invasive blood pressure (NIBP), invasive blood pressure (IBP), cardiac output (C.O.), carbon dioxide (CO₂), oxygen (O₂), and anesthetic gas (AG).

This patient monitor is intended for use only by clinical professionals or under their guidance. It must only be used by persons who have received adequate training in its use. Anyone unauthorized or untrained must not perform any operation on it.

The monitor can be used in pre-hospital and hospital environments, including but not limited to, ICU, general ward, outpatient department, emergency room, operating room, recovery room and preoperative observation ward, and etc.

2.2 Main Unit

2.2.1 Front View



1. Alarm lamp

When a physiological alarm or technical alarm occurs, this lamp will flash as defined below.

♦ High level alarms: the lamp quickly flashes red.

◆ Medium level alarms: the lamp slowly flashes yellow.

◆ Low level alarms: the lamp lights yellow without flashing.

2. Display Screen

3. AC power LED (iPM 12/iPM 10)

AC/DC power LED (iPM 8)

- 4. Power On/Off Switch
 - Pressing this switch turns the patient monitor on.
 - ♦ When the monitor is on, pressing and holding this switch turns the monitor off.

An indicator is built in this switch. It turns on when the patient monitor is on and turns off when the patient monitor is off.

5. Battery LED

- ♦ On: when the battery is installed and the AC source is connected.
- ◆ Off: when no battery is installed or the installed battery is malfunction, or no AC source is connected when the patient monitor is power off.
- ♦ Flash: when the patient monitor operates on battery power.
- 6. A Press to silence all alarm sounds.
- 7. A Press to pause or restore alarms.
- 8. Press to freeze or unfreeze waveforms.
- 9. Press to start or stop recordings.
- 10. Press to start or stop NIBP measurements.
- 11.


If no menu is displayed on the screen, pressing it will enter the main menu. If there is a menu displayed on the screen, pressing it will close that menu.

12. Knob

Rotate the Knob clockwise or anti-clockwise. With each click, the highlight jumps to the neighboring item. When you reach your desired item, press the Knob to select it.


2.2.2 Side View

2.2.2.1 iPM 12



- 1. Handle
- 3. Recorder
- 5. Connector for Temp probe 1
- 7 Connector for IBP cable
- 9. Connector for ECG cable
- 2. Battery compartment
- 4. Parameter module slot
- 6. Connector for Temp probe 2
- 8. Connector for SpO₂ cable
- 10. Connector for NIBP cuff

2.2.2.2 iPM 10/iPM 8




- 1. Handle
- 3. Recorder
- 5. Connector for Temp probe 1
- 7. Connector for SpO₂ cable
- 9. Connector for NIBP cuff

- 2. Battery compartment
- 4. Parameter module slot
- 6. Connector for Temp probe 2
- 8 Connector for ECG cable

2.2.3 Rear View

2.2.3.1 iPM 12/iPM 10



1. AC Power Input

2. Equipotential Grounding Terminal

When the patient monitor and other devices are to be used together, their equipotential grounding terminals should be connected together, eliminating the potential difference between them.

3. Parameter Module Slot

Used for connecting the parameter modules.

4. USB Connectors

It connects a USB disk for data or configuration transfer.

5. Multifunctional Connector

It outputs defibrillator synchronization signals, nurse call signals and analogy output signals.


6. Network Connector

It is a standard RJ45 connector which connects the patient monitor to CMS or other patient monitor for remote view. It also connects the patient monitor to PC for system upgrade.

7. VGA Connector

It connects a secondary display, which extends the display capability of your monitor. The contents displayed on the secondary display screen accords with those displayed on the monitor screen.

2.2.3.2 iPM 8



1. Equipotential Grounding Terminal

When the patient monitor and other devices are to be used together, their equipotential grounding terminals should be connected together, eliminating the potential difference between them.

2.AC Power Input

3. Parameter Module Slot

Used for connecting the parameter modules.

4. Multifunctional Connector

It outputs defibrillator synchronization signals, nurse call signals and analogy output signals.

5. DC Power Input

It is a standard DC power input connector with positive inside and negative outside. It can be connected to the external DC power supply.

6. USB Connectors

It connects a USB disk for data or configuration transfer.

7. Network Connector

It is a standard RJ45 connector which connects the patient monitor to CMS or other patient monitor for remote view. It also connects the patient monitor to PC for system upgrade.

8. VGA Connector

It connects a secondary display, which extends the display capability of your monitor. The contents displayed on the secondary display screen accords with those displayed on the monitor screen.

2.3 Modules




1.	IBP+C.O. Module:	Contains IBP cable connector and C.O. cable connector.
2.	IBP+C.O. + Sidestream CO ₂ Module:	Contains IBP cable connector, C.O. cable connector, CO ₂
	IDP+C.O. + Sidestream CO ₂ Module:	watertrap slot, and CO ₂ gas outlet.
3. IBP+C.	IDD CO Missostra and CO Madula	Contains IBP cable connector, C.O. cable connector, CO ₂
	IBP+C.O. + Microstream CO ₂ Module:	sampling line connector, and and CO ₂ gas outlet.
4	IBP+C.O.+ Mainstream CO ₂ Module:	Contains IBP cable connector, C.O. cable connector, and CO ₂
4.		transducer connector.
5.	IBP+C.O.+AG Module:	Contains IBP cable connector, C.O. cable connector, AG
		watertrap slot, and AG gas outlet.

NOTE

• The above modules support 2 invasive blood pressures through a dual-receptacle extended cable.

2.4 Display Screen

This patient monitor adopts a high-resolution LED display to display patient parameters and waveforms. A typical display screen is shown below.



1. Patient Information Area

This area shows the patient information such as department, bed number, patient name, patient category and paced status.

- indicates that no patient is admitted or the patient information is incomplete.
- indicates that the patient has a pacer.

If no patient is admitted, selecting this area will enter the [**Patient Setup**] menu. If a patient has been admitted, selecting this area will enter the [**Patient Demographics**] menu.

2. Alarm Symbols

- indicates alarms are paused.
- indicates alarm sounds are paused.
- indicates alarm sounds are turned off.
- indicates the system is in alarm off status.

3. Technical Alarm Area

This area shows technical alarm messages and prompt messages. When multiple messages come, they will be displayed circularly. Select this area and the technical alarm list will be displayed.

4. Physiological Alarm Area

This area shows physiological alarm messages. When multiple alarms occur, they will be displayed circularly. Select this area and the physiological alarm list will be displayed.

5. Waveform Area

This area shows measurement waveforms. The waveform name is displayed at the left upper corner of the waveform. Select this area and the corresponding measurement setup menu will be displayed.

6. Parameter Area A

This area shows measurement parameters. Each monitored parameter has a parameter window and the parameter name is displayed at the upper left corner. The corresponding waveform of each parameter is displayed in the same row in the waveform area. Select this area and the corresponding measurement setup menu will be displayed.

7. Parameter Area B

For the parameters displayed in this area, their corresponding waveforms are not displayed.

8. Prompt Message Area

This area shows the prompt messages, network status icons, battery status icons, date and time, etc. For details about battery status symbols, refer to the chapter *24 Batteries*.

- indicates patient monitor is connected to a wired network successfully.
- indicates the patient monitor has failed to connect a wired network.
- indicates the wireless function is working.
- indicates the wireless function is not working.
- ♦ indicates a USB disk is inserted.

9. QuickKeys Area

This area contains QuickKeys that give you fast access to functions.

2.5 QuickKeys

A QuickKey is a configurable graphical key, located at the bottom of the main screen. They give you fast access to functions. Their availability and the order in which they appear on your screen, depend on how your patient monitor is configured.

The following QuickKeys can be displayed on the screen:



Display more QuickKeys.



Hide the QuickKeys.



Enter the main menu



Enter standby mode



Change alarm settings



Review the patient's data



Enter the NIBP measurement menu



Stop all NIBP measurement



Zero IBP



Start the realtime print



Print Setup



Silence all alarm sounds



Pause or restore alarms



Change screen



Enter the patient setup menu



Trigger a manual event



Have a split-screen view of minitrends



Enter the volume setup menu



Default configurations

₩-	Start cardiac output procedure (not available in USA)
	Perform calculations
ήnή	Have a split-screen view of another patient's conditions
	Have a split-screen view of oxyCRG trends
	Enter the interpretation of resting 12-lead ECG screen (not available in USA)
~	Enter the full-screen 7-lead ECG screen
	Enter the [Parameters] menu
À	Start NIBP STAT measurement
	Enter the PAWP measurement screen
T	Enter the CPB mode
2	Enter the privacy mode
	Enter the night mode

You can also select your desired QuickKeys to display on the screen.

- Select [Main Menu]→[Maintenance >>]→[Manage Configuration >>]→enter the required password→[Ok].
- 2. In the [Manage Configuration] menu, select [Edit Config.>>].
- 3. In the pop-up menu, select the desired configuration and then select [**Edit**].
- 4. In the pop-up menu, select [Screen Setup >>].
- 5. In the [**Select QuickKeys**] screen, select your desired QuickKeys and the order of them.

3

Basic Operations

3.1 Installation



WARNING

- The equipment shall be installed by personnel authorized by us.
- The software copyright of the equipment is solely owned by us. No organization or individual shall resort to juggling, copying, or exchanging it or to any other infringement on it in any form or by any means without due permission.
- Devices connected to the equipment must meet the requirements of the applicable IEC standards (e.g. IEC 60950 safety standards for information technology equipment and IEC 60601-1 safety standards for medical electrical equipment). The system configuration must meet the requirements of the IEC 60601-1-1 medical electrical systems standard. Any personnel who connect devices to the equipment's signal input/output port is responsible for providing evidence that the safety certification of the devices has been performed in accordance to the IEC 60601-1-1. If you have any question, please contact us.
- If it is not evident from the equipment specifications whether a particular combination with other devices is hazardous, for example, due to summation of leakage currents, please consult the manufacturers or else an expert in the field, to ensure the necessary safety of patients and all devices concerned will not be impaired by the proposed combination.

3.1.1 Unpacking and Checking

Before unpacking, examine the packing case carefully for signs of damage. If any damage is detected, contact the carrier or us.

If the packing case is intact, open the package and remove the equipment and accessories carefully. Check all materials against the packing list and check for any mechanical damage. Contact us in case of any problem.

NOTE

Save the packing case and packaging material as they can be used if the equipment must be reshipped.



WARNING

- When disposing of the packaging material, be sure to observe the applicable waste control regulations and keep it out of children's reach.
- The equipment might be contaminated during storage and transport. Before use, please verify whether the packages are intact, especially the packages of single use accessories. In case of any damage, do not apply it to patients.

3.1.2 Environmental Requirements

The operating environment of the equipment must meet the requirements specified in this manual.

The environment where the equipment is used shall be reasonably free from noises, vibration, dust, corrosive, flammable and explosive substances. If the equipment is installed in a cabinet, sufficient space in front and behind shall be left for convenient operation, maintenance and repair. Moreover, to maintain good ventilation, the equipment shall be at least 2 inches (5cm) away from around the cabinet.

When the equipment is moved from one place to another, condensation may occur as a result of temperature or humidity difference. In this case, never start the system before the condensation disappears.



WARNING

 Make sure that the operating environment of the equipment meets the specific requirements. Otherwise unexpected consequences, e.g. damage to the equipment, could result.

3.2 Getting Started

3.2.1 Turning Power On

Once the patient monitor is installed, you can get ready for monitoring:

- 1. Before you start to make measurements, check the patient monitor and plug-in modules for any mechanical damage and make sure that all external cables, plug-ins and accessories are properly connected.
- 2. Plug the power cord into the AC power source. If you run the patient monitor on battery power, ensure that the battery is sufficiently charged.
- 3. Press the power on/off switch on the monitor's front.



WARNING

 Do not use the patient monitor for any monitoring procedure on a patient if you suspect it is not working properly, or if it is mechanically damaged. Contact your service personnel or us.

3.2.2 Starting Monitoring

- 1. Decide which measurements you want to make.
- 2. Connect the required modules, patient cables and sensors.
- 3. Check that the patient cables and sensors are correctly connected.
- 4. Check that the patient settings such as [Patient Cat.], [Paced], etc, are appropriate for your patient.
- 5. Refer to the appropriate measurement section for details of how to perform the measurements you require.

3.3 Disconnecting from Power

To disconnect the patient monitor from the AC power source, follow this procedure:

- 1. Confirm that the patient monitoring is finished.
- 2. Disconnect the patient cables and sensors from the patient monitor.
- 3. Make sure to save or clear the patient monitoring data as required.
- 4. Press and hold the power on/off switch. The patient monitor shuts down and you can unplug the power cable.



CAUTION

Although not recommended, you can press and hold the power on/off switch to forcibly shut down the
monitor when it could not be shut down normally or under some special situations. This may cause loss of
data of the patient monitor.

3.4 Using Keys

The monitor has three types of keys:


- Softkey: A softkey is a graphic key on the screen, giving you fast access to certain menus or functions. The monitor has three types of softkeys:
 - Parameter keys: Each parameter area can be seen as a softkey. You can enter a parameter setup menu by selecting its corresponding parameter area.
 - QuickKeys: QuickKeys are configurable graphical keys, located at the bottom of the main screen. For details, refer to the section *QuickKeys*.
- Hardkeys: A hardkey is a physical key on a monitoring device, such as the main menu hardkey on the monitor's front
- Pop-Up Keys: Pop-up keys are task-related keys that appear automatically on the monitor screen when required. For example, the confirm pop-up key appears only when you need to confirm a change.

3.5 Using the Touchscreen

Select screen items by pressing them directly on the patient monitor's screen. You can enable or disable touchscreen operation by pressing and holding the [Main Menu] QuickKey for 3 seconds. A padlock symbol is displayed if touchscreen operation is disabled.

3.6 Setting the Screen

You can enter the [Screen Setup] window as shown below by selecting [Main Menu] → [Screen Setup] → [Screen Layout >>]. In this window, you can allocate the positions of the parameters and waveforms. The parameters or waveforms whose positions are not allocated will not be displayed.



The ECG parameter and the first ECG waveform always display in the first row. The configurable areas can be classified as Area A, Area B, and Area C.

- In Area A, you can choose to display the parameters (having waveforms) and their waveforms. Each parameter and the associated waveform are displayed in the same row.
- In Area B, you can choose to display the parameters and their waveforms. When there is no parameter displayed in area C, both the parameters and their waveforms will be displayed in area B. Otherwise, only the parameters will be displayed.
- In Area C, you can choose to display all the parameters whose associated waveforms will not be displayed.

The screen can automatically adjust to ensure the best view based on your screen setup.

If no corresponding parameter or waveform is displayed after the module is inserted, you should perform the following inspections:

- Check the connection between the module and lead, cable, sensor, or external device.
- Enter the [**Screen Setup**] window for the desired display configuration.




WARNING

• The parameters whose positions are not allocated in the [Screen Setup] window will not be displayed. However, the monitor can still give alarms of these parameters.

3.7 Using the Main Menu

To enter the main menu, select the on-screen QuickKey or the hardkey on the monitor's front. Most of monitor operations and settings can be performed through the main menu.



Other menus are similar to the main menu and contain the following parts:

- 1. Heading: gives a sum-up for the current menu.
- 2. Main body: displays options, buttons, prompt messages, etc. The menu button with ">>" enlarges a secondary window to reveal more options or information.
- 3. Online help area: displays help information for the highlighted menu item.
- 4. X: select to exit the current menu.

3.8 Removing a Label Conflict

Every label is unique and is assigned only once. The measurement label is stored inside the module. The system will prompt module name conflict when two measurement modules with the same name are used.

For example, an IBP module (module X) is already loaded and the Art label is used for module X. Then another IBP module (module Y) is inserted and the Art label is also used for module Y. In this case, your patient monitor will prompt the message of label conflict and display the [**Label**] menu.

- To use module X for Art measurement, just modify the label of module Y on this channel in the [**Label**] menu. If the [**Label**] menu already exits inadvertently, you need to plug out and then plug in module Y.
- To use module Y for Art measurement, first exit the [**Label**] menu. Then select the Art parameter area on the screen and modify the label of module X on this channel in the popup menu. Finally, plug out and then plug in module Y.

3.9 Changing General Settings

This chapter covers only general settings such as language, brightness, date and time, etc. Measurement settings and other settings can be referred to in respective sections.

3.9.1 Setting up a Monitor

In situations where you install a patient monitor or change the patient monitor's application site, you need to setup the patient monitor as follows:

- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. In the [User Maintenance] menu, select, in turn, [Monitor Name], [Department] and [Bed No.], and then change their settings.

3.9.2 Changing Language

- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. In the [User Maintenance] menu, select [Language] and then select the desired language.
- 3. Restart the patient monitor.

3.9.3 Adjusting the Screen Brightness

- 1. Select the [Main Menu]→[Screen Setup >>]→[Brightness].
- 2. Select the appropriate setting for the screen brightness. 10 is the brightest, and 1 is the least bright.

If the patient monitor operates on battery power, you can set a less bright screen to prolong the operating time of the battery. When the patient monitor enters standby mode, the screen will change to the least brightness automatically.

3.9.4 Showing/Hiding the Help

The patient monitor provides online help information. The user can display or hide the help as required.

- 1. Select [Main Menu]→[Screen Setup >>].
- 2. Select [**Help**] and toggle between [**On**] and [**Off**].

3.9.5 Setting the Date and Time

- 1. Select [Main Menu] → [Maintenance >>] → [System Time >>].
- 2. Set the date and time.
- 3. Select [Date Format] and toggle between [yyyy-mm-dd], [mm-dd-yyyy] and [dd-mm-yyyy].
- 4. Select [**Time Format**] and toggle between [**24h**] and [**12h**].

If your patient monitor is connected to a central monitoring system (CMS), the date and time are automatically taken from that CMS. In that case, you cannot change the date and time settings on your patient monitor.



CAUTION

Changing date and time will affect the storage of trends and events and may cause data missing.

3.9.6 Adjusting Volume

Alarm Volume

- 1. Select the [Volume Setup] QuickKey, or [Main Menu]→[Alarm Setup >>]→[Others].
- 2. Select [**Alm Volume**] and then select the appropriate volume: X-10, in which X is the minimum volume, depending on the set minimum alarm volume (refer to the chapter Alarm), and 10 the maximum volume.

Key Volume

When you press the navigation knob or the touchscreen, or the hardkeys on the panel, the patient monitor prompts you by making a sound of the key volume you have set.

- 1. Select the [Volume Setup] QuickKey, or [Main Menu]→[Screen Setup >>].
- 2. Select [Key Volume] and then select the appropriate volume. 0 means off, and 10 the maximum volume.

QRS Volume

The QRS tone is derived from either the HR or PR, depending on which is currently selected as the alarm source in [**ECG Setup**] or [**SpO2 Setup**]. When monitoring SpO₂, there is a variable pitch tone which changes as the patient's saturation level changes. The pitch of the tone rises as the saturation level increases and falls as the saturation level decreases. The volume of this tone is user adjustable.

- Select the [Volume Setup] QuickKey, or the ECG parameter window→[Others >>], or the SpO₂ parameter window.
- 2. Select [**QRS Volume**] or [**Beat Vol**] and then select the appropriate volume. 0 means off, and 10 the maximum volume.

3.10 Operating Modes

Your monitor has different operating modes. Some are password protected. This section lists the major operating modes.

3.10.1 Monitoring Mode

This is the normal, everyday working mode that you use for monitoring patients. Your monitor automatically enters the monitoring mode after being turned on.

3.10.2 Night Mode

To avoid disturbing the patient, night mode may be used.

To activate the night mode:

- 1. Select the [Night Mode] QuickKey or [Main Menu]→[Screen Setup >>]→[Night Mode >>].
- 2. In the pop-up menu, set the desired brightness, alarm volume, QRS volume, key volume, or whether to stop NIBP measurement or not. When [**Stop NIBP**] is selected, all the NIBP measurements terminate after entering the night mode.
- 3. Select the [Enter Night Mode] button.

To cancel the night mode:

- Select the [Night Mode] QuickKey or [Main Menu]→[Screen Setup >>]→[Night Mode >>].
- 2. Select [**Ok**] in the popup.



WARNING

 Before entering night mode, confirm the settings of brightness, alarm volume, QRS volume, and key volume. Pay attention to the potential risk when the setting value is a bit low.

3.10.3 Privacy Mode

Privacy mode is only available when a patient who is admitted at a patient monitor is also monitored by the central station.

To activate the privacy mode, select [Main Menu]→[Screen Setup >>]→[Privacy Mode].

The patient monitor behaves as follows as soon as the privacy mode is activated:

- The screen turns blank and [Under monitoring. Press any key to exit the privacy mode.] is displayed.
- Monitoring and data storing continue but patient data is only visible at the central station.
- Alarms can still be triggered. But all audible alarms are suppressed and the alarm light is deactivated at the patient monitor.
- All system sounds are suppressed, including heart beat tone, pulse tone, all prompt tones, etc.

To cancel the privacy mode, press any key.

The patient monitor exits the privacy mode automatically in one of the following situations:

- The patient monitor disconnects from central station.
- Plug or unplug the storage card.

The touchscreen is locked automatically in the privacy mode.

WARNING

During privacy mode, all audible alarms are suppressed and the alarm light is deactivated at the patient monitor. Alarms sound only at the central station.

3.10.4 Demo Mode

In Demo mode, the monitor can demonstrate its major functions when patient or patient simulator is not connected. The Demo mode is password protected.

To enter the Demo mode:

- 1. Select [Main Menu]→[Maintenance >>].
- Select [Exit Demo].

To exit the Demo mode, select [Main Menu] \rightarrow [Maintenance >>] \rightarrow [Exit Demo].



⚠ WARNING

The Demo mode is for demonstration purpose only. To avoid that the simulated data are mistaken for the monitored patient's data, you must not change into Demo mode during monitoring. Otherwise, improper patient monitoring and delayed treatment could result.

3.10.5 Standby Mode

In standby mode, you can temperately stops patient monitoring without turning off the monitor. To enter the standby mode, Select [Main Menu]→[Standby].

FOR YOUR NOTES		

4 Managing Patients

4.1 Admitting a Patient

The patient monitor displays physiological data and stores them in the trends as soon as a patient is connected. This allows you to monitor a patient that is not admitted yet. However, it is recommended that you fully admit a patient so that you can clearly identify your patient, on recordings, reports and networking devices.

To admit a patient:

- 1. Select the [Patient Setup] QuickKey, or [Main Menu]→[Patient Setup >>].
- 2. Select [**Discharge Patient**] to clear any previous patient data. If you do not erase data from the previous patient, the new patient's data will be saved into the data of the previous patient. The monitor makes no distinction between the old and the new patient data.
- 3. If [Discharge Patient] button appears dimmed, directly select [Admit Patient] and then select:
 - ♦ [Yes] to apply the data saved in the patient monitor to the new patient, or
 - ◆ [No] to clear the data saved in the patient monitor.
- 4. In the [Patient Demographics] menu, enter the demographic details, of which:
 - [Patient Cat.] determines the way your patient monitor processes and calculates some measurements, and what safety and alarm limits are applied for your patient.
 - ◆ [Paced] determines whether to show pace pulse marks on the ECG waveform. When the [Paced] is set to [No], pace pulse marks are not shown in the ECG waveform.
- 5. Select [**Ok**].

WARNING

- [Patient Cat.] and [Paced] will always contain a value, regardless of whether the patient is fully admitted or not. If you do not specify settings for these fields, the patient monitor uses the default settings from the current configuration, which might not be correct for your patient.
- For paced patients, you must set [Paced] to [Yes]. If it is incorrectly set to [No], the patient monitor could mistake a pace pulse for a QRS and fail to alarm when the ECG signal is too weak.
- For non-paced patients, you must set [Paced] to [No].

4.2 Quick Admitting a Patient

Use [**Quick Admit**] only if you do not have the time or information to fully admit a patient. Complete the rest of the patient demographic details later. Otherwise, the ? symbol will always be displayed in the patient information area.

- 1. Select the [Patient Setup] QuickKey, or [Main Menu]→[Patient Setup >>].
- 2. Select [**Quick Admit**]. If a patient has been admitted at present, select [**OK**] to discharge the current patient. If .no patient is admitted, you can choose either:
 - ◆ [Yes] to apply the data in your patient monitor to the new patient, or
 - ◆ [No] to clear any previous patient data.
- 3. Enter the patient category and paced status for the new patient, and then select $[\mathbf{0k}]$.

4.3 Editing Patient Information

To edit the patient information after a patient has been admitted, or when the patient information is incomplete, or when you want to change the patient information:

- 1. Select the [Patient Setup] QuickKey, or [Main Menu]→[Patient Setup >>].
- 2. Select [Patient Demographics] and then make the required changes.
- 3. Select [**Ok**].

4.4 Discharging a Patient

To discharge a patient:

- 1. Select the [Patient Setup] QuickKey, or [Main Menu]→[Patient Setup >>].
- 2. Select [Discharge Patient]. In the popup menu, you can either:
 - ◆ Directly select [**Ok**] to discharge the current patient, or
 - ◆ Select [**Standby**] then [**Ok**]. The patient monitor enters the standby mode after discharging the current patient, or
 - ◆ Select [**Cancel**] to exit without discharging the patient.

NOTE

• Discharging a patient clears all history data in the monitor.

4.5 Transferring a Patient

You can transfer patient data between monitors with a USB drive without re-entering the patient demographic information. Transferring of patient data enables you to understand the patient's history condition. The patient data that can be transferred includes: patient demographics, trend data, alarm events and parameters alarm limits. USB drive can be used for patient data transfer of BeneView T Series patient monitors, iPM patient monitors, and iMEC patient monitors.

Select [Others >>] from [User Maintenance] menu. In the popup menu, you can set [Transferred Data Length]. The default is [4 h]. You can also set [Data Transfer Method]. The default is [Off].

4.5.1.1 Transferring Data from the Monitor to a USB Drive

- 1. Select [Main Menu]→[Patient Setup >>].
- 2. Select [**Transfer to Storage Medium**]. In the popup menu, you can:
 - ◆ Select [**Ok**] to transfer the patient data, or
 - ♦ Select [**Cancel**] to exit the menu.
- 3. Wait until the following message appears: [Transfer to storage medium successful. Please remove the USB drive.].
- 4. Remove the USB drive from patient monitor.

4.5.1.2 Transferring Data from a USB Drive

- 1. Connect the USB drive to the patient monitor.
- 2. In the popup menu, you can:
 - ◆ Select [**Transfer**] to transfer the patient data to the monitor, or
 - Select [Cancel Transfer] to cancel the operation of transferring patient data.
 - ♦ Select [**Unload USB Drive**] to not transfer the patient data and to unload the USB drive.
- 3. After you select [**Transfer**], in the popup menu you can further select the patient data contents that need to be transferred. [**Patient Demographics**] must be selected. After [**Ok**] is selected, the monitor compares the patient information stored in both the storage medium and monitor and deals with the patient data based on the following.
 - ◆ Different The monitor erases all the current patient data, transfers the patient data from the storage Patients: medium, and loads the configuration according to the patient category.
 - ◆ Same Patient: In the popup dialog box, you can:
 - ◆ Select [Yes] to merge the patient data in the monitor and storage medium.
 - ◆ Select [**No**] to erase all the current patient data in the monitor and to transfer the patient data from the storage medium.
- 4. Wait until the following message appears:[**Transfer from storage medium successful.**].

MARNING

- The USB drive you use may have write-protect function. In this case, please make sure the USB drive for data transfer is in read/write mode.
- Do not remove the storage medium during data transfer process. Otherwise, data files may be damaged.
- Do not discharge a patient before the patient is successfully transferred.
- After a patient is successfully transferred, check if the patient settings (especially patient category, paced status and alarm limits settings, etc) on the monitor are appropriate for this patient.

4.6 Connecting to a Central Monitoring System

If your patient monitor is connected to a central monitoring system (CMS):

- All patient information, measurement data and settings on the patient monitor can be transferred to the CMS.
- All patient information, measurement data and settings can be displayed simultaneously on the patient monitor and CMS. For some functions such as editing patient information, admitting a patient, discharging a patient, starting/stopping NIBP measurements, etc., bi-directional control can be achieved between your patient monitor and the CMS.

For details, refer to the CMS's instructions for use.

5

Managing Configurations

5.1 Introduction

When performing continuous monitoring on a patient, the clinical professional often needs to adjust the monitor's settings according to the patient's condition. The collection of all these settings is called a configuration. Allowing you to configure the monitor more efficiently, the monitor offers different sets of configuration to suit different patient categories and departments. You can change some settings from a certain set of configuration and then save the changed configuration as a user configuration.

The default configurations provided for your monitor are department-oriented. You can choose either from:

- General
- OR
- ICU
- NICU
- CCU

Each department has three different sets of configurations tailored for adult, pediatric and neonatal patients.



• The configuration management function is password protected. The configuration management tasks must be performed by clinical professionals.

The system configuration items can be classified as:

Parameter configuration items

These items relates to parameters, e.g., waveform gain, alarm switch, alarm limits..

■ Conventional configuration items

These items define how the monitor works, e.g., screen layout, record, print and alarm settings.

User maintenance items

These items relates to user maintenance settings, e.g., unit setup, time format and data format.

For the important configuration items and their default values and user maintenance items, see appendix *Configuration*Default Information.

5.2 Entering the [Manage Configuration] Menu

- 1. Press the hardkey on the monitor's front to enter the main menu.
- 2. Select [Maintenance >>]→[Manage Configuration >>]. Enter the required password and then select [Ok].



5.3 Changing Department

If the current department configuration is not the one you want to view, you can select [**Change Department >>**] in the [**Manage Configuration**] menu and then choose the one you want for viewing as shown below.



NOTE

• Changing the department will delete all current user configurations. Please act with caution.

5.4 Setting Default Configuration

The monitor will load the pre-set default configuration in the following cases.

- The patient monitor restarts after quitting over 120 seconds.
- A patient is admitted.
- A patient is discharged.
- Patient data is cleared.
- Patient category is changed.

To set default configuration:

- 1. Select [**Select Default Config.** >>] in the [Manage Configuration] menu.
- 2. In the [Select Default Config.] menu, select [Load the Latest Config.] or [Load Specified Config.].

When you select [**Load Specified Config.**], the configuration (adult, pediatric or neonate) to be restored is subject to the patient category. This configuration can be either factory configuration or saved user configuration. Take adult as an example, select [**Default Adu Config.**] and toggle between [**Defaults**] or user configuration(s).

NOTE

 To know what configuration is restored when the patient monitor starts, enter the main screen to check the prompt information at the lower part of the screen (displayed for about 10 seconds).


5.5 Saving Current Settings

Current settings can be saved as user configuration. Up to 5 user configurations can be saved. To save current settings:

- 1. Select [Save Current Settings As >>] in the [Manage Configuration] menu.
- 2. In the popup dialog box, enter the configuration name and then select [**Ok**].

5.6 Editing Configuration

1. Select [**Edit Config. >>**] in the [**Manage Configuration**] menu. The following menu appears.



2. The popup menu shows the existing configurations on the monitor. Selecting [**Config. on USB drive >>**] will show the existing configurations on the USB drive. Select the desired configuration and then select the [**Edit**] button. The following menu appears.



- 3. Select [**Alarm Setup >>**], [**Screen Setup >>**] or [**Parameter >>**] to enter the corresponding menu in which settings can be changed. The changed items of alarm setup will be marked in red.
- 4. You can select [**Save**] or [**Save as**] to save the changed configuration. Select [**Save**] to overwrite the original configuration. Select [**Save as**] to save the changed configuration in another name.

5.7 Deleting a Configuration

- 1. Select [**Delete Config. >>**] in the [**Manage Configuration**] menu.
- 2. The popup menu shows the existing user configurations on the monitor. Selecting [**Config. on USB drive** >>] will show the existing user configurations on the USB drive. Select the user configurations you want to delete and then select [**Delete**].
- 3. Select [**Yes**] in the popup.

5.8 Transferring a Configuration

When installing several monitors with identical user configuration it is not necessary to set each unit separately. An USB drive may be used to transfer the configuration from monitor to monitor.

To export the current monitor's configuration:

- 1. Connect the USB drive to the monitor's USB port.
- Select [Export Config. >>] in the [Manage Configuration] menu.
 In the [Export Config.] menu, select the configurations and [User Maintenance Settings] to export. Then select the [Export] button. A status message will report completion of the transfer.

To import the configuration on the USB drive to the monitor:

- 1. Connect the USB drive to the monitor's USB port.
- 2. Select [Import Config. >>] in the [Manage Configuration] menu.
- 3. In the [Import Config.] menu, select the configurations and [User Maintenance Settings] to import. Then select the [Import] button. A status message will report completion of the transfer.

5.9 Loading a Configuration

You may make changes to some settings during operation. However, these changes or the pre-selected configuration may not be appropriate for the newly admitted patient. Therefore, the monitor allows you to load a desired configuration so as to ensure that all the settings are appropriate for your patient.

To load a configuration,

- 1. Select [**Load Configuration** >>] from the main menu.
- 2. The popup menu shows the existing configurations on the monitor. Selecting [**Config. on USB drive >>**] will show the existing configurations on the USB drive.
- 3. Select a desired configuration.
- 4. Select [View] to view the configuration details. In the popup menu, you can select [**Alarm Setup >>**], [**Screen Setup >>**] or [**Parameter >>**] to view the corresponding contents. The alarm setup items which are different than those currently used are marked in red.
- 5. Select [**Load**] to load this configuration.

5.10 Restoring the Latest Configuration Automatically

During operation, you may make changes to some settings. However, these changes may not be saved as user configuration. To prevent the changes from losing in case of a sudden power failure, the patient monitor stores the configuration in real time. The saved configuration is the latest configuration.

The monitor restore the latest configuration if restarts within 60 seconds after the power failure. And it will restore the default configuration rather than the latest configuration if restarts 120 seconds later after the power failure. The monitor may load either the latest configuration or the default configuration if restarts from 60-120 seconds after the power failure.

5.11 Modifying Password

To modify the password for accessing the [Manage Configuration] menu,

- 1. Select [Modify Password >>] in the [Manage Configuration] menu.
- 2. Input a new password in the popup menu.
- 3. Select [**Ok**].

6 User Screens

6.1 Tailoring Your Screens

You can tailor your patient monitor's screens by setting:

- Wave line size
- The color in which each measurement's numerics and waveform are displayed
- The parameter to be monitored.

Changing some settings may be hazardous. Therefore, those setting are password-protected and can be modified by authorized personnel only. Once change is made, those who use the patient monitor should be notified.

6.1.1 Changing the Wave Line Size

- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. Select [Others >>].
- 3. Select [Wave Line] and toggle between [Thick], [Mediate] and [Thin].

6.1.2 Changing Measurement Colors

- 1. Select [Main Menu]→[Screen Setup >>]→[Measurement Color Setup >>].
- 2. Select the color box next to your desired measurement and then select a color from the popup menu.

6.1.3 Changing Screen Layout

Select the [Screens] QuickKey, or [Main Menu]→[Screen Setup >>]→[Screen Layout >>] to enter the [Screens] menu.

- You can choose the desired screen type in the [**Choose Screen**] window.
- You can select the parameters and waveforms you want to view in the [Screen Setup] window. For details, please refer to the section Setting the Screen.
- You can select the parameters you want to view on big numerics screen in the [**Big Numerics Screen Setup**] window.


6.2 Viewing Minitrends

6.2.1 Having a Split-Screen View of Minitrends

You can split the normal screen so that one part of the screen, on the left hand side, continuously shows graphic minitrends beside waveforms as shown in the figure below.

To have a split-screen view of minitrends, you can:

- Select [Minitrends] QuickKey, or
- Select [Screens] QuickKey→[Minitrends Screen]→X, or
- Select [Main Menu]→[Screen Setup >>]→[Screen Layout >>]→[Minitrends Screen]→X.



The split-screen view provides minitrends for multiple parameters. In each field, the label, scale and time are respectively displayed at the top, left, and bottom as shown below.



6.2.2 Setting Minitrends


Select the minitrends area. From the pop-up [Minitrend Setup] menu, you can:

- Select the parameters to be displayed, or
- Select [Minitrend Length] and then select the appropriate setting.

6.3 Viewing oxyCRG

To have a split screen view of oxyCRG, you can:

- Select [**oxyCRG**] QuickKey, or
- Select [Screens] QuickKey→[OxyCRG Screen]→X, or.
- Select [Main Menu]→[Screen Setup >>]→[Screen Layout >>]→[OxyCRG Screen]→X.



The split-screen view covers the lower part of the waveform area and shows HR trend, SpO_2 trend and RR trend (or Resp wave). At the bottom, there are controls:

1. Trend length list box

In the trend length list box, you can select [1 min], [2 min], [4 min], or [8 min].

2. Resp Wave (or RR Trend) list box

From this list box, you can select either [**Resp Wave**] or [**RR Trend**] for display.

3. Record

Through this button, you can print out the currently displayed oxyCRG trends by the recorder.

6.4 Viewing Other Patients

6.4.1 Care Group

You can select up to 10 patient monitors (including telemetry) connected to the same central monitoring system into a Care Group. This lets you:

- View information on the monitor screen from another bed in the same Care Group.
- Be notified of physiological and technical alarm conditions at the other beds in the same Care Group.

To have a Care Group:

- 1. Open the [View Other Patient] window by:
 - ◆ Selecting [**Others**] QuickKey, or
 - ◆ Selecting [Screens] QuickKey→[View Others Screen]→ X, or
 - ◆ Selecting [Main Menu]→[Screen Setup >>]→[Screen Layout >>]→[View Others Screen]→X.
- 2. Select [Setup] in the [View Other Patient] window.
- 3. Select the desired patient monitors from the [**Connected Monitor List**], and then select the button. The selected patient monitors constitute a Care Group.

6.4.2 Viewing the Care Group Overview Bar



The Care Group overview bar locates at the bottom of the [**View Other Patient**] window. In the overview bar, the department and bed label for any Care Group beds are displayed. For telemetry, # is displayed before the department label. The color in which a Care Group bed appears matches its status:


- Red: indicates the bed is giving high-level physiological alarms or the telemetry is giving alarm, such as nurse call or event.
- Yellow: indicates the bed is giving medium- or low-level physiological alarms, or medium-level technical alarms.
- Blue: indicates the bed is giving low-level technical alarms.
- Grey: indicates the bed fails to be networked or stays in the standby mode.

You can view a Care Group bed's alarms by selecting it from the care group, and as well you can select the [**View This Patient**] button to view this bed in the [**View Other Patient**] window.

For more details about Care Group alarms, refer to the *Alarms* chapter.

6.4.3 Understanding the View Other Patient Window

When you first open the [**View Other Patient**] window, the patient monitor automatically selects a monitor from the network to display in the [**View Other Patient**] window.



The [View Other Patient] window covers the lower part of the waveform area and consists of:

- 1. Information Area: shows the patient information (including department, bed number, patient name, etc.), network status symbol.
- View Area: shows physiological waveforms and parameters. You can switch a waveform area to a parameter area by selecting your desired waveform area and then selecting [Switch to Parameter Area], or switch a parameter area to a waveform area by selecting your desired parameter area and then selecting [Switch to Waveform Area].
- 3. Care Group Overview Bar.
- 4. Message Area: shows physiological, technical and prompt messages from the currently viewed patient monitor. It also shows the alarm given by the telemetry such as nurse call or event. By selecting this area, you can enter the [Alarm Information List] to view all physiological, technical and prompt messages coming from the currently viewed patient.

Additionally, you can change a waveform or parameter for viewing

- To change a waveform for viewing, select the waveform segment where you want a new waveform to appear and then select the waveform you want from the popup menu.
- To change a parameter for viewing, select the parameter window where you want a new parameter to appear and then select the parameter you want from the popup menu.




WARNING

• The data presented in the [View Other Patient] window have delay. Do not rely on this window for realtime data.

6.5 Understanding the Big Numerics Screen

To enter the big numerics screen:

- 1. Select the [Screens] QuickKey, or [Main Menu]→[Screen Setup >>]→[Screen Layout >>].
- 2. Select [**Big Numerics**]→ X.



You can select your desired parameters to display in this screen: select the [**Screens**] QuickKey→[**Big Numerics Screen Setup**] and then select the parameters you want. For parameters having a waveform, the waveform will also be displayed.

7 Alarms

Alarms, triggered by a vital sign that appears abnormal or by technical problems of the patient monitor, are indicated to the user by visual and audible alarm indications.

NWARNING

- A potential hazard can exist if different alarm presets are used for the same or similar equipment in any single area, e.g. an intensive care unit or cardiac operating room.
- If your patient monitor is connected to a CMS, remote suspension, inhibition, silence and reset of monitor alarms via the CMS may cause a potential hazard.

7.1 Alarm Categories

By nature, the patient monitor's alarms can be classified into three categories: physiological alarms, technical alarms and prompt messages.

1. Physiological alarms

Physiological alarms, also called patient status alarms, are triggered by a monitored parameter value that violates set alarm limits or an abnormal patient condition. Physiological alarm messages are displayed in the physiological alarm area.

2. Technical alarms

Technical alarms, also called system status alarms, are triggered by a device malfunction or a patient data distortion due to improper operation or mechanical problems. Technical alarm messages are displayed in the technical alarm area.

3. Prompt messages

As a matter of fact, prompt messages are not alarm messages. Apart from the physiological and technical alarm messages, the patient monitor will show some messages telling the system status or patient status. Messages of this kind are included into the prompt message category and usually displayed in the prompt information area. Some prompt messages that indicate the arrhythmia events are displayed in the physiological alarm area. For some measurements, their related prompt messages are displayed in their respective parameter windows.

7.2 Alarm Levels

By severity, the patient monitor's alarms can be classified into three categories: high level, medium level and low level..

	Physiological alarms	Technical alarms
High level	Indicate that your patient is in a life	Indicate a severe device malfunction or an improper operation,
	threatening situation, such as Asystole,	which could make it possible that the monitor cannot detect
	Vfib/Vtac and so forth, and an	critical patient status and thus threaten the patient's life, such as
	emergency treatment is demanded.	low battery.
Medium	Indicate that your patient's vital signs	Indicate a device malfunction or an improper operation, which
level	appear abnormal and an immediate	may not threaten the patient's life but may compromise the
	treatment is required.	monitoring of vital physiological parameters.
Low level	Indicate that you patient's vital signs	Indicate a device malfunction or an improper operation, which
	appear abnormal and an immediate	may compromise a certain monitoring function but will not
	treatment may be required.	threaten the patient's life.

7.3 Alarm Indicators

When an alarm occurs, the patient monitor will indicate it to the user through visual or audible alarm indications.

- Alarm lamp
- Alarm message
- Flashing numeric
- Audible alarm tones

7.3.1 Alarm Lamp

If a technical alarm or physiological alarm occurs, the alarm lamp will flash. The flashing color and frequency match the alarm level as follows:

High level alarms: the lamp quickly flashes red.Medium level alarms: the lamp slowly flashes yellow.

■ Low level alarms: the lamp turns yellow without flashing.

7.3.2 Alarm Message

When an alarm occurs, an alarm message will appear in the technical or physiological alarm area. For physiological alarms, the asterisk symbols (*) before the alarm message match the alarm level as follows:

■ High level alarms: ***

■ Medium level alarms: **

■ Low level alarms: *

Additionally, the alarm message uses different background color to match the alarm level:

High level alarms: redMedium level alarms: yellowLow level alarms: yellow

You can view the alarm messages by selecting the physiological or technical alarm area.

7.3.3 Flashing Numeric

If an alarm triggered by an alarm limit violation occurs, the numeric of the measurement in alarm will flash every second, and the corresponding alarm limit will also flash at the same frequency indicating the high or low alarm limit is violated.

7.3.4 Audible Alarm Tones

The alarm tone is distinct from heart beat tone, keystroke tone and pulse tone in frequency. This monitor has three choices of alarm tones and patterns: ISO, Mode 1 and Mode 2. For each pattern, the alarm tones identify the alarm levels as follows:

■ ISO pattern:

♦ High level alarms: triple+double+triple+double beep.

Medium level alarms: triple beep.Low level alarms: single beep.

■ Mode 1:

♦ High level alarms: high-pitched single beep.

◆ Medium level alarms: double beep.

♦ Low level alarms: low-pitched single beep.

■ Mode 2:

♦ High level alarms: high-pitched triple beep.

Medium level alarms: double beep.

◆ Low level alarms: low-pitched single beep.

NOTE

 When multiple alarms of different levels occur simultaneously, the patient monitor will select the alarm of the highest level and give visual and audible alarm indications accordingly.

7.3.5 Alarm Status Symbols

Apart from the aforementioned alarm indicators, the patient monitor still uses the following symbols telling the alarm status:

- indicates alarms are paused.
- indicates alarm sound is silenced.
- indicates the alarm sound is turned off.
- indicates individual measurement alarms are turned off or the system is in alarm off status.

7.4 Alarm Tone Configuration

7.4.1 Setting the Minimum Alarm Volume

- Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. Select [Alarm Setup >>] to enter the [Alarm Setup] menu.
- 3. Select [Minimum Alarm Volume] and toggle between 0 and 10.

The minimum alarm volume refers to the minimum value you can set for the alarm volume, which is not affected by user or factory default configurations. The setting of minimum alarm volume remains unchanged when the patient monitor shuts down and restarts.

7.4.2 Changing the Alarm Volume

- Select the [Volume Setup] QuickKey or the [Alarm Setup] QuickKey→[Others], or [Main Menu]→[Alarm Setup >>]→[Others].
- 2. Select the appropriate volume from [**Alm Volume**]: X to 10, in which X is the minimum volume, depending on the set minimum alarm volume, and 10 the maximum volume.
- 3. Select [**High Alarm Volume**] to set the volume of the high priority alarm as [**Alm Volume+0**], [**Alm Volume+1**] or [**Alm Volume+2**].
- 4. Select [**Reminder Vol**] to set the volume of the reminder tone as [**High**], [**Med**] or [**Low**].

When alarm volume is set to 0, the alarm sound is turned off and a prears on the screen.

7.4.3 Setting the Interval between Alarm Sounds

You cannot change the interval between alarm tones if you choose mode 1 or 2 as your desired alarm tone pattern. For these two patterns, the interval between alarm tones identifies the alarm levels as follows:

- Mode 1:
 - ♦ Interval between high level alarm tones: continuously.
 - ♦ Interval between medium level alarm tones: 5 s.
 - ♦ Interval between low level alarm tones: 20 s.
- Mode 2:
 - ◆ Interval between high level alarm tones: 1 s.
 - ◆ Interval between medium level alarm tones: 5 s.
 - Interval between low level alarm tones:
 20 s.

If you choose the ISO pattern, you can change the interval between alarm tones. To change the interval between alarm tones:

- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. Select [Alarm Setup >>] to enter the [Alarm Setup] menu.
- 3. Select [**High Alarm Interval (s)**], [**Med Alarm Interval (s)**] and [**Low Alarm Interval (s)**] in turn and then select the appropriate settings.



🖺 WARNING

- When the alarm sound is switched off, the patient monitor will give no audible alarm tones even if a new alarm occurs. Therefore the user should be very careful about whether to switch off the alarm sound or not.
- Do not rely exclusively on the audible alarm system for patient monitoring. Adjustment of alarm volume to a low level may result in a hazard to the patient. Always keep the patient under close surveillance.

7.4.4 Changing the Alarm Tone Pattern

To change the alarm tone pattern:

- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. Select [Alarm Setup >>] to enter the [Alarm Setup] menu.
- 3. Select [Alarm Sound] and toggle between [ISO], [Mode 1] and [Mode 2].

User or factory default configurations exert no impact on the setup of alarm tone pattern. The alarm tone pattern remains unchanged after the monitor restarts.

7.4.5 Setting the Reminder Tones

When the alarm volume is set to zero, or the alarm tone is silenced or turned off, the patient monitor issues a periodical reminder tone. To set the reminder tones:


- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. Select [Alarm Setup >>] to enter the [Alarm Setup] menu.
 - ◆ To switch the reminder tones on or off, select [Reminder Tones] and toggle between [On] and [Off].
 - ◆ To set the interval between reminder tones, select [**Reminder Interval**] and toggle between []**1min**], [**2min**] and [**3min**].

In addition, you can set the volume of alarm reminder tones. To set the volume of alarm reminder tones, select [Main Menu]→[Alarm Setup >>]→[Others] or the [Alarm Setup] QuickKey→[Others]. Then, select [Reminder Vol] and toggle between [High], [Medium] and [Low].

7.5 Understanding the Alarm Setup Menu

Select the [**Alarm Setup**] QuickKey or [**Main Menu**]→[**Alarm Setup** >>] to enter the [**Alarm Setup**], where you can:

- Set alarm properties for all parameters.
- Change ST alarm settings.
- Change arrhythmia alarm settings.
- Set the threshold for some arrhythmia alarms.
- Change other settings.



Please refer to the *ECG* section for how to change ST alarm settings, how to change arrhythmia alarm settings and how to set the threshold for some arrhythmia alarms.

7.5.1 Setting Alarm Properties for All Parameters

In the main menu, select [**Alarm Setup >>**]→[**Parameters**]. You can review and set alarm limits, alarm switches, alarm level and alarm recordings for all parameters.

When a measurement alarm occurs, automatic recording of all the measurement numerics and related waveforms is possible when the measurement's [**On/Off**] and [**Record**] are set on.

- Make sure that the alarm limits settings are appropriate for your patient before monitoring.
- Setting alarm limits to extreme values may cause the alarm system to become ineffective. For example,
 High oxygen levels may predispose a premature infant to retrolental fibroplasia. If this is a consideration do NOT set the high alarm limit to 100%, which is equivalent to switching the alarm off.

7.5.2 Adjusting Alarm Limits Automatically

The monitor can automatically adjust alarm limits according to the measured vital signs, using the auto limits function. When auto limits are selected, the monitor calculates safe auto limits based on the latest measured values.

To get accurate auto alarm limits, you need to collect a set of measured vital signs as a baseline. Then, in the main menu, select [**Alarm Setup >>**] \rightarrow [**Parameters**] \rightarrow [**Auto Limits**] \rightarrow [**Ok**]. The monitor will create new alarm limits based on the measured values.

Before applying these automatically created alarm limits, confirm if they are appropriate for your patient in the mass alarm setup menu. If not, you can adjust them manually. These alarm limits will remain unchanged until you select auto limits again or adjust them manually.

The monitor calculates the auto limits based on the following rules.

		Low alarm limit		High alarm limit		
Module Pa	Parameter	Adult/ pediatric	Neonate	Adult/ pediatric	Neonate	Auto alarm limits range
ECG	HR/PR	(HR × 0.8) or 40bpm (whichever is greater)	(HR – 30) or 90bpm (whichever is greater)	(HR × 1.25) or 240bpm (whichever is smalle)	(HR + 40) or 200bpm (whichever is smaller)	Adult/pediatric: 35 to 240 Neonate: 55 to 225
Resp	RR	(RR × 0.5) or 6 rpm (whichever is greater)	(RR – 10) or 30 rpm (whichever is greater)	(RR × 1.5) or 30 rpm (whichever is smaller)		Adult/pediatric: 6 to 55 Neonate: 10 to 90
SpO₂	SpO ₂	Same as the default alarm limit	Same as the default alarm limit	Same as the default alarm limit	Same as the default alarm limit	Same as the measurement range
NIBP-S	(SYS × 0.68 + 10) mmHg	(SYS – 15) or 45mmHg (whichever is greater)	(SYS × 0.86 + 38) mmHg	(SYS + 15) or 105mmHg (whichever is smaller)	Adult: 45 to 270 Pediatric: 45 to 185 Neonate: 35 to 115	
NIBP	NIBP-D	(Dia × 0.68 + 6) mmHg	(Dia – 15) or 20mmHg (whichever is greater)	(Dia × 0.86 + 32) mmHg	(Dia + 15) or 80mmHg (whichever is smaller)	Adult: 25 to 225 Pediatric: 25 to 150 Neonate: 20 to 90
	NIBP-M	(Mean × 0.68 + 8) mmHg	(Mean – 15) or 35mmHg (whichever is greater)	(Mean × 0.86 + 35) mmHg	(Mean + 15) or 95 mmHg (whichever is smaller)	Adult: 30 to 245 Pediatric: 30 to 180 Neonate: 25 to 105
Temp	T1 T2	(T1 − 0.5) °C (T2 − 0.5) °C	(T1 − 0.5) °C (T2 − 0.5) °C	(T1 + 0.5)℃ (T2 + 0.5)℃	(T1 + 0.5) °C (T2 + 0.5) °C	1 to 49℃ 1 to 49℃

		Low alarm limit		High alarm limit		Auto alarm limits
Module	Parameter	Adult/ pediatric	Neonate	Adult/ pediatric	Neonate	range
	TD	Same as the default alarm limit	Same as the default alarm limit		Same as the default alarm limit	Same as the measurement range
IBP: ART/ Ao/ UAP/ BAP/	IBP-S	(SYS × 0.68+10) mmHg	(SYS – 15) or 45mmHg (whichever is greater)	(SYS × 0.86 + 38) mmHg	(SYS + 15) or 105mmHg (whichever is smaller)	Adult: 45 to 270 Pediatric: 45 to 185 Neonate: 35 to 115
	IBP-D	(Dia × 0.68+ 6)mmHg	(Dia – 15) or 20mmHg (whichever is greater)	(Dia × 0.86+ 32)mmHg	(Dia + 15) or 80mmHg (whichever is smaller)	Adult: 25 to 225 Pediatric: 25 to 150 Neonate: 20 to 90
P1-P4 (Arterial pressure)	ІВР-М	(Mean × 0.68 + 8)mmHg	(Mean – 15) or 35mmHg (whichever is greater)	(Mean × 0.86 + 35)mmHg	(Mean + 15) or 95mmHg (whichever is smaller)	Adult: 30 to 245 Pediatric: 30 to 180 Neonate: 25 to 105
	IBP-S	SYS × 0.75	SYS × 0.75	SYS × 1.25	SYS × 1.25	
PA	IBP-D	Dia × 0.75	Dia × 0.75	Dia × 1.25	Dia × 1.25	3 to 120mmHg
	IBP-M	Mean × 0.75	Mean × 0.75	Mean × 1.25	Mean × 1.25	
IBP: CVP/ ICP/ LAP/ RAP/ UVP/ P1-P4 (Venous pressure)	IBP-M	Mean × 0.75	Mean × 0.75	Mean × 1.25	Mean × 1.25	3 to 40mmHg
CO ₂		0 to 32mmHg: remains the same	0 to 32mmHg: remains the same	0 to 32mmHg: remains the same	0 to 32mmHg: remains the same	
		32 to 35mmHg: 29mmHg	32 to 35mmHg: 29mmHg	32 to 35mmHg: 41mmHg	32 to 35mmHg: 41mmHg	
	EtCO ₂	_	35 to 45mmHg: (etCO ₂ -6) mmHg	35 to 45mmHg: (etCO ₂ +6) mmHg	35 to 45mmHg: (etCO ₂ +6) mmHg	Same as the measurement range
		45 to 48mmHg:39 mmHg	45 to 48mmHg:39 mmHg	45 to 48mmHg:51 mmHg	45 to 48mmHg:51 mmHg	
		>48mmHg: remains the same	>48mmHg: remains the same	>48mmHg: remains the same	>48mmHg: remains the same	

		Low alarm limit		High alarm limit		Auto alarm limits		
Module	Parameter	Adult/ pediatric	Neonate	Adult/ pediatric	Neonate	range		
	FiCO ₂	N/A	N/A	Same as the default alarm limit	Same as the default alarm limit	Same as the measurement range		
	awRR	awRR × 0.5 or 6 rpm (whichever is greater)	(awRR – 10) or 30 rpm (whichever is greater)	awRR × 1.5 or 30 rpm (whichever is smaller)	rpm (whichever is	Adult/pediatric: 6 to 55 Neonate:10 to 90		
	EtCO ₂ (AG)							
	FiCO₂ (AG)	Same as CO₂ mod	dule					
AG	awRR	awRR × 0.5 or 6 rpm (whichever is greater)	awRR – 10 or 30 rpm (whichever is greater)	rpm (whichever is	rpm (whichever is	Adult/pediatric: 6 to 55 Neonate: 10 to 90		
	FiAA/ EtAA	Same as the default alarm limit	Same as the default alarm limit	Same as the default alarm limit	Same as the default alarm limit	Same as the measurement range		
	FiO ₂ / EtCO ₂	Same as the default alarm limit	Same as the default alarm limit	Same as the default alarm limit	Same as the default alarm limit	Same as the measurement range		
	FiN ₂ O/ EtN ₂ O	Same as the default alarm limit	Same as the default alarm limit	Same as the default alarm limit	Same as the default alarm limit	Same as the measurement range		
C.O.	ВТ	Adult: (BT – 1)°C	N/A	Adult: (BT – 1)℃	N/A	Same as the measurement range		

7.5.3 Setting Alarm Delay Time

You can set the alarm delay time for over-limit alarms of continuously measured parameters. If the alarm-triggered condition disappears within the delay time, the patient monitor will not give the alarm. You can set the [Alarm Delay], [Apnea Delay] and [ST Alarm Delay] in the [Others] window of [Alarm Setup] menu.

7.5.4 Setting Recording Length

You can change the length of the recorded waveforms. In the [Others] window of the [Alarm Setup] menu, select [Recording Length] and toggle between [8 s], [16 s] and [32 s]:

- [8 s]: 4 seconds respectively before and after the alarm or manual event trigger moment.
- [16 s]: 8 seconds respectively before and after the alarm or manual event trigger moment.
- [32 s]: 16 seconds respectively before and after the alarm or manual event trigger moment.

7.5.5 Entering CPB Mode

When performing Cardiopulmonary bypass (CPB), you can set the patient monitor to enter CPB mode in order to reduce unnecessary alarms. In the CPB mode, all the physiological alarms are switched off except for the following alarms.

- FiCO₂/EtCO₂ too high (for CO₂ module and AG module)
- FiO₂/EtO₂ too high or too low
- FiAA/EtAA too high (AA represents the anaesthetic gas)
- FiN₂O/EtN₂O too high

In CPB mode, [CPB Mode] is displayed in the physiological alarm area with red background color.

To enter CPB mode:

Select the [CPB Mode] Quickkey or select [Enter CPB Mode] in the [Others] window of the [Alarm Setup] menu. Then select [Ok] in the popup dialog box.

7.6 Pausing Alarms

If you want to temporarily prevent alarms from sounding, you can pause alarms by pressing the A hardkey on the monitor's front. When alarms are paused:

- No alarm lamps flash and no alarms are sounded.
- No numeric and alarm limit flash.
- No alarm messages are shown.
- The remaining pause time is displayed in the physiological alarm area.
- The 💹 alarms paused symbol is displayed in the sound symbol area.

If the time interval of the monitor's last shutdown from this starting-up is greater that 2 minutes, the patient monitor enters into the alarm paused status as soon as it is turned on. The alarm pause time is fixed to be 2 minutes.

When the alarm pause time expires, the alarm paused status is automatically cancelled and the alarm tone will sound. You can also cancel the alarm paused status by pressing the 🖄 hardkey.

You can set the alarm pause time to [1 min], [2 min], [3 min], [5 min], [10 min], [15 min], or [Permanent]. The default alarm pause time is 2 minutes.

- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. Select [Alarm Setup >>]→[Alarm Pause Time] and then select the appropriate setting from the popup list.

7.7 Switching Off All Alarms

If [**Alarm Pause Time**] is set to [**Permanent**], the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters into the alarm off status after the Alarm part in the patient monitor enters in the patient mo

- As for physiological alarms, no alarm lamps flash and no alarms are sounded.
- As for physiological alarms, no numeric and alarm limit flash.
- No physiological alarm messages are shown.
- [Alarm Off] is displayed in the physiological alarm area with red background.
- As for technical alarms, no alarms are sounded.
- The XX alarm off symbol is displayed in the sound symbol area.

You can cancel the alarm off status by pressing the 🖄 hardkey.



Pausing or switching off alarms may result in a hazard to the patient. Please be very careful.

7.8 Silencing the Alarm Sound

You can silence all alarm sounds by pressing the 🌋 hardkey on the monitor's front. In that case, the alarm lamp flashing and alarm tones are cleared and 🔯 appears in the sound symbol area. After the physiological alarm is silenced, ✓ appears before the alarm message and the numeric and alarm limit still flash. For the performance after the technical alarm is silenced, please refer to the *Silencing Technical Alarms* section.

The alarm silenced status will be automatically cancelled if you switch the patient monitor to other alarm statues or when a new physiological or technical alarm occurs.

7.9 Latching Alarms

The alarm latching setting for your patient monitor defines how the alarm indicators behave when you do not acknowledge them. When alarms are set to non-latching, their alarm indications end when the alarm condition ends. If you switch alarm latching on, all visual and audible alarm indications last until you acknowledge the alarms, except that the measurement numeric and violated alarm limit stop flashing as soon as the initial alarm condition goes away.

To set alarms to latching or non-latching:

- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. Select [Alarm Setup >>].
- 3. Select [Latching Alarms] and toggle between [High only], [Hi&Med], [All] and [Off]. If you select [High only], only high priority alarm are latched; if you select [Hi&Med], both high priority alarms and mediate priority alarms are latched; if you select [All], all alarms are latched; if you select [Off], the alarm latching is turned off.

NOTE

Changing of alarm priority may affect the latching status of corresponding alarm. Please determine if you
need to reset the latching status for the specific alarm when you have changed its alarm priority.

7.10 Silencing Technical Alarms

For some technical alarms, their alarm lamp flashing and alarm tones are cleared and the alarm messages change to prompt messages after the Alarm hardkey is pressed. After the patient monitor restores the normal status, the patient monitor can give alarm indications correctly when these alarms are triggered again.

For some technical alarms, all their alarm indications are cleared after the A hardkey is pressed. After the patient monitor restores the normal status, the patient monitor can give alarm indications correctly when these alarms are triggered again.

For some other technical alarms, their alarm lamp flashing and alarm tones are cleared and √ appears before the alarm message after the Ardkey is pressed. After the patient monitor restores the normal status, the patient monitor can give alarm indications correctly when these alarms are triggered again.

7.11 Testing Alarms

When the monitor starts up, a selftest is performed. In this case the alarm lamp is lit in yellow and red respectively, and the system gives a beep. This indicates that the visible and audible alarm indicators are functioning correctly.

For further testing of individual measurement alarms, perform the measurement on yourself (for example SpO_2 or CO_2) or use a simulator. Adjust alarm limits and check that appropriate alarm behaviour is observed.

7.12 When an Alarm Occurs

When an alarm occurs, observe the following steps and take proper actions:

- 1. Check the patient's condition.
- 2. Confirm the alarming parameter or alarm category.
- 3. Identify the source of the alarm.
- 4. Take proper action to eliminate the alarm condition.
- 5. Make sure the alarm condition is corrected.

For troubleshooting specific alarms, see appendix Alarm Messages.

7.13 Using Care Group Alarms

7.13.1 Care Group Auto Alarms

When auto alarm is set on for viewing other patient and a Care Group is set up on your monitor, a flashing symbol will appear beside the QuickKeys area if any monitor in your Care Group, which is not currently viewed by your monitor, is alarming. The alarm symbol is shown as below.



The department and bed label of the alarming monitor appear on the symbol. You can enter the view other patient window by pressing this symbol.

To switch auto alarm on or off for viewing other patient:

- 1. In the main menu, select [Screen Setup >>]→[Screen Layout >>]→[View Others Screen]→X.
- 2. In the view other patient window, select the care group setup button. Then, select [**Auto Alarm**] and toggle between [**On**] and [**Off**].

7.13.2 Silencing Care Group Alarms

You can silence the alarm sound of the currently viewed bed in the view other patient window. This function can be set in the [**Alarm Setup**] menu from the [**User Maintenance**] menu only.

When the alarm silence function for other patients is active and the currently viewed bed is in normal alarm status or alarm sound off status, press the [**Silence**] button in the view other patient window. The currently viewed bed will then enter into the alarm silenced status.

Note that this button is disabled when the currently viewed bed is in alarms off or paused status.



• Silencing care group alarms may cause a potential hazard. Please act with caution.

FOR YOUR NOTES		

8 Monitoring ECG

8.1 Introduction

The electrocardiogram (ECG) measures the electrical activity of the heart and displays it on the patient monitor as a waveform and a numeric. This section also tells you about ST monitoring and arrhythmia monitoring. 12-lead monitoring is available for iPM 12 and iPM 10 patient monitors only. ECG monitoring provides two algorithms:

- Mindray algorithm
- Mortara algorithm

The patient monitor incorporating Mortara algorithm is labelled with the logo of Mortara.

8.2 Safety



WARNING

- Use only ECG electrodes and cables specified by Mindray.
- When connecting electrodes and/or patient cables, make sure that the connectors never come into contact
 with other conductive parts, or with earth. In particular, make sure that all of the ECG electrodes are
 attached to the patient, to prevent them from contacting conductive parts or earth.
- Periodically inspect the electrode application site to ensure skin quality. If the skin quality changes,
 replace the electrodes or change the application site.
- Use defibrillator-proof ECG cables during defibrillation.
- Do not touch the patient, or table, or instruments during defibrillation.
- After defibrillation, the screen display recovers within 10 seconds if the correct electrodes are used and applied in accordance with the manufacturer's instructions for use.
- Interference from a non-grounded instrument near the patient and electrosurgery interference can cause problems with the waveform.

8.3 Preparing to Monitor ECG

8.3.1 Preparing the Patient and Placing the Electrodes

- 1. Prepare the patient's skin. Proper skin preparation is necessary for good signal quality at the electrode, as the skin is a poor conductor of electricity. To properly prepare the skin, choose flat areas and then follow this procedure:
 - ◆ Shave hair from skin at chosen sites.
 - ◆ Gently rub skin surface at sites to remove dead skin cells.
 - ◆ Thoroughly cleanse the site with a mild soap and water solution. We do not recommend using ether or pure alcohol, because this dries the skin and increases the resistance.
 - Dry the skin completely before applying the electrodes.
- 2. Attach the clips or snaps to the electrodes before placing them.
- 3. Place the electrodes on the patient.
- 4. Attach the electrode cable to the patient cable and then plug the patient cable into the ECG connector.

8.3.2 Choosing AHA or IEC Lead Placement

- 1. Select the ECG parameter window or waveform area to enter the [**ECG Setup**] menu.
- 2. Select [Lead Set] and then select [3-lead], [5-lead], [12-lead] or [Auto] according to the applied electrodes.
- Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password
- Select [Others >>]→[ECG Standard] and then select [AHA] or [IEC] according to the standard that is applied for your hospital.

8.3.3 ECG Lead Placements

The electrode placement illustrations in this chapter adopt the AHA standard.

3-Leadwire Electrode Placement


Following is an electrode configuration when using 3 leadwires:

- RA placement: directly below the clavicle and near the right shoulder.
- LA placement: directly below the clavicle and near the left shoulder.
- LL placement: on the left lower abdomen.

5-Leadwire Electrode Placement

Following is an electrode configuration when using 5 leadwires:

- RA placement: directly below the clavicle and near the right shoulder.
- LA placement: directly below the clavicle and near the left shoulder.
- RL placement: on the right lower abdomen.
- LL placement: on the left lower abdomen.
- V placement: on the chest.



The chest (V) electrode can be placed on one of the following positions:


- V1 placement: on the fourth intercostal space at the right sternal border.
- V2 placement: on the fourth intercostal space at the left sternal border.
- V3 placement: midway between the V2 and V4 electrode positions.
- V4 placement: on the fifth intercostal space at the left midclavicular line.
- V5 placement: on the left anterior axillary line, horizontal with the V4 electrode position.
- V6 placement: on the left midaxillary line, horizontal with the V4 electrode position.
- V3R-V6R placement: on the right side of the chest in positions corresponding to those on the left.
- VE placement: over the xiphoid process.
- V7 placement: on posterior chest at the left posterior axillary line in the fifth intercostal space.
- V7R placement: on posterior chest at the right posterior axillary line in the fifth intercostal space.



12-lead ECG uses 10 electrodes, which are placed on the patient's four limbs and chest. The limb electrodes should be placed on the soft skin and the chest electrodes placed according to the physician's preference.

Lead Placement for Surgical Patients


The surgical site should be taken into consideration when placing electrodes on a surgical patient. e.g. for open-chest surgery, the chest electrodes can be placed on the lateral chest or back. To reduce artifacts and interference from electrosurgical units, you can place the limb electrodes close to the shoulders and lower abdomen and the chest electrodes on the left side of the mid-chest. Do not place the electrodes on the upper arm. Otherwise, the ECG waveform will be very small.





WARNING

- When using electrosurgical units (ESU), patient leads should be placed in a position that is equal distance from the Electrosurgery electrotome and the grounding plate to avoid burns to the patient. Never entangle the ESU cable and the ECG cable together.
- When using electrosurgical units (ESU), never place ECG electrodes near to the grounding plate of the ESU, as this can cause a lot of interference on the ECG signal.



8.3.4 Checking Paced Status

It is important to set the paced status correctly when you start monitoring ECG. The paced symbol is displayed when the [**Paced**] status is set to [**Yes**]. The pace pulse markers "|" are shown on the ECG wave when the patient has a paced signal.

To change the paced status, you can select either:

- the patient information area, or
- [Main Menu]→[Patient Setup]→[Patient Demographics], or,
- the ECG parameter window or waveform area → [Others >>],

and then, select [Paced] from the popup menu and toggle between [Yes] and [No].

If you do not set the paced status, the patient monitor issues a prompt tone when pace pulse is detected. At the same time, the paced symbol flashes and the message "Please confirm the pace of patient" appears in the ECG waveform area. Then, please check and set the paced status of the patient.



WARNING

- For paced patients, you must set [Paced] to [Yes]. If it is incorrectly set to [No], the patient monitor could
 mistake a pace pulse for a QRS and fail to alarm when the ECG signal is too weak. Do not rely entirely on
 rate meter alarms when monitoring patients with pacemakers. Always keep these patients under close
 surveillance.
- For non-paced patients, you must set [Paced] to [No].
- The auto pacer recognization fucntion is not applicable to neonatal patients.

8.4 Understanding the ECG Display

Your display may be configured to look slightly different.



- 1. Lead label of the displayed wave
- 2. ECG gain
- 3. ECG filter label
- 4. Notch filter status

Besides, when a paced signal has been detected, the pace pulse marks "|" are shown on the ECG wave if the [**Paced**] has been set to [**Yes**].



- 1. Current heart rate alarm limits
- 2. Current heart rate
- 3. Heart beat symbol

For 12-lead ECG display screen, refer to the section 12-Lead ECG Monitoring.

8.5 Changing ECG Settings

8.5.1 Accessing ECG Menus

By selecting the ECG parameter window or waveform area, you can access the [**ECG Setup**] menu.

8.5.2 Setting Pacemaker Rate (For Mortara only)

Some pacemaker pulses can be difficult to reject. When this happens, the pulses are counted as a QRS complex and could result in an incorrect HR and failure to detect some arrhythmias. You can set [**Pacemaker Rate**] to the pacemaker's rate in the [**ECG Setup**] menu. In this way, the patient monitor can calculate HR and detect arrhythmias more accurately. When [**Paced**] is set to [**No**], the pacemaker rate cannot be set.

8.5.3 Choosing the Alarm Source

In most cases the HR and PR numerics are identical. In order to avoid simultaneous alarms on HR and PR, the monitor uses either HR or PR as its active alarm source. To change the alarm source, select [**Alm Source**] in the [**ECG Setup**] menu and then select either:

- [**HR**]: if you want the HR to be the alarm source for HR/PR.
- [**PR**]: if you want the PR to be the alarm source for HR/PR.
- [Auto]: If the [Alm Source] is set to [Auto], the patient monitor will use the heart rate from the ECG measurements as the alarm source whenever a valid heart rate is available. If the heart rate becomes unavailable, for example the ECG module is turned off or becomes disconnected, the patient monitor will automatically switch to PR as the alarm source.

8.5.4 Setting the ECG Lead Set

You can set the [**Lead Set**] by selecting [**ECG Setup**]→[**Others>>**]. You can set the [**Lead Set**] as [**Auto**] if the auto lead detection function is available.

8.5.5 Choosing an ECG Display Screen

When monitoring with a 5-lead or 12-lead set, you can select the [**Screens**] Quickkey. In the [**Choose Screen**] window, choose the screen type as:

- [Normal Screen]: The ECG waveform area shows 2 ECG waveforms.
- [ECG 7-Lead Full-Screen]: The whole waveform area shows 7 ECG waveforms only.
- [ECG 7-Lead Half-Screen]: The upper half part of the whole waveform area displays 7 ECG waveforms.

When monitoring with a 12-lead set, you can also choose the screen type as [**ECG 12-Lead Full-Screen**]. When the screen type is set to [**Normal Screen**], cascaded ECG waveforms can be displayed. To cascade ECG waveforms:

- 1. Select the [Screens] Quickkey→[Screen Setup].
- 2. Select [**ECG1 Casc.**] in the second row. A cascaded waveform is displayed in two waveform positions.

8.5.6 Changing the ECG Filter Settings

The ECG filter setting defines how ECG waves are smoothed. To change the filter setting, select [**Filter**] from [**ECG Setup**] and then select the appropriate setting.

- [Monitor]: Use under normal measurement conditions.
- [**Diagnostic**]: Use when diagnostic quality is required. The unfiltered ECG wave is displayed so that changes such as R-wave notching or discrete elevation or depression of the ST segment are visible.
- [Surgery]: Use when the signal is distorted by high frequency or low frequency interference. High frequency interference usually results in large amplitude spikes making the ECG signal look irregular. Low frequency interference usually leads to wandering or rough baseline. In the operating room, the surgery filter reduces artifacts and interference from electrosurgical units. Under normal measurement conditions, selecting [Surgery] may suppress the QRS complexes too much and then interfere with ECG analysis.
- [**ST**]: Use when ST monitoring is applied.



• The [Diagnostic] filter is recommended when monitoring a patient in an environment with slight interference only.

8.5.7 Switching the Notch Filter On or Off

The notch filter removes the line frequency interference. When [**Filter**] is not set to [**Diagnostic**], the notch filter always stays on. When [**Filter**] is set to [**Diagnostic**], you can switch the notch filter on or off as required.

- 1. Select the ECG parameter window or waveform area to enter its setup menu. Then select [Others >>].
- 2. Select [**Notch Filter**] and toggle between [**On**] and [**Off**]. Switching the notch filter on is recommended when there is interference (such as spikes) with the waveform.
- When [Notch Filter] is set on, select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 4. Select [Others >>]→[Notch Freq.] and then select [50Hz] or [60Hz] according to the power line frequency.

8.5.8 Changing the Pacer Reject Settings

Select [**ECG Setup**]→[**Others>>**]→[**Pacer Reject**], and toggle between [**On**] and [**Off**]. When [**Paced**] is set to [**Yes**]:

- When [Pacer Reject] is switched on, the pace pulses are not counted as extra QRS complexes.
- The pace pulse marks "|" are shown on the ECG wave when pace pulses are detected.

When [Paced] is set to [No], the pace markers are not shown on the ECG wave, and the options of [Pacer Reject] are invalid.

8.5.9 About the Defibrillator Synchronization

If a defibrillator is connected, a defibrillator synchronization pulse (100 ms, +5V) is outputted through the Multifunctional Connector every time when the patient monitor detects an R-wave. The defibrillator synchronization function is always enabled.

WARNING

- Improper use of a defibrillator may cause injury to the patient. The user should determine whether to perform defibrillation or not according to the patient's condition.
- Before defibrillation, the user must ensure both defibrillator and monitor has passed the system test and can be safely used jointly.
- Before defibrillation, make sure that the [Filter] is set to [Diagnostic].
- After defibrillation is finished, select the filter mode as required.

8.5.10 Changing ECG Wave Settings

In the [**ECG Setup**] menu:

- If the wave is too small or clipped, you can change its size by selecting an appropriate [**Gain**] setting. If you select [**Auto**] from [**Gain**], the patient monitor will automatically adjust the size of the ECG waves. In normal screen, only the selected ECG wave's size is adjusted. In other screens, all ECG waves' size is adjusted simultaneously.
- You can change the wave sweep speed by selecting [**Sweep**] and then selecting the appropriate setting.

8.5.11 Enabling Smart Lead Off

When the smart lead off function is set on and there is a "lead off" in the lead that has an ECG waveform in filter mode and notch status, if another lead is available, this available lead automatically becomes that lead. The system will re-calculate HR and analyze and detect arrhythmia. When the "lead off" condition is corrected, the leads are automatically switched back.

To switch on/off the smart lead off function, select [Others >>] from the [**ECG Setup**] menu; select [**Smart Lead Off**] and toggle between [**On**] and [**Off**] from the popup menu.

8.5.12 Setting the Alarm Level for ECG Lead Off Alarms

Select [Alarm Setup >>] from the [User Maintenance] menu. You can set [ECGLeadOff Lev.] from the popup menu.

8.5.13 Adjusting QRS Volume

QRS sounds are produced based on the alarm source. To adjust the QRS volume, select [Others >>] from the [ECG Setup] menu; select [QRS Volume] from the popup menu and select the appropriate setting. When valid SpO₂ measured value is available, the system adjusts the pitch tone of QRS sound based on the SpO₂ value.

8.6 About ST Monitoring

- ST segment analysis is not intended for neonatal patients.
- ST segment analysis calculates ST segment elevations and depressions for individual leads and then displays them as numerics in the ST1 and ST2 areas.
- A positive value indicates ST segment elevation; a negative value indicates ST segment depression.
- Measurement unit of the ST segment: mV or mm. You can set the unit in the [Unit Setup] menu from the [User Maintenance] menu.
- Measurement range of the ST segment: -2.0 mV to +2.0 mV.



∕!\ WARNING

The ST algorithm has been tested for accuracy of the ST segment data. The significance of the ST segment changes need to be determined by a clinician.

8.6.1 Switching ST On and Off

To switch ST monitoring on or off:

- 1. In the [ECG Setup] menu, select [ST Analysis >>].
- 2. Select [**ST Analysis**] to toggle between [**On**] and [**Off**].

Reliable ST monitoring can hardly be ensured if:

- You are unable to get a lead that is not noisy.
- Arrhythmias such as atrial fib/flutter cause irregular baseline.
- The patient is continuously ventricularly paced.
- The patient has left bundle branch block.

In these cases, you may consider switching ST monitoring off.

8.6.2 Changing ST Filter Settings

ST-segment analysis can be carried out only when the filter mode is set to [**Diagnostic**] or [**ST**]. When ST-segment analysis is switched on, [**Filter**] will automatically switch to [**ST**] if it is not [**Diagnostic**] or [**ST**]. When ST-segment analysis is switched off, the filter mode automatically switches to previous manual setting.

However, if you switch [**Filter**] to [**Monitor**] or [**Surgery**], ST-segment analysis will turn off automatically. In case that you change [**Monitor**] or [**Surgery**] to [**Diagnostic**] or [**ST**], ST-segment analysis keeps off, you can turn it on manually.

NOTE

- When the filter mode is switched to [Diagnostic], the notch filter automatically switches to [Off]. In this case, you can still set the notch filter to [On] manually.
- When the filter mode is [Monitor], [Surgery] or [ST], the notch filter is fixed to [On], and can not be changed.

8.6.3 Understanding the ST Display

8.6.3.1 ST Numerics

This example shows ST numerics with 5-lead ECG. Your monitor screen may look slightly different from the illustration.

8.6.3.2 ST Segment


ST segment shows a QRS complex segment for each measured ST lead. The current ST segment is drawn in the same color as the ECG wave, usually green, superimposed over the stored reference segment, drawn in a different color. The information is updated once every ten seconds.

To display the ST segment on normal screen:

- 1. Enter the [ST Analysis] menu. Set [ST Analysis] to [On].
- 2. Enter the [Screen Setup] window of [Screens] menu. Set [ST Segment] to be displayed.



Select the ST parameter window or ST segment area and you can enter the [ST Analysis] menu.



8.6.4 Saving the Current ST Segment as Reference

Select [**Save Ref**.] in the [**ST Analysis**] menu to save the current segment as reference. Up to 20 reference segment groups can be saved.

NOTE

 If the memory is full and you do not delete a group before saving a new one, the oldest saved group is deleted automatically.

8.6.5 Changing the Reference Segment

Select the **I** and **I** arrow keys beside the [**Change Ref**.] to switch between different reference segment groups.

8.6.6 Deleting a Reference Segment

To delete the current ST reference segment, select [**Delete Ref.**] in the [**ST Analysis**] menu and then select [**Ok**] in the popup.

8.6.7 Recording the ST Segment

To record the current ST segment and reference segment, select [Record] in the [ST Analysis] menu.

8.6.8 Changing the ST Alarm Limits


High and low ST alarm limits can be set individually for each ECG lead. Alarm limits can also be set separately for single-lead and multi-lead ST monitoring. You can select [**ST Alarm Setup >>**] from [**ST Analysis**] menu and then change ST alarm settings for each lead.

8.6.9 Setting the ST Alarm Delay Time

You can set the ST alarm delay time from the [Others] window of [Alarm Setup] menu.

8.6.10 Adjusting ST Measurement Points

As shown in the figure below, the ST measured for each beat complex is the vertical difference between two measurement points with the R-wave peak as the baseline for the measurement.



The ISO and ST points need to be adjusted when you start monitoring and if the patient's heart rate or ECG morphology changes significantly. Exceptional QRS complexes are not considered for ST-segment analysis.



WARNING

Always make sure that the positions of ST measurement points are appropriate for your patient.

To adjust the ST measurement points:

- 1. In the [ST Analysis] menu, select [Adjust ST Point >>]. In the [Adjust ST Point] window, three vertical lines represent the ISO, J and ST point positions respectively.
- 2. Select [View Leads] and use the Knob to select an ECG lead with obvious J point and R wave.
- 3. Select [ISO], [J] or [ST Point] and then use the Knob to adjust the position of each point.
 - ◆ The ISO-point (isoelectric) position is given relative to the R-wave peak. Position the ISO-point in the middle of the flattest part of the baseline (between the P and Q waves of in front of the P wave).
 - ◆ The J-point position is given relative to the R-wave peak and helps locating the ST-point. Position the J-point at the end of the QRS complex and the beginning of the ST segment.
 - ◆ The ST-point is positioned a fixed distance from the J-point. Move the J-point to position the ST-point at the midpoint of the ST segment. Position the ST-point relative to the J-point at either [J+60/80ms], [J+40ms], [J+60ms] or [J+80ms]. When [J+60/80ms] is selected, the ST-point will be positioned 80 ms (heart rate 120 bpm or less) or 60 ms (heart rate more than 120 bpm) from the J-point.

8.7 About Arrhythmia Monitoring

Arrhythmia analysis provides information about your patient's condition, including heart rate, PVC rate, rhythm and ectopics.



riangle warning

- Arrhythmia analysis program is intended to detect ventricular arrhythmias. It is not designed to detect atrial or supraventricular arrhythmias. It may incorrectly identify the presence or absence of an arrhythmia. Therefore, a physician must analyze the arrhythmia information with other clinical findings.
- The arrhythmia detection is not intended for neonatal patients.

8.7.1 Understanding the Arrhythmia Events

Mindray algorithm

Arrhythmia message	Description	Category		
Asystole	No QRS detected within the set time threshold in absence of ventricular fibrillation			
	or chaotic signal.			
Vfib/Vtac	A fibrillatory wave for 6 consecutive seconds.			
VIID/VLaC	A dominant rhythm of adjacent Vs and a HR > the V-Tac HR limit.			
Vtac	The consecutive PVCs > Vtac PVCs limit, and the HR > the Vtac HR limit.	Lethal arrhythmia		
Vent. Brady	The consecutive PVCs ≥ the Vbrd threshold and the ventricular HR < the Vbrd Rate			
vent. brady	threshold.			
Extreme Tachy	The heart rate is greater than the extreme tachycardia limit.			
Extreme Brady	The heart rate is less than the extreme bradycardia limit.			
PVCs	PVCs/min exceeds high limit	Nonlethal		
Da a a a Nat Da aire a	No pace pulse detected for 1.75 x average R-to-R intervals following a QRS complex	arrhythmia		
Pacer Not Pacing	(for paced patients only).			
De seu Net Continue	No QRS complex detected for 300 milliseconds following a pace pulse (for paced			
Pacer Not Capture	patients only).			
PVC	One PVC detected in normal heartbeats.			
Couplet	Paired PVCs detected in normal heartbeats.			
VT > 2	More than 2 consecutive PVCs within the last minute.			
Bigeminy	A dominant rhythm of N, V, N, V, N, V.			
Trigeminy	A dominant rhythm of N, N, V,N, N, V, N, N, V.			
R on T	R on T detected in normal heartbeats.			
	No beat detected for 1.75 x average R-R interval for HR <120, or			
Missed Beats	No beat for 1 second with HR > 120 (for non-paced patients only), or			
	No beat detected for more than the set pause threshold.			
Brady	The average heart rate is less than the bradycardia limit.			
Tachy	The average heart rate is greater than the tachycardia limit.			
** - **	The consecutive PVCs > the Vbrd PVCs limit, and the HR is between			
Vent. Rhythm	Vbrd Rate limit and the Vtac Rate limit.			

Arrhythmia message	Description	Category
Multif. PVC	Multiform PVCs detected in Multif. PVC's Window (which is adjustable).	
Nonsus. Vtac	The consecutive PVCs < the Vtac PVCs limit but > 2, and HR > the Vtac Rate limit.	
Pause	No QRS detected within the set time threshold of pause.	
Irr. Rhythm	Consistently irregular rhythm.	

Mortara algorithm

Arrhythmia Message	Description	Category
Asystole	No QRS complex detected within the set time threshold (in absence of	Lethal
	ventricular fibrillation or chaotic signals).	arrhythmia
Vfib	Ventricular fibrillation occurs and persists for 6 seconds.	
Vtac	Ventricular HR is greater or equal to the preset threshold and the number of	
	consecutive PVCs is greater than the preset threshold.	
PVCs	PVCs/min exceeds high limit	Nonlethal
PNP	No pace pulse detected for (60*1000/pace rate +90) milliseconds following a QRS	arrhythmia
	complex or a pacer pulse (for paced patients only).	
PNC	No QRS complex detected for 300 milliseconds following a pace pulse (for paced	
	patients only).	
Multif. PVC	More than 2 PVCs of different forms occur in the predefined search window	
	(3-31).	
Couplet	Paired PVCs are detected.	
VT > 2	Ventricular HR is greater than or equal to the preset threshold and the number of	
	PVCs is greater than or equal to 3 but less than the preset threshold.	
Vent. Rhythm	Ventricular HR is less than the preset threshold and the number of PVCs is greater	
	than or equal to 3.	
Bigeminy	A dominant rhythm of N, V,N, V, N, V.	
Trigeminy	A dominant rhythm of N, N, V,N, N, V, N, N, V.	
R on T	R on T is detected.	
Irr. Rhythm	Consistently irregular rhythm	
Missed Beats	No beat detected for 1.75x average R-R interval for HR <120, or	
	No beat for 1 second with HR >120 (for non-paced patients only), or	
	No beat detected for more than the set pause threshold.	
Brady	The HR is less than the set bradycardia low limit.	
Tachy	The HR is greater than the set tachycardia high limit.	

8.7.2 Changing Arrhythmia Alarm Settings

To change arrhythmia alarm settings, select the ECG parameter area or waveform area \rightarrow [ECG Setup] \rightarrow [Arrh. Analysis >>]. In the pop-up menu, you can set the [Alm Lev] to [High], [Med], [Low] or [Message], or switch on lethal arrhythmia analysis alarms only or switch on/off all arrhythmia analysis alarms. In the [Alarm Setup] menu from the [User Maintenance] menu, you can enable/disable turning off lethal arrhythmia analysis alarms.



WARNING

If you switch off all arrhythmia analysis alarms, the monitor cannot give any arrhythmia analysis alarm. Always keep the patient under close surveillance.

8.7.3 Changing Arrhythmia Threshold Settings

Select the ECG parameter window or waveform area → [Arrh. Analysis >>] → [Arrh. Threshold], and you can then change threshold settings for some arrhythmia alarms. In case an arrhythmia violates its threshold, an alarm will be triggered. The asystole delay time relates to ECG relearning. When HR is less than 30 bpm, it is recommended to set the asystole delay time to 10 seconds.

Mindray algorithm

Arrh. event	Range	Default	Step	Unit
PVCs High	1 to 10	10	1	/min
Asys. Delay	3 to 10	5	1	S
Ta about that	50.4.200	Adult: 120	г	h m m
Tachy High	60 to 300	Pediatric: 160	5	bpm
Product our	15 to 120	Adult: 50	5	bpm
Brady Low	15 to 120	Pediatric: 75	3	
Future and Talahu	120 to 300	Adult: 160		bpm
Extreme Tachy		Pediatric: 180	5	
Extreme Brady	15 to 60	Adult: 35	5	bpm
		Pediatric: 50		
Multif. PVC's Window	3 to 31	15	1	/min
Vtac Rate	100 to 200	130	5	bpm
Vtac PVCs	3 to 99	6	1	/min
Pause Time	1.5, 2.0,2.5	2	/	S
Vbrd PVCs	3 to 99	5	1	/min
Vbrd Rate	15 to 60	40	5	bpm

Mortara algorithm

Arrh. event	Range	Default	Step	Unit
PVCs High	1 to 10	10	1	/min
Asys. Delay	2 to 10	5	1	S
Vtac Rate	100 to 200	130	5	bpm
Vtac PVC	3 to 12	6	1	beats
Multif. PVC	3 to 31	15	1	beats
Tasky Hisk	Adult: 100 to 300	Adult: 100	5	bpm
Tachy High	Pediatric: 160 to 300	Pediatric: 160	٠	
Brady Low	Adult: 15 to 60	Adult: 60	5	bpm
brady Low	Pediatric: 15 to 80	Pediatric: 80	3	

8.7.4 Setting the Extended Arrhythmia (For Mindray Algorithm only)

The following arrhythmia events are defined as extended arrhythmia:

- Extreme Tachy
- Extreme Brady
- Vent. Brady
- Nonsus. Vtac
- Multif. PVC
- Irr. Rhythm
- Pause

You can select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password→select [Alarm Setup >>], and set [Extended Arrh.] to [Enable] or [Disable]. When [Extended Arrh.] is set to [Disable], the patient monitor does not analysis the extended arrhythmia events and corresponding alarms are not given.

ACAUTION

Set [Extended Arrh.] to [Disable] when the patient monitor is connected to the Central Monitoring
System of version prior to 06.01.00. Failure to do so may cause the Central Monitoring System unable to
display extended arrhythmia related alarms normally when extended arrhythmia occurs.

8.7.5 Reviewing Arrhythmia Events

Please refer to the *Review* chapter.

8.8 ECG Relearning

8.8.1 Initiating an ECG Relearning Manually

During ECG monitoring, you may need to initiate an ECG relearning when the patient's ECG template changes dramatically. A change in the ECG template could result in:

- incorrect arrhythmia alarms
- loss of ST measurement, and/or
- inaccurate heart rate

ECG relearning allows the monitor to learn the new ECG template so as to correct arrhythmia alarms and HR value, and restore ST measurements. To initiate relearning manually, select the ECG parameter window or waveform area → [Relearn]. When the patient monitor is learning, the message [ECG Learning] is displayed in the technical alarm area.

ACAUTION

Take care to initiate ECG relearning only during periods of normal rhythm and when the ECG signal is
relatively noise-free. If ECG learning takes place during ventricular rhythm, the ectopics may be incorrectly
learned as the normal QRS complex. This may result in missed detection of subsequent events of V-Tach
and V-Fib.


8.8.2 Automatic ECG Relearning

ECG relearning is initiated automatically whenever:

- The ECG lead or lead label is changed
- The ECG lead is re-connected
- A new patient is admitted
- After the calibration is completed, select [Stop Calibrating ECG]
- A switch happens between the options of screen type during 5/12-lead ECG monitoring.
- The paced status of the patient is changed.

8.9 12-Lead ECG Monitoring (for iPM 12 and iPM 10 patient monitors only)

- 1. Refer to the section **8.3.3 ECG Lead Placements** for placing the electrodes.
- In the [ECG Setup] menu, select [Lead Set]→[12-Lead]. Select [Screens] Quickkey→[ECG 12-Lead Full-Screen].



There are totally 12 ECG waves and 1 rhythm wave displayed on the screen. The rhythm lead is ECG I before entering the 12-lead ECG monitoring screen. The ST numerics are displayed in three groups:

■ ST Ant (anterior): V1, V2, V3, V4

■ ST Inf (inferior): II, III, aVF, (aVR)

■ ST Lat (lateral): I, aVL, V5, V6

Although aVR is displayed in the ST Inf group, it is not an inferior lead.

Additionally, the 12-lead ECG monitoring has the following features:

- The [Filter] mode is fixed to [Diagnostic] and cannot be changed.
- In the adult mode, the 🖼 hardkey on the monitor's front is disabled

FOR YOUR NOTES

9

Monitoring Respiration (Resp)

9.1 Introduction


Impedance respiration is measured across the thorax. When the patient is breathing or ventilated, the volume of air changes in the lungs, resulting in impedance changes between the electrodes. Respiration rate (RR) is calculated from these impedance changes, and a respiration waveform appears on the patient monitor screen.

9.2 Safety Information



- When monitoring the patient's respiration, do not use ESU-proof ECG cables.
- If you do not set the detection level for the respiration correctly in manual detection mode, it may not be
 possible for the monitor to detect apnea. If you set the detection level too low, the monitor is more likely
 to detect cardiac activity, and to falsely interpret cardiac activity as respiratory activity in the case of
 apnea.
- The respiration measurement does not recognize the cause of apneas. It only indicates an alarm if no
 breath is detected when a pre-adjusted time has elapsed since the last detected breath. Therefore, it
 cannot be used for diagnostic purpose.
- If operating under conditions according to the EMC Standard EN 60601-1-2 (Radiated Immunity 3V/m), field strengths above 1V/m may cause erroneous measurements at various frequencies. Therefore it is recommended to avoid the use of electrically radiating equipment in close proximity to the respiration measurement unit.

9.3 Understanding the Resp Display



By selecting the waveform area or parameter area, you can enter the [**Resp Waveform**] menu. By selecting the Resp parameter window, you can enter the [**Resp Setup**] menu.

NOTE

• Respiration monitoring is not for use on the patients who are very active, as this will cause false alarms.


9.4 Placing Resp Electrodes

As the skin is a poor conductor of electricity, preparing the skin is necessary for a good Respiration signal. You can refer to the ECG section for how to prepare the skin.

As the Respiration measurement adopts the standard ECG electrode placement, you can use different ECG cables (3-lead, 5-lead or 12-lead). Since the respiration signal is measured between two ECG electrodes, if a standard ECG electrode placement is applied, the two electrodes should be RA and LA of ECG Lead I, or RA and LL of ECG Lead II.

NOTE

To optimize the respiration waveform, place the RA and LA electrodes horizontally when monitoring respiration with ECG Lead I; place the RA and LL electrodes diagonally when monitoring respiration with ECG Lead II.



9.4.1 Optimizing Lead Placement for Resp

If you want to measure Resp and you are already measuring ECG, you may need to optimize the placement of the two electrodes between which Resp will be measured. Repositioning ECG electrodes from standard positions results in changes in the ECG waveform and may influence ST and arrhythmia interpretation.

9.4.2 Cardiac Overlay

Cardiac activity that affects the Resp waveform is called cardiac overlay. It happens when the Resp electrodes pick up impedance changes caused by the rhythmic blood flow. Correct electrodes placement can help to reduce cardiac overlay: avoid the liver area and the ventricles of the heart in the line between the respiratory electrodes. This is particularly important for neonates.

9.4.3 Abdominal Breathing

Some patients with restricted movement breathe mainly abdominally. In these cases, you may need to place the left leg electrode on the left abdomen at the point of maximum abdominal expansion to optimise the respiratory wave.

9.4.4 Lateral Chest Expansion

In clinical applications, some patients (especially neonates) expand their chests laterally, causing a negative intrathoracic pressure. In these cases, it is better to place the two respiration electrodes in the right midaxillary and the left lateral chest areas at the patient's maximum point of the breathing movement to optimise the respiratory waveform.

9.5 Choosing the Respiration Lead

In the [Resp Setup] menu, select [Resp Lead] and toggle between [I] and [II].

9.6 Changing the Apnea Alarm Delay

The apnea alarm is a high-level alarm used to detect apneas. You can set the apnea alarm delay time after which the patient monitor alarms if the patient stops breathing. In the [**Resp Setup**] menu, select [**Apnea Delay**] and then select the appropriate setting. The [**Apnea Delay**] of Resp, CO₂, and AG module keeps consistent with each other.

9.7 Changing Resp Detection Mode

In the [Resp Setup] menu, select [Detection Mode] and toggle between [Auto] and [Manual].

In auto detection mode, the patient monitor adjusts the detection level automatically, depending on the wave height and the presence of cardiac artifact. Note that in auto detection mode, the detection level (a dotted line) is not displayed on the waveform.

Use auto detection mode for situations where:

- ◆ The respiration rate is not close to the heart rate.
- Breathing is spontaneous, with or without continuous positive airway pressure (CPAP).
- Patients are ventilated, except patients with intermittent mandatory ventilation (IMV).
- In manual detection mode, you adjust the dotted detection level line to the desired level by selecting [**Upper Line**] or [**Lower Line**] and then selecting or beside them. Once set, the detection level will not adapt automatically to different respiration depths. It is important to remember that if the depth of breathing changes, you may need to change the detection level.

Use manual detection mode for situations where:

- The respiration rate and the heart rate are close.
- Patients have intermittent mandatory ventilation.
- Respiration is weak. Try repositioning the electrodes to improve the signal.

In Auto Detection Mode, if you are monitoring Resp and ECG is switched off, the monitor cannot compare the ECG and Resp rates to detect cardiac overlay. The respiration detection level is automatically set higher to prevent the detection of cardiac overlay as respiration.

In Manual Detection Mode, cardiac overlay can in certain situations trigger the respiration counter. This may lead to a false indication of a high respiration or an undetected apnea condition. If you suspect that cardiac overlay is being registered as breathing activity, raise the detection level above the zone of cardiac overlay. If the Resp wave is so small that raising the detection level is not possible, you may need to optimize the electrode placement as described in the section "Lateral Chest Expansion".

9.8 Changing Resp Wave Settings



• When monitoring in manual detection mode, make sure to check the respiration detection level after you have increased or decreased the size of the respiration wave.

In the [Resp Setup] menu, you can:

- Select [Gain] and then select an appropriate setting. The bigger the gain is, the larger the wave amplitude is.
- Select [Sweep] and then select an appropriate setting. The faster the wave sweeps, the wider the wave is.

9.9 Setting RR Source

To set RR source:

- 1. Enter the [**Resp Setup**] menu.
- 2. Select [RR Source] and then select a source or [Auto] from the dropdown list.

The dropdown list displays the currently available RR source. When you select [**Auto**], the system will automatically select the RR source according to the priority. When the current RR source does not have valid measurement, the system will automatically switch the [**RR Source**] to [**Auto**]. RR source switches back to impedance respiration if you press the silence hardkey on the monitor's front during an apnea alarm.

The priority of RR source is (from high to low): CO₂ measurement, and impedance respiration measurement. The [**RR Source**] settings of Resp, CO₂, and AG module are linked.

The RR source options and description are shown in the table below.

Option	Description
Auto	RR source is automatically selected according to the priority.
CO ₂	RR source is from CO ₂ measurement.
ECG	RR source is from impedance respiration measurement.

9.10 Setting alarm properties

Select [**Alarm Setup >>**] from the [**Resp Setup**] menu. In the popup menu, you can set alarm properties for this parameter.

9.11 Switching Resp Measurement On/Off

To switch Resp measurement on, select [Imped. Resp Measure. ON] from the [Resp Setup] menu. To switch Resp measurement off, select [Imped. Resp Measure. OFF] from the [Resp Setup] menu and then select [Yes] from the popup dialog box. Then, a line is displayed in the waveform area and no numeric but [Measurement OFF] message is displayed in the parameter area.

FOR YOUR NOTES		

10 Monitoring PR

10.1 Introduction

The pulse numeric counts the arterial pulsations that result from the mechanical activity of the heart. You can display a pulse from any measured SpO_2 or any arterial pressure (see the IBP section). The displayed pulse numeric is color-coded to match its source.



- 1. PR: detected beats per minute.
- 2. PR Source

10.2 Setting the PR Source

The current pulse source is displayed in the PR parameter area. The pulse rate chosen as pulse source:

- is monitored as system pulse and generates alarms when you select PR as the active alarm source;
- is stored in the monitor's database and reviewed in the graphic/tabular trends; in trend graphs, as the PR curve is in the same color with the PR source, it is unlikely to distinguish the PR source;
- is sent via the network to the central monitoring system, if available.

To set which pulse rate as PR source:

- 1. Enter the [**SpO**₂ **Setup**] menu.
- 2. Select [**PR Source**] and then select a label or [**Auto**] from the popup menu.

The popup menu displays the currently available PR sources from top to bottom by priority. When you select [**Auto**], the system will automatically select the first option as the PR source from the popup menu. When the current PR source is unavailable, the system will automatically switch [**PR Source**] to [**Auto**]. When you select [**IBP**], the system will automatically select the first pressure label as the PR source from the popup menu.

10.3 Selecting the Active Alarm Source

In most cases the HR and pulse numerics are identical. In order to avoid simultaneous alarms on HR and Pulse, the monitor uses either HR or Pulse as its active alarm source. To change the alarm source, select [Alm Source] in the [ECG Setup] or [SpO₂ Setup] menu and then select either:

- [**HR**]: The monitor will use the HR as the alarm source for HR/pulse.
- [**PR**]: The monitor will use the PR as the alarm source for HR/pulse.
- [Auto]: If the [Alm Source] is set to [Auto], the monitor will use the heart rate from the ECG measurement as the alarm source whenever the ECG measurement is switched on and a valid heart rate is available. If the heart rate becomes unavailable, for example if leads becomes disconnected, and a pulse source is switch on and available, the monitor will automatically switch to Pulse as the alarm source. When the Leads Off condition is corrected, the monitor will automatically switch back to the heart rate as the alarm source.


10.4 QRS Tone

When PR is used as the alarm source, the PR source will be used as a source for the QRS tone. You can change the QRS volume by adjusting [**Beat Vol**] in the [**SpO**₂ **Setup**] menu. When a valid SpO_2 value exists, the system will adjust the pitch tone of QRS volume according to the SpO_2 value.

11.1 Introduction

 SpO_2 monitoring is a non-invasive technique used to measure the amount of oxygenated haemoglobin and pulse rate by measuring the absorption of selected wavelengths of light. The light generated in the probe passes through the tissue and is converted into electrical signals by the photodetector in the probe. The SpO_2 module processes the electrical signal and displays a waveform and digital values for SpO_2 and pulse rate.

This device is calibrated to display functional oxygen saturation. It provides four measurements:



- 1. Pleth waveform (Pleth): visual indication of patient's pulse. The waveform is not normalized.
- 2. Oxygen saturation of arterial blood (SpO₂): percentage of oxygenated hemoglobin in relation to the sum of oxyhemoglobin and deoxyhemoglobin.
- 3. Perfusion index (PI): gives the numerical value for the pulsatile portion of the measured signal caused by arterial pulsation. PI is an indicator of the pulsatile strength. You can also use it to assess the quality of SpO₂ measurement. Above 1 is optimal, between 0.3 and 1 is acceptable. Below 0.3 indicates low perfusion; reposition the SpO₂ sensor or find a better site. If low perfusion persists, choose another method to measure oxygen saturation if possible. PI is available for Mindray SpO₂ module or Masimo SpO₂ module.
- 4. Perfusion indicator: the pulsatile portion of the measured signal caused by arterial pulsation.
- 5. Pulse rate (derived from pleth wave): detected pulsations per minute.

11.2 Safety



WARNING

- Use only SpO₂ sensors specified in this manual. Follow the SpO₂ sensor's instructions for use and adhere to all warnings and cautions.
- When a trend toward patient deoxygenation is indicated, blood samples should be analyzed by a laboratory co-oximeter to completely understand the patient's condition.
- Do not use SpO₂ sensors during magnetic resonance imaging (MRI). Induced current could potentially cause burns. The sensor may affect the MRI image, and the MRI unit may affect the accuracy of the oximetry measurements.
- Prolonged continuous monitoring may increase the risk of undesirable changes in skin characteristics, such as irritation, reddening, blistering or burns. Inspect the sensor site every two hours and move the sensor if the skin quality changes. Change the application site every four hours. For neonates, or patients with poor peripheral blood circulation or sensitive skin, inspect the sensor site more frequently.

11.3 Identifying SpO₂ Modules

To identify which SpO_2 module is incorporated into your patient monitor, see the color of the SpO_2 connector and the company logo located at the patient monitor. The color of the cable connector matches the company as shown below:

- Mindray SpO₂ module: a blue connector without logo.
- Masimo SpO₂ module: a purple connector with a logo of Masimo SET.
- Nellcor SpO₂ module: a grey connector with a logo of Nellcor.

The connectors for these three SpO₂ sensors are mutually exclusive.

11.4 Applying the Sensor

- 1. Select an appropriate sensor according to the module type, patient category and weight.
- 2. Remove colored nail polish from the application site.
- 3. Apply the sensor to the patient.
- 4. Select an appropriate adapter cable according to the connector type and plug this cable into the SpO_2 connector.
- 5. Connect the sensor cable to the adapter cable.

11.5 Changing SpO₂ Settings

11.5.1 Accessing SpO₂ Menus

By selecting the SpO₂ parameter window or waveform area, you can access the [**SpO**₂ **Setup**] menu.

11.5.2 Adjusting the Desat Alarm

The desat alarm is a high level alarm notifying you of potentially life threatening drops in oxygen saturation. Select [Alarm Setup >>] from the [SpO₂ Setup] menu. From the popup menu, you can set low alarm limit, alarm switch, and alarm recording for [Desat]. When the SpO₂ value is below the desat alarm limit and desat alarm switch is set on, the message [SpO₂ Desat] is displayed.

11.5.3 Setting SpO₂ Sensitivity

For Masimo SpO₂ module, you can set [**Sensitivity**] to [**Normal**] or [**Maximum**] in the [**SpO₂ Setup**] menu. When the [**Sensitivity**] is set to [**Maximum**], the patient monitor is more sensitive to minor signals. When monitoring critically ill patients whose pulsations are very weak, it is strongly recommended that the sensitivity is set to [**Maximum**]. When monitoring neonatal or non-critically ill patients who tend to move a lot, noise or invalid signals may be caused. In this case, it is recommended that the sensitivity is set to [**Normal**] so that the interference caused by motion can be filtered and therefore the measurement stability can be ensured.

11.5.4 Changing Averaging Time

The SpO_2 value displayed on the monitor screen is the average of data collected within a specific time. The shorter the averaging time is, the quicker the patient monitor responds to changes in the patient's oxygen saturation level. Contrarily, the longer the averaging time is, the slower the patient monitor responds to changes in the patient's oxygen saturation level, but the measurement accuracy will be improved. For critically ill patients, selecting shorter averaging time will help understanding the patient's state.

To set the averaging time:

- For Mindray SpO₂ module, select [**Sensitivity**] in the [**SpO**₂ **Setup**] menu and then toggle between [**High**], [**Med**] and [**Low**], which respectively correspond to 7 s, 9 s and 11 s.
- For Masimo SpO₂ module, select [**Averaging**] in the [**SpO**₂ **Setup**] menu and then toggle between [**2-4 s**], [**4-6 s**], [**8 s**], [**10 s**], [**12 s**], [**14 s**] and [**16 s**].

11.5.5 Monitoring SpO₂ and NIBP Simultaneously

When monitoring SpO₂ and NIBP on the same limb simultaneously, you can switch [**NIBP Simul**] on in the [**SpO**₂ **Setup**] menu to lock the SpO₂ alarm status until the NIBP measurement ends. If you switch [**NIBP Simul**] off, low perfusion caused by NIBP measurement may lead to inaccurate SpO₂ readings and therefore cause false physiological alarms.

11.5.6 Sat-Seconds Alarm Management

With traditional alarm management, high and low alarm limits are set for monitoring oxygen saturation. During monitoring, as soon as an alarm limit is violated, an audible alarm immediately sounds. When the patient % SpO₂ fluctuates near an alarm limit, the alarm sounds each time the limit is violated. Such frequent alarm can be distracting. Nellcor's Sat-Seconds alarm management technique is used to reduce these nuisance alarms.

The Sat-Seconds feature is available with the Nellcor SpO₂ module to decrease the likelihood of false alarms caused by motion artifacts. To set the Sat-Seconds limit, select [**Sat-Seconds**] in the [**SpO**₂ **Setup**] menu and then select the appropriate setting.


With Sat-Seconds alarm management, high and low alarm limits are set in the same way as traditional alarm management. A Sat-Seconds limit is also set. The Sat-Seconds limit controls the amount of time that SpO_2 saturation may be outside the set limits before an alarm sounds. The method of calculation is as follows: the number of percentage points that the SpO_2 saturation falls outside the alarm limit is multiplied by the number of seconds that it remains outside the limit. This can be stated as the equation:

Sat-Seconds = Points × Seconds

Only when the Sat-Seconds limit is reached, the monitor gives a Sat-Seconds alarm. For example, the figure below demonstrates the alarm response time with a Sat-Seconds limit set at 50 and a low SpO_2 limit set at 90%. In this example, the patient % SpO_2 drops to 88% (2 points) and remains there for 2 seconds. Then it drops to 86% (4 points) for 3 seconds, and then to 84% (6 points) for 6 seconds. The resulting Sat-Seconds are:

% SpO₂	Seconds	Sat-Seconds	
2×	2=	4	
4×	3=	12	
6×	6=	36	
Total Sat-Seconds=		52	

After approximately 10.9 seconds, a Sat-Second alarm would sound, because the limit of 50 Sat-Seconds would have been exceeded.



Saturation levels may fluctuate rather than remaining steady for a period of several seconds. Often, the patient % SpO₂ may fluctuate above and below an alarm limit, re-entering the non-alarm range several times. During such fluctuation, the monitor integrates the number of %SpO₂ points, both positive and negative, until either the Sat-Seconds limit is reached, or the patient %SpO₂ re-enters the non-alarm range and remains there.

11.5.7 Changing the Speed of the Pleth Wave

In the [\mathbf{SpO}_2 \mathbf{Setup}] menu, select [\mathbf{Sweep}] and then select the appropriate setting. The faster the waveform sweeps, the wider the waveform is.

11.5.8 Setting the Alarm Level for SpO₂ Sensor Off Alarm

Select [**Alarm Setup >>**] from the [**User Maintenance**] menu. You can set the [**SpO2SensorOff Lev.**] in the popup menu.

11.5.9 Setting the SpO₂ Tone Mode

Select [Others >>] from the [User Maintenance] menu. In the popup menu, you can set [SpO₂ Tone] as [Mode 1] or [Mode 2].



• The same SpO₂ tone mode shall be used for the same patient monitors in a single area.

11.6 Measurement Limitations

If you doubt the measured SpO_2 , check patient vital signs first. Then check the patient monitor and SpO_2 sensor. The following factors may influence the accuracy of measurement:

- Ambient light
- Physical movement (patient and imposed motion)
- Diagnostic testing
- Low perfusion
- Electromagnetic interference, such as MRI environment
- Electrosurgical units
- Dysfunctional haemoglobin, such as carboxyhemoglobin (COHb)and methemoglobin (MetHb)
- Presence of certain dyes, such as methylene and indigo carmine
- Inappropriate positioning of the SpO₂ sensor, or use of incorrect SpO₂
- Drop of arterial blood flow to immeaurable level caused by shock, anemia, low temperature or vasoconstrictor.

11.7 Masimo Information



Masimo Patents

This device is covered under one or more the following U.S.A. patents: 5,758,644, 6,011,986, 6,699,194, 7,215,986, 7,254,433, 7,530,955 and other applicable patents listed at: www.masimo.com/patents.htm.

■ No Implied License

Possession or purchase of this device does not convey any express or implied license to use the device with unauthorized sensors or cables which would, alone, or in combination with this device, fall within the scope of one or more of the patents relating to this device.

11.8 Nellcor Information



■ Nellcor Patents

This device may be covered by one or more of the following US patents and foreign equivalents: 5,485,847, 5,676,141, 5,743,263, 6,035,223, 6,226,539, 6,411,833, 6,463,310, 6,591,123, 6,708,049, 7,016,715, 7,039,538, 7,120,479, 7,120,480, 7,142,142, 7,162,288, 7,190,985, 7,194,293, 7,209,774, 7,212,847, 7,400,919.

No Implied License

Possession or purchase of this device does not convey any express or implied license to use the device with unauthorized replacement parts which would, alone, or in combination with this device, fall within the scope of one or more of the patents relating to this device.

12 Monitoring NIBP

12.1 Introduction

The patient monitor uses the oscillometric method for measuring the non-invasive blood pressure (NIBP). This measurement can be used for adults, pediatrics and neonates.

Automatic non-invasive blood pressure monitoring uses the oscillometric method of measurement. To understand how this method works, we'll compare it to the auscultative method. With auscultation, the clinician listens to the blood pressure and determines the systolic and diastolic pressures. The mean pressure can then be calculated with reference to these pressures as long as the arterial pressure curve is normal.

Since the monitor cannot hear the blood pressure, it measures cuff pressure oscillation amplitudes. Oscillations are caused by blood pressure pulses against the cuff. The oscillation with the greatest amplitude is the mean pressure. This is the most accurate parameter measured by the oscillometric method. Once the mean pressure is determined, the systolic and diastolic pressures are calculated with reference to the mean.

Simply stated, auscultation measures systolic and diastolic pressures and the mean pressure is calculated. The oscillometric method measures the mean pressure and determines the systolic and diastolic pressures.

As specified by IEC 60601-2-30/EN60601-2-30, NIBP measurement can be performed during electro-surgery and discharge of defibrillator.

NIBP diagnostic significance must be decided by the doctor who performs the measurement.

NOTE

 Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method or an intra-arterial blood pressure measurement device, within the limits prescribed by the American National Standard, Manual, electronic, or automated sphygmomanometers.

12.2 Safety



WARNING

- Be sure to select the correct patient category setting for your patient before measurement. Do not apply
 the higher adult settings for pediatric or neonatal patients. Otherwise it may present a safety hazard.
- Do not measure NIBP on patients with sickle-cell disease or any condition where skin damage has occurred
 or is expected.
- Use clinical judgement to determine whether to perform frequent unattended blood pressure measurements on patients with severe blood clotting disorders because of the risk of hematoma in the limb fitted with the cuff.
- Do not use the NIBP cuff on a limb with an intravenous infusion or arterial catheter in place. This could cause tissue damage around the catheter when the infusion is slowed or blocked during cuff inflation.
- If you doubt the NIBP readings, determines the patient's vital signs by alternative means and then verify that the monitor is working correctly.

12.3 Measurement Limitations

Measurements are impossible with heart rate extremes of less than 40bpm or greater than 240bpm, or if the patient is on a heart-lung machine.

The measurement may be inaccurate or impossible:

- If a regular arterial pressure pulse is hard to detect
- With excessive and continuous patient movement such as shivering or convulsions
- With cardiac arrhythmias
- Rapid blood pressure changes
- Severe shock or hypothermia that reduces blood flow to the peripheries
- Obesity, where a thick layer of fat surrounding a limb dampens the oscillations coming from the artery

12.4 Measurement Methods

There are three methods of measuring NIBP:

- Manual: measurement on demand.
- Auto: continually repeated measurements at set intervals.
- STAT: continually rapid series of measurements over a five minute period, then return to the previous mode.

12.5 Setting Up the NIBP Measurement

12.5.1 Preparing to Measure NIBP

- 1. Power on the monitor.
- 2. Verify that the patient category is correct. Change it if necessary.
- 3. Plug the air tubing into the NIBP connector on the patient monitor.
- 4. Select a correct sized cuff and then apply it as follows:
 - Determine the patient's limb circumference.
 - Select an appropriate cuff by referring to the limb circumference marked on the cuff. The width of the cuff should be 40% (50% for neonates) of the limb circumference, or 2/3 of the upper arm's length. The inflatable part of the cuff should be long enough to encircle at least 50% to 80% of the limb.
 - lackApply the cuff to an upper arm or leg of the patient and make sure the Φ marking on the cuff matches the artery location. Do not wrap the cuff too tightly around the limb. It may cause discoloration, and ischemia of the extremities. Make sure that the cuff edge falls within the marked range. If it does not, use a larger or smaller cuff that will fit better.
- 5. Connect the cuff to the air tubing and make sure that the bladder inside the cover is not folded and twisted.

NOTE

• The use of the equipment is restricted to one patient at a time.

12.5.2 Starting and Stopping Measurements

Select the [**NIBP Measure**] QuickKey and you can start the desired measurement from the popup menu. You can select [**Stop All**] QuickKey to stop all NIBP measurements. You can also start and stop measurements by using the hardkey on the monitor's front panel.

12.5.3 Correcting the Measurement if Limb is not at Heart Level

The cuff should be applied to a limb at the same level as the patient's heart. If the limb is not at the heart level, to the displayed value:

- Add 0.75 mmHg (0.10 kPa) for each centimetre higher, or
- Deduct 0.75 mmHg (0.10 kPa) for each centimeter lower.

12.5.4 Enabling NIBP Auto Cycling and Setting the Interval

- 1. Select the NIBP parameter window to enter the [**NIBP Setup**] menu.
- 2. Select [**Interval**] and then select a desired time interval. Selecting [**Manual**] switches to manual mode.
- 3. Start a measurement manually. The monitor will then automatically repeat NIBP measurements at the set time interval.

12.5.5 Starting a STAT Measurement


- 1. Select the NIBP parameter window to enter the [**NIBP Setup**] menu.
- 2. Select [NIBP STAT]. The STAT mode initiates 5 minutes of continuous, sequential, automatic NIBP measurements.



Continuous non-invasive blood pressure measurements may cause purpura, ischemia and neuropathy in the limb with the cuff. Inspect the application site regularly to ensure skin quality and inspect the extremity of the cuffed limb for normal color, warmth and sensitivity. If any abnormity occurs, move the cuff to another site or stop the blood pressure measurements immediately.

12.6 Understanding the NIBP Numerics

The NIBP display shows numerics only as below. Your display may be configured to look slightly different.



- 1. Time of last measurement
- 2. Time remaining to next measurement
- 3. Measurement mode
- 4. Unit of pressure: mmHg or kPa
- 5. Prompt message area: shows NIBP-related prompt messages
- 6. Systolic pressure
- 7. Diastolic pressure
- 8. Mean pressure obtained after the measurement and cuff pressure obtained during the measurement

12.7 Changing NIBP Settings

By selecting the NIBP parameter window, you can enter the [NIBP Setup] menu.

12.7.1 Setting the Initial Cuff Inflation Pressure


You can set the initial cuff inflation pressure manually. In the [**NIBP Setup**] menu, select [**Initial Pressure**] and then select the appropriate setting.

12.7.2 Setting NIBP Alarm Properties

Select [**Alarm Setup >>**] from the [**NIBP Setup**] menu. You can set the alarm properties for this parameter in the popup menu.

12.7.3 Displaying NIBP List

Select [**Screens**] QuickKey→[**Screen Setup**]. You can set [**NIBP List**] to be displayed at the bottom area of the screen. Then, multiple sets of most recent NIBP measurements will be displayed. And PR displayed is derived from NIBP.



You can not display NIBP list in some screens such as the big numerics screen.

12.7.4 Setting the Pressure Unit

Select [**Unit Setup >>**] from the [**User Maintenance**] menu. In the popup menu, select [**Press. Unit**] and toggle between [**mmHg**] and [**kPa**].

12.8 Assisting Venous Puncture

You can use the NIBP cuff to cause sub-diastolic pressure to block the venous blood vessel and therefore help venous puncture.

- 1. Select [**VeniPuncture** >>] from the [**NIBP Setup**] menu. In the popup menu, verify that the [**Cuff Press.**] value is appropriate. Change it if necessary.
- 2. Select [VeniPuncture].
- 3. Puncture vein and draw blood sample.
- 4. Select the shardkey on the monitor's front, or the [**Stop All**] QuickKey to deflate the cuff. The cuff deflates automatically after a set time if you do not deflate it.

During measurement, the NIBP display shows the inflation pressure of the cuff and the remaining time in venous puncture mode.

12.9 NIBP Leakage Test


The NIBP leakage test checks the integrity of the system and of the valve. It is required at least once every two years or when you doubt the measured NIBP. If the test failed, corresponding prompt messages will be given. If no message is displayed, it means no leakage is detected.

Tools required:

- An adult cuff
- An air tubing
- A correct sized cylinder

Follow this procedure to perform the leakage test:

- 1. Set the patient category to [**Adu**].
- 2. Connect the cuff to the NIBP connector on the monitor.
- 3. Wrap the cuff around the cylinder as shown below.



- Select [Main Menu]→[Maintenance >>]→[NIBP Leakage Test]. The NIBP display shows [Leakage Testing...].
- 5. After about 20 seconds, the monitor will automatically deflate. This means the test is completed.
- 6. If the message [**NIBP Pneumatic Leak**] is displayed, it indicates that the NIBP airway may have leakages. Check the tubing and connections for leakages. If you ensure that the tubing and connections are all correct, perform a leakage test again.

If the problem persists, contact your service personnel.

NOTE

The leakage test is intended for use to simply determine whether there are leakages in the NIBP airway. It
is not the same as that specified in the EN 1060-3 standard.

12.10 NIBP Accuracy Test


The NIBP accuracy test is required at least once every two years or when you doubt the measured NIBP.

Tools required:

- T-piece connector
- Approprating tubing
- Balloon pump
- Metal Vessel (volume 500±25 ml)
- Reference manometer (calibrated with accuracy higher than 1 mmHg)

Follow this procedure to perform the accuracy test:

1. Connect the equipment as shown.



- 2. Before inflation, the reading of the manometer should be 0. If not, disconnect the airway and reconnect it until the readings is 0.
- 3. Select [Main Menu]→[Maintenance >>]→[NIBP Accuracy Test].
- 4. Compare the manometer values with the displayed values. The difference between the manometer and displayed values should be no greater than 3 mmHg.
- 5. Raise the pressure in the metal vessel to 50 mmHg with the balloon pump. Repeat step 3 and 4.
- 6. Raise the pressure in the metal vessel to 200 mmHg with the balloon pump. Repeat step 3 and 4.

If the difference between the manometer and displayed values is greater than 3 mmHg, contact your service personnel.

12.11 Calibrating NIBP

NIBP is not user-calibrated. Cuff-pressure transducers must be verified and calibrated once every two years by a qualified service professional. Contact your service personnel when a calibration is necessary.

FOR YOUR NOTES		

13 Monitoring Temp

13.1 Introduction

You can simultaneously monitor two temperature sites using the patient monitor.

13.2 Safety



WARNING

Verify that the probe detection program works correctly before monitoring. Plug out the temperature
probe cable from the T1 or T2 connector, and the monitor can display the message [T1 Sensor Off] or [T2
Sensor Off] and give alarm tones correctly.

13.3 Making a Temp Measurement

- 1. Select an appropriate probe for your patient.
- 2. If you are using a disposable probe, connect the probe to the temperature cable.
- 3 Plug the probe or temperature cable to the temperature connector.
- 4. Attach the probe to the patient correctly.
- 5. Check that the alarm settings are appropriate for this patient.

13.4 Understanding the Temp Display

The temperature monitoring is displayed on the monitor as three numerics: T1, T2 and TD. By selecting this area, you can enter the [**Alarm Setup**] menu.



13.5 Setting the Temperature Unit

Select [**Unit Setup >>**] from the [**User Maintenance**] menu. In the popup menu, select [**Temp Unit]** and toggle between [°**C**] and [°**F**].

14 Monitoring IBP

14.1 Introduction

The iPM 12 patient monitor can monitor up to 4 invasive blood pressures and iPM 10 and iPM 8 patient monitors can monitor up to 2 invasive blood pressures. The patient monitor can display the systolic, diastolic and mean pressures and a waveform for each pressure.

14.2 Safety




WARNING

- Use only pressure transducers specified in this manual. Never reuse disposable pressure transducers.
- Make sure that the applied parts never contact other conductive parts.
- To reduce the hazard of burns during high-frequency surgical procedure, ensure that the monitor's cables and transducers never come into contact with the high-frequency surgical units.
- When using accessories, their operating temperature should be taken into consideration. For details, refer to instructions for use of accessories.

14.3 Zeroing the Transducer

To avoid inaccurate pressure readings, the monitor requires a valid zero. Zero the transducer in accordance with your hospital policy (at least once per day). Zero whenever:

- A new transducer or adapter cable is used.
- You reconnect the transducer cable to the monitor.
- The monitor restarts.
- You doubt the readings.
- 1. Turn off the stopcock to the patient.



- 2. Vent the transducer to the atmospheric pressure by turning on the stopcock to the air.
- 3. In the setup menu for the pressure (e.g. Art), select [**Art Zero** >>]→[**Zero**]. During zero calibration, the [**Zero**] button appears dimmed. It recovers after the zero calibration is completed. To zero all IBP channels, select [**Zero** IBP] hotkey, and then select [**Zero All Channels**] in the popup menu.
- 4. After the zero calibration is completed, close the stopcock to the air and open the stopcock to the patient.

NOTE

 Your hospital policy may recommend that the ICP transducer is zeroed less frequently than other transducers.


14.4 Setting Up the Pressure Measurement

- Plug the pressure cable into the IBP connector.
- Prepare the flush solution. 2.
- 3. Flush the system to exhaust all air from the tubing. Ensure that the transducer and stopcocks are free of air bubbles.



WARNING

- If air bubbles appear in the tubing system, flush the system with the infusion solution again. Air bubble may lead to wrong pressure reading.
- Connect the pressure line to the patient catheter.
- Position the transducer so that it is level with the heart, approximately at the level of the midaxillary line. 5.
- Select the appropriate label.
- 7. Zero the transducer. After a successful zeroing, turn off the stopcock to the atmosphere and turn on the stopcock to the patient.






🗥 WARNING

If measuring intracranial pressure (ICP) with a sitting patient, level the transducer with the top of the patient's ear. Incorrect leveling may give incorrect values.

14.5 Understanding the IBP Display

The IBP measurement is displayed on the monitor as a waveform and numeric pressures. The figure below shows the waveform and numerics for the Art pressure. For different pressures, this display may be slightly different.




- 1. Waveform
- 2. Pressure unit
- 3. Systolic pressure
- 4. Mean pressure
- 5. Diastolic pressure

For some pressures, the parameter window may show the mean pressure only. For different pressures, their defaults unit may be different. If the Art and ICP pressures are measured simultaneously, the ICP parameter area will display numeric CPP, which is obtained by subtracting ICP from the Art mean.

14.6 Changing IBP Settings

14.6.1 Changing a Pressure for Monitoring

1. Select the pressure you want to change to enter its setup menu. In the menu, there is a figure showing the current IBP measurement connector.



2. Select [**Label**] and then select your desired label from the list. The already displayed labels cannot be selected.

Label	Description	Label	Description
PA	Pulmonary artery pressure	CVP	Central venous pressure
Ao	Aortic pressure	LAP	Left atrial pressure
UAP	Umbilical arterial pressure	RAP	Right atrial pressure
BAP	Brachial arterial pressure	ICP	Intracranial pressure
FAP	Femoral arterial pressure	UVP	Umbilical venous pressure
Art	Arterial blood pressure	LV	Left ventricular pressure
P1 to P4 Non-specific pressure label		ssure label	

14.6.2 Setting Alarm Properties

Select [**Alarm Setup >>**] from the parameter setup menu. You can set alarm properties for this parameter in the popup menu.

14.6.3 Changing Averaging Time

The IBP value displayed on the monitor screen is the average of data collected within a specific time. The shorter the averaging time is, the quicker the patient monitor responds to changes in the patient's blood pressure. Contrarily, the longer the averaging time is, the slower the patient monitor responds to changes in the patient's blood pressure, but the measurement accuracy will be improved. For critically ill patients, selecting shorter averaging time will help understanding the patient's state.

To set the averaging time, in the parameter setup menu, select [**Sensitivity**] and toggle between [**High**], [**Med**] and [**Low**], the corresponding averaging time is about 1 s, 8 s and 12 s respectively.

14.6.4 Setting the Pressure Unit

Select [Unit Setup >>] from the [User Maintenance] menu. In the popup menu, select [Press. Unit] and toggle between [mmHg] and [kPa]. Select [CVP Unit] and toggle between [mmHg], [cmH2O] and [kPa].

14.6.5 Setting Up the IBP Wave

In the setup menu for the pressure, you can:

- Select [**Sweep**] and then select the appropriate setting. The faster the wave sweeps, the wider the wave is.
- Select [Scale] and then select the appropriate setting. If [Auto] is selected, the size of the pressure's waveform will be adjusted automatically.
- Select [**Filter**] and then select the desired option.

14.7 Measuring PAWP


Pulmonary Artery Wedge Pressure (PAWP) values, used to assess cardiac function, are affected by fluid status, myocardial contractility, and valve and pulmonary circulation integrity.

Obtain the measurement by introducing a balloon-tipped pulmonary artery flotation catheter into the pulmonary artery. When the catheter is in one of the smaller pulmonary arteries, the inflated balloon occludes the artery allowing the monitor to record changes in the intrathoracic pressures that occur throughout the respiration cycle.

The pulmonary wedge pressure is the left ventricular end diastolic pressure when the airway pressure and valve function are normal. The most accurate PAWP values are obtained at the end of the respiration cycle when the intrathoracic pressure is fairly constant and the artifact caused by respiration is minimal.

14.7.1 Preparing to Measure PAWP

- 1. Prepare the same accessories as in the C.O. measurement. Connect the parts such as catheter, syringe, etc. following the C.O. measurement steps and use the balloon inflation port.
- 2. Connect the PAWP cable into the IBP connector on the monitor. Since PAWP is measured on PA, selecting [**PA**] as the IBP label is recommended.
- Select the PA parameter window or waveform area to enter its setup menu. Then, select [PAWP] to enter the PAWP measurement window. You can also enter the PAWP measurement window from the P1-P4 parameter window.




NOTE

After entering the PAWP measurement window, the monitor will turn off the PA alarm automatically.

14.7.2 Setting Up the PAWP Measurement

1. Wedge the flotation catheter into the pulmonary artery. Then inflate the balloon and pay attention to PA waveform changes on the screen.



- 2. After obtaining a stable PAWP waveform, press the [Freeze] key to freeze the waveform and deflate the balloon.
- 3. You can adjust the PAWP scale to an appropriate position by adjusting or beside the [**Adjust**] button. Press the [**Confirm**] key to save one PAWP measurement.
- 4. If you need to start a new measurement, select [**Next Measure**].

WARNING

- Prolonged inflation can cause pulmonary hemorrhage, infarction or both. Inflate the balloon for the minimum time necessary to get an accurate measurement.
- If the PAWP is greater than the PA (systolic), deflate the balloon and report the incident in accordance with hospital policy. Because the pulmonary artery could be accidentally ruptured, and the PAWP value derived will not reflect the patient's hemodynamic state, but will merely reflect the pressure in the catheter or balloon.

14.7.3 Understanding the PAWP Setup Menu

Select [**Setup**] to enter the [**PAWP Setup**] menu. In this menu, you can:

- Select a ECG lead wave as the first reference wave.
- Select a respiration wave as the second reference wave.
- Select a sweep speed for the displayed waveform.
- Change the size of the PA waveform by adjusting the scale height.

FOR YOUR NOTES		

15 Monitoring Cardiac Output

15.1 Introduction

The cardiac output (C.O.) measurement invasively measures cardiac output and other hemodynamic parameters using the right heart (atria) thermodilution method. A cold solution of known volume and temperature is injected into the right atrium through the proximal port of a pulmonary artery (PA) catheter. The cold solution mixes with the blood in the right ventricle and the change in blood temperature is measured with a thermistor at the distal end of the catheter in the pulmonary artery. The temperature change is displayed as a curve in the C.O. split screen, and the monitor calculates the C.O. value from this curve. The C.O. value is inversely proportional to the area under the curve. As cardiac output varies continuously, a series of measurements must be carried out to achieve a reliable C.O. average value. Always use the average of multiple thermodilution measurements for therapy decisions. The monitor is capable of storing 6 measurements.



• C.O. monitoring is restricted to adult patients only.

15.2 Understanding the C.O. Display

The C.O. measurement is displayed on the monitor as numeric C.O., C.I. and TB in the C.O. parameter window as shown below. To enter the **[C.O. Setup**] menu, select the C.O. parameter window.



- 1. Cardiac output
- 2. Time at which the C.O. average is calculated
- 3. Cardiac index
- 4. Blood temperature

15.3 Influencing Factors

The factors that affect cardiac output are:

- temperature of injectate solution,
- volume of injectate solution,
- patient's baseline blood temperature,
- patient's inspiratory/expiratory cycle,
- placement of catheter with relation to proximity of lung field,
- the catheter itself,
- the patient rhythm and hemodynamic status, and
- any other rapid IV solutions which are infused while the C.O. measurement is being performed.


Followings are some technique suggestions to obtain accurate C.O.:

- Injectate solution must be cooler than the patient's blood.
- Inject solution rapidly and smoothly.
- Inject at end expiration.


15.4 Setting Up the C.O. Measurement



- Use only accessories specified in this manual. Make sure that the accessories never come into contact with conductive parts.
- 1. Connect the C.O. cable to the C.O. connector on the monitor.
- 2. Interconnect the C.O. module, catheter and syringe as shown below. Make sure that:
 - ◆ The module is securely inserted.
 - ◆ The PA catheter is in place in the patient.
 - ◆ The C.O. cable is properly connected to the module.



- 3. Select the C.O. parameter window to enter the [**C.O. Setup**] menu. Check if the height and weight are appropriate for your patient. Change if necessary.
- 4. In the [**C.O. Setup**] menu:
 - Check that the correct computation constant is entered. To change the computation constant, select [**Comp. Const**] and then enter the correct value. When a new catheter is used, the computation constant should be adjusted in accordance with the manufacturer's instructions for use.
 - ◆ Set the [Auto TI] to [Manual] or [Auto]. If you select [Auto], the system automatically detects the injectate temperature, and the [Manual TI] is disabled. If you select [Manual], you need to enter the injectate temperature at [Manual TI] manually.
- 5. Select [**C.O. Measure**] to enter the C.O. measurements window.



- A. Currently measured numeric
- B. Currently measured C.O. curve
- C. Prompt message area
- D. Buttons
- E. Averaged values
- F. Measurement windows
- 6. Proceed as follows.

When you see the message [**Ready for new measurement**], select the [**Start**] button and then inject the solution quickly when the message [**Inject now!**] and prompt tone appear. As shown in the figure above, during the measurement, the currently measured thermodilution curve and message [**Measuring...**] are displayed. At the end of the measurement, the thermodilution curve is transferred to one of the 6 measurement windows and the monitor prompts you to wait for a certain period of time before starting a new measurement.

A maximum of 6 measurements can be stored. If you perform more than six measurements without rejecting any, the oldest will automatically be deleted when a seventh curve is stored. Select from the 6 measurement curves and the system will automatically calculate and display the averaged C.O. and C.I. values. Then select the [**Accept Average**] button to accept and store the averaged values.

When injecting, the stopcock to the PA catheter is open and the stopcock to the injectate solution is closed. After the measurement is completed, turn off the stopcock to the PA catheter and turn on the stopcock to the injectate solution, and then draw the injectate solution into the injectate syringe.

In the buttons area, you can:

- Select [Cancel] during a measurement to cancel the measurement. Selecting it after a measurement deletes the measured results.
- Select [Record] to print out the curves selected for average calculation, numerics and averaged values by the recorder.
- Select [**Setup** >>] to access the [**C.O. Setup**] menu.
- Select [Calc >>]→[Hemodynamic >>] to access the [Hemodynamic Calculation] menu.


The system can automatically adjust the X-axis scale range to 30 s or 60 s and Y-axis scale range to 0.5 ℃, 1.0 ℃, or 2.0 ℃.

NOTE

- Starting measuring without blood temperature being stable yet may cause measuring failure.
- During the cardiac output measurement, blood temperature alarms are inactive.

15.5 Measuring the Blood Temperature

As shown below, the blood temperature is measured with a temperature sensor at the distal end of the catheter in the pulmonary artery. During C.O. measurements, blood temperature alarms are suppressed to avoid false alarms. They will automatically recover as soon as the C.O. measurements are completed.



15.6 Changing C.O. Settings

15.6.1 Setting the Temperature Unit

Select [**Unit Setup** >>] from the [**User Maintenance**] menu. In the popup menu, select [**Temp Unit**] to toggle between [$^{\circ}$ C] and [$^{\circ}$ F].

15.6.2 Setting Alarm Properties

Select [**Alarm Setup >>**] from the [**C.O. Setup**] menu. You can set alarm properties for this parameter in the popup menu.

FOR YOUR NOTES		

16 Monitoring Carbon Dioxide

16.1 Introduction


 CO_2 monitoring is a continuous, non-invasive technique for determining the concentration of CO_2 in the patient' airway by measuring the absorption of infrared (IR) light of specific wavelengths. The CO_2 has its own absorption characteristic and the amount of light passing the gas probe depends on the concentration of the measured CO_2 . When a specific band of IR light is passed through respiratory gas samples, some of IR light will be absorbed by the CO_2 molecules. The amount of IR light transmitted after it has been passed through the respiratory gas sample is measured with a photodetector. From the amount of IR light measured, the concentration of CO_2 is calculated.

There are two methods for measuring CO₂ in the patient's airway:

- Mainstream measurement uses a CO₂ sensor attached to an airway adapter directly inserted into the patient's breathing system.
- 2. Sidestream/Microstream measurement samples expired patient gas at a constant sample flow from the patient's airway and analyzes it with a CO₂ sensor built into the CO₂ module.

The measurement provides:


- 1. A CO₂ waveform
- 2. End tidal CO₂ value (EtCO₂): the CO₂ value measured at the end of the expiration phase.
- 3. Fraction of inspired CO₂ (FiCO₂): the smallest CO₂ value measured during inspiration.
- 4. Airway respiration rate (awRR): the number of breaths per minute, calculated from the CO₂ waveform.



16.2 Measuring CO₂

16.2.1 Making a Sidestream CO₂ Measurement

1. Attach the watertrap to the module and then connect the CO₂ components as shown below.



- The CO₂ module needs time to warm up to reach the operating temperature. The message [CO₂ Sensor Warmup] is displayed during warm-up. If you perform CO₂ measurements during warm-up, the measurement accuracy may be compromised.
- 3. After warm-up is finished, you can perform CO₂ measurements.

ACAUTION


- The watertrap collects water drops condensed in the sampling line and therefore prevents them from
 entering the module. If the collected water reaches a certain amount, you should drain it to avoid blocking
 the airway.
- The watertrap has a filter preventing bacterium, water and secretions from entering the module. After a long-term use, dust or other substances may compromise the performance of the filter or even block the airway. In this case, replace the watertrap. It is recommended to replace the watertrap once every two months, or when the watertrap is found leaky, damaged or contaminated.

NOTE

• To extend the lifetime of the watertrap and module, disconnect the watertrap and set the operating mode to standby mode when CO₂ monitoring is not required.


16.2.2 Making a Microstream CO₂ Measurement

Connect the sampling line to the module and then connect the CO2 components as shown below. After warm-up is finished, you can perform CO2 measurements.



16.2.3 Making a Mainstream CO₂ Measurement

- 1. Connect the sensor to the module. The message [CO₂ Sensor Warmup] appears on the screen when the CO₂ module is plugged.
- 2. After warm-up is finished, connect the transducer to the airway adapter.
- 3. Perform a zero calibration per the **Zeroing the Sensor** section.
- 4. After the zero calibration is finished, connect the airway as shown below.



5. Make sure there are no leakages in the airway and then start a measurement.

NOTE

 Always position the sensor with the adapter in an upright position to avoid collection of fluids on the windows of the adapter. Large concentrations of fluids at this point will obstruct gas analysis.

16.3 Changing CO₂ Settings

16.3.1 Accessing CO₂ Menus

By selecting the CO₂ parameter window, you can access the [CO₂ Setup] menu.

16.3.2 Entering the Standby Mode

The standby mode of the CO₂ module relates to the standby mode of the monitor as follows:

- If the monitor enters the standby mode, the CO₂ module will also enter the standby mode.
- If the monitor exits the standby mode, the CO₂ module will also exit the standby mode.
- If the CO₂ module enters or exits the standby mode, it will not affect the monitor.

To enter or exit the standby mode manually, select [**Operating Mode**] in the [**CO**₂ **Setup**] menu and then toggle between [**Standby**] and [**Measure**].

When you set the sidestream CO_2 module to the strandby mode, the CO_2 gas sample intake pump automatically sets the sample flow rate to zero. When exiting the standby mode, the CO_2 module continues to work at the preset sample flow rate.

For the sidestream CO_2 module, you can set the delay time. After the delay time the CO_2 module enters the standby mode if no breath is detected.

For the microstream CO_2 module, you can also set a period of time after which the CO_2 module enters the standby mode if no breath is detected since the CO_2 module is powered on or the CO_2 module switches to the measuring mode or the automatic standby time is re-set. To set the standby time, in the [CO_2 Setup] menu, select [Auto Standby] and then select the appropriate setting.

16.3.3 Setting the CO₂ Unit

Select [Unit Setup >>] from the [User Maintenance] menu. In the popup menu, select [CO₂. Unit] and toggle between [mmHg], [%] and [kPa].

16.3.4 Setting up Gas Compensations

NWARNING

 Make sure that the appropriate compensations are used. Inappropriate compensations may cause inaccurate measurement values and result in misdiagnosis.

For the sidestream CO₂ module:

- 1. Select [**CO**₂ **Setup**].
- 2. According to the actual condition, set the concentration required for the following compensations:
 - ♦ [O₂ Compen]
 - ♦ [N₂O Compen]
 - ♦ [Des Compen]

For the microstream CO₂ module, gas compensations are not required.

For the mainstream CO_2 module, in the [CO_2 Setup] menu, respectively select:

- **Balance Gas**] and toggle between [**Room Air**] and [N_2O]. Select [**Room Air**] when air predominates in the ventilation gas mixture and select [N_2O] when N_2O predominates in the ventilation gas mixture and select [N_2O] when N_2O predominates in the ventilation gas mixture.
- **O**₂ **Compen**] and then select [**Off**] or an appropriate setting according to the amount of O_2 in the ventilation gas mixture. When the amount of O_2 is less than 30%, you'd better switch this compensation off.
- [AG Compen] and enter the concentration of anesthetic gas present in the ventilation gas mixture. This could compensate for the effect of AG on the readings.

16.3.5 Setting up Humidity Compensation

Sidestream and microstream CO_2 modules are configured to compensate CO_2 readings for either Body Temperature and Pressure, Saturated Gas (BTPS), to account for humidity in the patient's breath, or Ambient Temperature and Pressure, Dry Gas (ATPD).

- 1. ATPD: $P_{co2}(mmHg) = CO_2(vol\%) \times P_{amb} / 100$
- 2. BTPS: $P_{CO2}(mmHg) = CO_2(vol\%) \times (P_{amb} 47)/100$

Where, P_{CO2} = partial pressure, vol% = CO₂ concentration, P_{amb} = ambient pressure, and unit is mmHg.

As the mainstream CO₂ module has a built-in heating component to prevent water vapour from condensing, setting humidity compensation is not needed. For the sidestream and microstream CO₂ module, you can set the humidity compensation on or off according to the actual condition. To set the humidity compensation:

- 1. In the [CO₂ Setup] menu, select [BTPS Compen].
- 2. Select either [**On**] for BTPS or [**Off**] for ATPD, depending on which compensation applies.

16.3.6 Setting the Apnea Alarm Delay

In the [CO₂ Setup] menu, select [Apnea Delay] and then select the appropriate setting. The monitor will alarm if the patient has stopped breathing for longer than the preset apnea time. The [Apnea Delay] of Resp, CO₂, and AG module keeps consistent with each other.

WARNING

The respiration measurement does not recognize the cause of apneas. It only indicates an alarm if no
breath is detected when a pre-adjusted time has elapsed since the last detected breath. Therefore, it
cannot be used for diagnostic purpose.

16.3.7 Choosing a Time Interval for Peak-Picking

For microstream and mainstream CO_2 modules, you can select a time interval for picking the highest CO_2 as the $EtCO_2$ and the lowest as the $FiCO_2$.

In the [CO₂ Setup] menu, select [Max Hold] and toggle between [Single Breath], [10 s], [20 s] and [30 s] (for microstream CO₂ module only)..

- [Single Breath]: EtCO₂ and FiCO₂ are calculated for every breath.
- [10 s], [20 s] or [30 s]: EtCO₂ and FiCO₂ are calculated using 10, 20 or 30 seconds of data.

16.3.8 Setting the Flow Rate

For the sidestream CO_2 module, you can change the sampling rate of respiratory gas in the patient's airway by setting the flow rate. To set the flow rate, enter the [CO_2 Setup] menu and select an appropriate setting from [Flow Rate].



 Please consider the patient's actual bearing capability and select the appropriate flow rate when setting the flow rate.

16.3.9 Setting up the CO₂ Wave

In the [CO₂ Setup] menu, you can:

- Select [Wave Type] and toggle between [Draw] and [Fill]:
 - ◆ [**Draw**]: The CO₂ wave is displayed as a curved line.
 - ◆ [Fill]: The CO₂ wave is displayed as a filled area.
- Select [**Sweep**] and then select the appropriate setting. The faster the wave sweeps, the wider the wave is.
- Change the size of the CO₂ waveform by adjusting the wave [**Scale**].

16.4 Setting RR Source

To set RR source:

- 1. Enter the [CO₂ Setup] menu.
- 2. Select [RR Source] and then select a source or [Auto] from the dropdown list.

The [RR Source] settings of Resp, CO₂, and AG module are linked. For details, please refer to the section **Setting RR Source** of chapter **Resp**.

16.5 Setting Barometric Pressure Compensation

Both sidestream and microstream CO_2 modules have the function of automatic barometric pressure compensation (the system automatically measures the barometric pressure which the patient monitor is exposed to). However, the mainstream CO_2 module does not have such function. For the mainstream CO_2 module, the default barometric pressure is 760 mmHg. You must modify the barometric pressure based on the actual situation as follows:

- Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password→[Maintain CO₂ >>]→[Calibrate CO₂ >>].
- Select [Barometric Pressure] and then enter the value of barometric pressure to which the patient monitor is exposed to.



 Be sure to set the barometric pressure properly before using the mainstream CO₂ module. Improper settings will result in erroneous CO₂ reading.

16.6 Measurement Limitations

The following factors may influence the accuracy of measurement:

- Leaks or internal venting of sampled gas
- Mechanical shock
- Cyclic pressure up to 10 kPa (100 cmH₂O)
- Other sources of interference, if any

16.7 Leakage test

When the sidestream CO_2 module needs maintenance, the monitor prompts on the CO_2 waveform area: [**Need** maintenance. Enter CO_2 setup menu.] Then, you can access [CO_2 Setup] \rightarrow [Maintain CO_2], and perform leakage test according to the prompt messages on the menu.

16.8 Troubleshooting the Sidestream CO₂ Sampling System

When the sampling system of the sidestream CO_2 module works incorrectly, check if the sampling line is kinked. If not, remove it from the watertrap. If the monitor gives a message indicating the airway still works incorrectly, it indicates that the watertrap must have been blocked, and you should replace with a new one. Otherwise, you can determine that the sampling line must have been blocked. Replace with a new sampling line.

16.9 Removing Exhaust Gases from the System



Anesthetics: When using the Sidestream or Microstream CO₂ measurement on patients who are receiving
or have recently received anesthetics, connect the outlet to a scavenging system, or to the anesthesia
machine/ventilator, to avoid exposing medical staff to anesthetics.

To remove the sample gas to a scavenging system, connect an exhaust tube to the gas outlet connector of the module.

16.10 Zeroing the Sensor

The zero calibration eliminates the effect of baseline drift during CO_2 measurement exerted on the readings and therefore maintains the accuracy of the CO_2 measurements.

16.10.1 For Sidestream and Microstream CO₂ Modules

For sidestream and microstream CO_2 modules, a zero calibration is carried out automatically when necessary. You can also start a manual zero calibration if necessary. To manually start a zero calibration, select [**Maintain CO₂** >>] from the [**User Maintenance**] menu. Then select [**Calibrate CO₂** >>] \rightarrow [**Start Zero Cal.**]. Disconnecting the patient airway is not required when performing a zero calibration.

16.10.2 For Mainstream CO₂ Modules

For mainstream CO₂ modules, zero the sensor whenever:

- A new adapter is used;
- You reconnect the sensor to the module;
- You see the message [CO₂ Zero Required]. In this case, check the airway adapter for any blockage, e.g. mucus, etc. If a blockage is detected, clear or replace the adapter.

To zero the sensor, follow this procedure:

- 1. Connect the sensor to the module.
- 2. In the [CO₂ Setup] menu, set the [Operating Mode] to [Measure]. The message [CO₂ Sensor Warmup] is displayed.
- 3. After warm-up is finished, connect the sensor to a clean, dry airway adapter. The adapter should be vented to the air and isolated from CO_2 sources, such as ventilator, the patient's breathing, your own breathing, etc.
- 4. Select [Start Zero Cal.] in the [CO₂ Setup] menu. The message [CO₂ Zero Running] is displayed.
- 5. It takes about 15 to 20 seconds. The message disappears when the zero calibration is completed.



- When perform a zero calibration during the measurement, disconnect the transducer from the patient's airway first.
- Please do not rely on the readings during zeroing.

16.11 Calibrating the Sensor

For sidestream or microstream CO_2 modules, a calibration should be performed once every year or when the readings go far beyond the range. For mainstream CO_2 modules, no calibration is required. For details, refer to the chapter **26** *Maintenance*.

16.12 Oridion Information

Microstream

This trademark is registered in Israel, Japan, German and America.

Oridion Patents

The capnography component of this product is covered by one or more of the following US patents: 6,428,483; 6,997,880; 6,437,316; 7,488,229; 7,726,954 and their foreign equivalents. Additional patent applications pending.

No Implied License

Possession or purchase of this device does not convey any express or implied license to use the device with unauthorized CO_2 sampling consumables which would, alone, or in combination with this device, fall within the scope of one or more of the patents relating to this device and/or CO_2 sampling consumable.

17 Monitoring AG

17.1 Introduction

The anaesthetic gas (AG) module measures the patient's anesthetic and respiratory gases, and incorporates the features of the O₂ module as well. AG monitoring is for iPM 12 and iPM 10 patient monitors only.

The AG module determines the concentration of certain gases using the infrared (IR) light absorption measurement. The gases that can be measured by the AG module absorb IR light. Each gas has its own absorption characteristic. The gas is transported into a sample cell, and an optical IR filter selects a specific band of IR light to pass through the gas. For multiple gas measurement, there are multiple IR filters. The higher the concentration of gas in a given volume the more IR light is absorbed. This means that higher concentration of IR absorbing gas cause a lower transmission of IR light. The amount of IR light transmitted after it has been passed though an IR absorbing gas is measured. From the amount of IR light measured, the concentration of gas present can be calculated.


Oxygen does not absorb IR light as other breathing gases and is therefore measured relying on its paramagnetic properties. Inside the O_2 sensor are two nitrogen-filled glass spheres mounted on a strong rare metal taut-band suspension. This assembly is suspended in a symmetrical non-uniform magnetic field. In the presence of paramagnetic oxygen, the glass spheres are pushed further away from the strongest part of the magnetic field. The strength of the torque acting on the suspension is proportional to the oxygen concentration. From the strength of the torque, the concentration of oxygen is calculated.

AG module can identify two anesthetic gases in a mixture automatically and distinguish between them according to their contributions to the MAC value for display as the primary and secondary anesthetis agent.

NOTE

The AG module is configured with automatic barometric pressure compensation function.

17.2 Understanding the AG Display



The AG module can send waves and numerics for all measured anesthetic gases for display on the monitor, including:

- CO₂, O₂, N₂O and AA waves
- awRR: airway respiratory rate
- MAC: minimum alveolar concentration
- End tidal (Et) and fraction of inspired (Fi) numerics for CO₂, O₂, N₂O and AA

Where AA represents Des (desflurane), Iso (isoflurane), Enf (enflurane), Sev (sevoflurane), or Hal (halothane).

The AA waveform area displays the primary anesthetic gas's waveform. When O_2 module does not exist, no O_2 waveform will be displayed. When O_2 module exists, the O_2 waveform will be displayed only when the O_2 waveform is currently switched on.



• To avoid explosion hazard, do not use flammable anesthetic agent such as ether and cyclopropane for this equipment.

17.3 MAC Values

Minimum alveolar concentration (MAC) is the minimum concentration of the agent in the alveoli. It is a basic index to indicate the depth of anesthesia. The standard ISO 21647 defines MAC as this: alveolar concentration of an inhaled anesthetic agent that, in the absence of other anesthetic agents and at equilibrium, prevents 50% of patients from moving in response to a standard surgical stimulus.

Minimum alveolar concentration (MAC) values are listed below:

Agent	Des	Iso	Enf	Sev	Hal	N2O
1 MAC	7.3%*	1.15%	1.7%	2.1%	0.77%	105%**

^{*} The data is taken from a patient of 25 years old.

^{**} indicates 1 MAC nitrous oxide can only be reached in hyperbaric chamber.

NOTE

- The MAC values shown in the table above are those published by the U.S. Food and Drug Administration for a healthy 40-year-old adult male patient.
- In actual applications, the MAC value may be affected by age, weight and other factors.

The formula to calculate the MAC value is as follows:

$$MAC = \sum_{i=0}^{N-1} \frac{EtAgent_i}{AgentVol_i}$$

Where N is the number of all agents (including N_2O) that the AG module can measure, EtAgenti is the concentration of each agent and AgentVoli is the concentration of each agent at 1 MAC.

For example, the AG module measures there are 4% of Des, 0.5% of Hal and 50% of N₂O in the patient's end-tidal gas:


$$MAC = \frac{4.0\%}{7.3\%} + \frac{0.5\%}{0.77\%} + \frac{50\%}{105\%} = 1.67$$

NOTE

• The formula mentioned above is intended for adult patients only.

17.4 Preparing to Measure AG

- 1. Select an appropriate watertrap according to patient category and attach it to the module.
- 2. Connect the gas sample line to the connector of the watertrap.
- 3. Connect the other end of the gas sampling line to the patient via the airway adapter.
- 4. Connect the gas outlet to a scavenging system using an exhaust tube.



5. Insert the AG module to the patient monitor and the patient monitor will prompt [**AG Startup**]. Then the AG module starts to warmup and at the same time the patient monitor prompts [**AG Warmup**]. After 45 seconds, the AG module enters the iso accuracy mode. After 10 minutes, the module enters the full accuracy mode.

Awarning

- Make sure that the connections are tight. Any leak in the system can result in erroneous readings due to ambient air mixing with patient gases.
- Do not apply adult watertrap to the neonate patient. Otherwise, patient injury could result.
- Using high-frequency electrosurgical units may increase the risk of skin burn. In this case, do not use antistatic or conductive respiratory tubing.

ACAUTION

- Position the airway adapter so that the part connecting to the gas sample line is pointing upwards. This
 prevents condensed water from passing into the gas sample line and causing an occlusion.
- The watertrap collects water drops condensed in the sampling line and therefore prevents them from entering the module. If the collected water reaches to a certain amount, you should drain it to avoid blocking the airway.
- The watertrap has a filter preventing bacterium, water and secretions from entering the module. After a
 long-term use, dust or other substances may compromise the performance of the filter or even block the
 airway. In this case, replace the watertrap. Replacing the watertrap once a month is recommended.

17.5 Changing AG Settings

17.5.1 Setting Gas Unit

For N2O and AA, the unit of the measured gas is fixed to "%".

Select [Unit Setup >>] from the [User Maintenance] menu. In the popup menu, you can select [CO₂ Unit] or [O₂ Unit] and toggle between [mmHg], [%] and [kPa].

17.5.2 Setting the Apnea Alarm Delay

In the [AG Setup] menu, select [Apnea Delay] and select the appropriate setting. The monitor will alarm if the patient has stopped breathing for longer than the preset apnea time.

The [Apnea Delay] of Resp, CO₂, AG, and RM module keeps consistent with each other.

!WARNING

The respiration measurement does not recognize the cause of apneas. It only indicates an alarm if no
breath is detected when a pre-adjusted time has elapsed since the last detected breath. Therefore, it
cannot be used for diagnostic purpose.

17.5.3 Changing the Sample Flow Rate

In the setup menu for any gas, select [Flow Rate] and then choose either:

- [High]: 200 ml/min for adult and pediatric patients, and 120 ml/min for neonatal patients.
- [Med]: 150 ml/min for adult and pediatric patients, and 90 ml/min for neonatal patients.
- [Low]: 120 ml/min for adult and pediatric patients, and 70 ml/min for neonatal patients.

17.5.4 Setting up the O₂ Compensation

If the AG module does not incorporate the O_2 module, you need to manually select [\mathbf{O}_2 **Compen**] and then select [\mathbf{Off}] or an appropriate setting according to the amount of O_2 in the ventilation gas mixture. When the amount of O_2 is less than 30%, you'd better switch this compensation off.

If the AG module incorporates the O_2 module, the system will directly use the O_2 concentration detected by the O_2 module to make compensation. At this time, the [O_2 Compen] in the setup menu for any gas is fixed to [Off].

17.5.5 Entering the Standby Mode

For the AG module, the default operating mode is measure. When you set the AG module to the standby mode, the AG gas sample intake pump automatically sets the sample flow rate to zero. When exiting the standby mode, the AG module continues to work at preset sample flow rate with no need to warm up again. After nearly 1 minute, the module enters the full accuracy mode. The standby mode of the AG module relates to the standby mode of the monitor as follows:

- If the monitor enters the standby mode, the AG module will also enter the standby mode.
- If the monitor exits the standby mode, the AG module will also exit the standby mode.
- If the AG module enters or exits the standby mode, it will not affect the monitor.

To enter or exit the standby mode manually, in the agent's setup menu, select [**Operating Mode**] and then toggle between [**Standby**] and [**Measure**]. You can also set a period of time after which the AG module enters the standby mode automatically if no breath is detected since the last detected breath. To set the standby time, in the agent's setup menu, select [**Auto Standby (min)**] and then select the appropriate setting.

17.5.6 Setting up the AG Wave

In the [**AG Setup**] menu, you can:

- Select [CO₂ Wave Type] and toggle between [Draw] and [Fill]:
 - ♦ [**Draw**]: The CO₂ wave is displayed as a curved line.
 - lacktriangle [**Fill**]: The CO₂ wave is displayed as a filled area.
- Select [**Sweep**] and then select the appropriate setting. The faster the wave sweeps, the wider the wave is.
- Change the size of the waveform by adjusting the scale.

17.5.7 Setting RR Source

To set RR source:

- 1. Enter the [AG Setup] menu.
- 2. Select [RR Source] and then select a source or [Auto] from the dropdown list.

The [RR Source] settings of Resp, CO₂, and AG module are linked. For details, please refer to the section **Setting RR Source** of chapter **Resp**.

17.6 Changing the Anesthetic Agent

When the anesthetic agent used on the patient is changed, the AG module can detect the mixed anesthetic gas during the transition of two anesthetic agents. The time required for completing the replacement of anesthetic agent depends on anesthesia type (low flow or high flow) and the characteristics of anesthetic agents (pharmacokinetics). During the transition of two anesthetic agents, the patient monitor gives no prompt messages and the MAC value displayed may be inaccurate.

The AG module can identify two anesthetic agents automatically. When the proportion of the primary and secondary anesthetic agents in the mixture changes, the AG module can distinguish between them according to their contributions to the MAC value. Then the primary and secondary anesthetic agents will be exchanged for display.

17.7 Measurement Limitations

The following factors may influence the accuracy of measurement:

- Leaks or internal venting of sampled gas
- Mechanical shock
- Cyclic pressure up to 10 kPa (100 cmH₂O)
- Other sources of interference, if any

17.8 Troubleshooting

17.8.1 When the Gas Inlet is Blocked

If the gas inlet (including watertrap, sampling line and airway adapter) is occluded by condensed water, the message [**AG Airway Occluded**] will appear.

To remove the occlusion:

- Check the airway adapter for an occlusion and replace if necessary.
- Check the sampling line for an occlusion or kinking and replace if necessary.
- Check the watertrap for a build up of water. Empty the watertrap. If the problem persists, replace the watertrap.

17.8.2 When an Internal Occlusion Occurs

Condensed water may enter the module and cause contamination and/or internal occlusions. In this case, the message [**AG Airway Occluded**] will be displayed.

To remove the occlusion:

- Check for any occlusion in the gas inlet and/or outlet system.
- If the problem persists, internal occlusions may exist. Contact your service personnel.

17.9 Removing Exhaust Gases from the System



Anesthetics: When using the AG measurement on patients who are receiving or have recently received
anesthetics, connect the outlet to a scavenging system, or to the anesthesia machine/ventilator, to avoid
exposing medical staff to anesthetics.

To remove the sample gas to a scavenging system, connect an exhaust tube to the gas outlet connector of the module.


FOR YOUR NOTES		

18 Freezing Waveforms

During patient monitoring, the freeze feature allows you to freeze the currently displayed waveforms on the screen so that you can have a close examination of the patient's status. Besides, you can select any frozen waveform for recording.

18.1 Freezing Waveforms

- 1. To freeze waveforms, select the hardkey on the monitor's front.
- 2. The system closes the displayed menu (if any), and opens the [Freeze] menu.



3. All displayed waveforms are frozen, i.e. the waveforms stop being refreshed or scrolling.

The freeze feature exerts no effect on the split-screen view of minitrends, oxyCRG and other patients.

18.2 Viewing Frozen Waveforms

To view the frozen waveforms, you can either:

- Select the [Scroll] button and then rotate the Knob clockwise or counter-clockwise, or
- Directly select the ◀ or ▶ beside the [**Scroll**] button using the touchscreen.

The frozen waveforms will scroll left or right accordingly. And meanwhile, at the lower right corner of the bottommost waveform, there is an upward arrow. The freeze time is displayed below the arrow and the initial frozen time is [0 s]. With the waveforms scrolling, the freeze time changes at intervals of 1 second. This change will be applied for all waveforms on the screen.

18.3 Unfreezing Waveforms

To unfreeze the frozen waveforms, you can either:

- Select the × button at the upper right corner of the [Freeze] menu,
- Select the M hardkey on the monitor's front, or
- Perform any other action that causes the screen to be readjusted or opens a menu, such as plugging in or out a module, pressing the hardkey, etc.

18.4 Recording Frozen Waveforms

1.	In the [Freeze] menu, select, in turn, [Wave 1], [Wave 2] and [Wave 3] and then select your desired
	waveforms.

2.	Select the [Record] button. The selected waveforms and all numerics at the frozen time are printed out by the
	recorder.


19 Review

19.1 Accessing Respective Review Windows

- Select the [Review] QuickKey, or [Main Menu]→[Review >>].
- 2. Select [**Graphic Trends**], [**Tabular Trends**], [**Events**], [**Full Disclosure**] or [**12-lead ECG**] to access their respective review windows.

19.2 Reviewing Graphic Trends

In the [Review] menu, select [Graphic Trends] to access the following window.



- 1. Event mark area
- 2. Time axis
- 3. Graphic trends area

- 4. Parameter area
- 5. Cursor

Events are marked with colors in the event mark area. Red represents high level alarm event. Yellow represents medium/low level alarm event. Green represents manual event.

In this review window:

- Select [Trend Group] and you can select a trend group for viewing in the popup menu. If [Custom 1] or [Custom 2] is selected, you can further select [Define Trend Group]. Then you can select the parameters for viewing in the popup menu.
- You can set the time length of the review window by selecting [Zoom].
- You can set the number of waves displayed in one page by selecting [Waves].

- To browse the graphic trends, you can either:
 - ◆ Select or beside [**Scroll**] to move the cursor one step to the left or right to navigate through the graphic trends, or
 - Select or both to move the cursor one page to the left or right to navigate through the graphic trends.

A time indicating your current position is displayed above the parameter area. Numeric measurement values corresponding to the cursor location change as the cursor is moved. The measurement value that triggered high level alarm has red background. The one that triggered medium/low level alarm has yellow background.

- By selecting or beside [**Event**], you can position the cursor to different event time.
- By selecting the [**Record**] button, you can print out the currently displayed graphic trends by the recorder.
- By selecting the [**Print**] button, you can set and print out the graphic trends report by the printer. For how to set the graphic trends report, please refer to the **Print** chapter.

19.3 Reviewing Tabular Trends

In the [Review] menu, select [Tabular Trends] to access the following window.



Events are marked with colors in window's top area. Red represents high level alarm event. Yellow represents medium/low level alarm event. Green represents manual event.

In this review window:

- Select [Trend Group] and you can select a trend group for viewing in the popup menu. If [Custom 1] or [Custom 2] is selected, you can further select [Define Trend Group]. Then you can select the parameters for viewing in the popup menu.
- You can change the resolution of the trend data by selecting [Interval] and then selecting the appropriate setting:
 - ◆ [5 s] or [30 s]: select to view up to 4 hours of tabular trends at 5- or 30-second resolution.

- ♦ [1 min], [5 min], [10 min], [15 min], [30 min], [1 h], [2 h] or [3 h]: select to view up to 120 hours of tabular trends at your selected resolution.
- [NIBP]: select to view the tabular trends when NIBP measurements were acquired.
- To browse the tabular trends, you can either:
 - ♦ Select or beside [**Scroll**] to drag the scrollbar left or right to navigate through the trend database, or
 - Select or to scroll left or right to navigate through the trend database.

The measurement value that triggered high level alarm has red background. The one that triggered medium/low level alarm has yellow background.

- By selecting or beside [**Event**], you can position the cursor to different event time.
- By selecting the [**Record**] button, you can access the [**Record Setup**] menu and set the start and end time of the tabular trends you want to record. This feature is not available when reviewing a history patient. By further selecting [**Record**], you can print out the currently displayed tabular trends by the recorder.
- By selecting the [**Print**] button, you can set and print out the tabular trends report by the printer. For how to set the tabular trends report, please refer to the **Print** chapter.

19.4 Events

19.4.1 Marking Events

During patient monitoring, some events may exert effects on the patient and as a result change the waveforms or numerics displayed on the monitor. To help analysing the waveforms or numerics at that time, you can mark these events.

Select [Main Menu]→[Mark Event >>]. In the popup menu, you can select the waves to be stored when a manual event is triggered. You can select [Trigger Manual Event] from the [Mark Event] menu or the [Manual Event] QuickKey to trigger a manual event and store it at the same time.

When you are reviewing graphic trends, tabular trends or full-disclosure waveforms, the manual event symbol is displayed at the time the event is triggered.

19.4.2 Reviewing Events

In the [Review] menu, select [Events] to access the following window.


The events that can be reviewed include parameter alarm events, arrhythmia alarm events and manual events. When an event occurs, all the measurement numerics at the event trigger time and related waveforms 4 seconds respectively before and after the event trigger time are stored.



In this window:

- You can view the desired events by selecting [**Event**].
- You can view the desired events according to the level by selecting [**Level**].

After selecting the desired event, you can select [**Details**] to access the following window. In this window, the waveform area displays the waveforms related to the event, and the parameter area displays the parameter values happened at the event trigger time.




- 1. Waveform area
- 2. Parameter area

In this window:

- You can select or to navigate through the waveforms.
- You can select or beside the [**Event**] button to switch between events.
- You can set the desired [**Gain**] for ECG waveform.
- You can set the desired [**Sweep**].
- By selecting the [**Record**] button, you can print out the currently displayed alarm events by the recorder.
- By selecting the [**Print**] button, you can print out the currently displayed alarm events by the printer.
- By selecting the [**Events List**] button, you can view the events list.

19.5 Reviewing Waveforms

In the [Review] menu, select [Full Disclosure] to access the following window.



In this review window:

- To review full-disclosure waveforms, you need to save waveforms first. Select [**Save Waves** >>] and then select the parameters whose waveforms you want to view. To save full-disclosure waveform, your monitor must be equipped with a storage card.
- To view the waveforms, you can either:
 - ◆ Select or beside the [**Scroll**] button to move the cursor one step left or right to navigate through the waveforms, or
 - Select or right to navigate through the waveforms.

A time indicating your current position is displayed at the top of the waveform area. Numeric measurement values corresponding to the cursor location are displayed in the parameter area, and change as the cursor is moved.

- You can change the ECG wave gain by selecting [**Gain**] and then selecting the appropriate setting.
- You can change the waveform sweep speed by selecting [**Sweep**] and then selecting the appropriate setting.

By selecting the [Record] button, you can print out the first three waveforms and measurement numerics by the recorder.
By selecting or beside the [Event] button, you can position the cursor between events.

20 Calculations

20.1 Introduction

The calculation feature is available with your patient monitor. The calculated values, which are not directly measured, are computed based on the values you provide.

Your can perform the following calculations:

- Dose calculations
- Oxygenation calculations
- Ventilation calculations
- Hemodynamic calculations
- Renal calculations

To perform a calculation, select [**Main Menu**] → [**Calc >>**], or the [**Calculations**] QuickKey and then select the calculation you want to perform.



∆ WARNING

 After the calculation is finished, verify the entered values are correct and the calculated values are appropriate. We assume no responsibility for any consequences caused by wrong entries and improper operations.

NOTE

The calculation feature is independent of other monitoring functions and can be therefore used for
patients being monitored by other monitors. Any operation in a calculation window does not affect the
patient monitoring by the local patient monitor.

20.2 Dose Calculations

20.2.1 Performing Calculations

To perform a dose calculation:

- Select [Main Menu]→[Calculations >>]→[Dose >>], or select [Calculations] QuickKey→[Dose >>].
- 2. Select, in turn, [**Patient Cat.**] and [**Drug Name**] and then select the appropriate settings. The dose calculation program has a library of commonly used drugs, of which Drug A through Drug E are for those not specified in this library.
 - ♦ Drug A, B, C, D, E
 - ♦ Aminophylline
 - ♦ Dobutamine
 - ◆ Dopamine
 - ◆ Epinephrine
 - ♦ Heparin

- ◆ Isuprel
- ◆ Lidocaine
- ♦ Nipride
- ♦ Nltroglycerin
- ◆ Pitocin
- 3. The system gives a set of default values when the above steps are finished. However, these values cannot be used as the calculated values. The user must enter values following the doctor's instructions, and then the calculated values can only be used
- 4. Enter the patient's weight.
- 5. Enter other values.
- 6. Verify if the calculated values are correct.

20.2.2 Selecting the Proper Drug Unit

Each drug has its fixed unit or unit series. Among a unit series, one unit may change to another automatically depending on the entered value.

The units for each drug are as follows:

- Drug A, B, C, Aminophylline, Dobutamine, Dopamine, Epinephrine, Isuprel, Lidocaine, Nipride and Nitroglycerin use the unit series: g, mg and mcg.
- Drug D, Heparin and Pitocin use the unit series: Unit, KU (kilo units) and MU (million units).
- Drug E uses the unit: mEq (milli-equivalents).

You must select the proper drug name (A, B, C, D or E) according to the units when you define a drug not listed in this library.

NOTE

• For neonate patients, [Drip Rate] and [Drop Size] are disabled.

20.2.3 Titration Table

To open the titration table, select [**Titration Table >>**] in the [**Dose Calculation**] window after the dose calculation is finished.

In the titration table, when you change:

- [Reference]
- [Interval]
- [Dose Type]

The titrated values change accordingly.

You can also:

- Select or , or or beside the vertical scrollbar to view more values.
- Select [**Record**] to print out the currently displayed titrated values by the recorder.

20.3 Oxygenation Calculations

20.3.1 Performing Calculations

To perform an oxygenation calculation:

- Select [Main Menu]→[Calculations >>]→[Oxygenation >>], or select [Calculations]
 QuickKey→[Oxygenation >>].
- 2. Enter values for calculation.
- 3. Select the [**Calculate**] button. The system performs a calculation per the current settings and displays the calculated values.
 - If a calculated value is outside the range, its background will highlight in yellow. You can select [**Range**] to view its normal range in the unit field.
 - ◆ Invalid values are displayed as [---].

In the [Oxygenation Calculation] window, you can:

- Change the pressure unit, Hb unit and oxygen content unit by selecting [Press. Unit], [Hb Unit] and [OxyCont Unit] and then selecting the appropriate settings. The changes take effect automatically.
- Trigger a recording by selecting the [**Record**] button. The currently displayed oxygenation calculations are printed out by the recorder.
- Review the previously performed calculations by selecting [Review].

20.3.2 Entered Parameters

Abbreviation	Unit	Full spelling
C.O.	L/min	cardiac output
FiO ₂	%	percentage fraction of inspired oxygen
PaO ₂	mmHg	partial pressure of oxygen in the arteries
PaCO ₂	mmHg	partial pressure of carbon dioxide in the arteries
SaO ₂	%	arterial oxygen saturation
PvO ₂	mmHg	partial pressure of oxygen in venous blood
SvO ₂	%	venous oxygen saturation
Hb	g/L	hemoglobin
CaO ₂	ml/L	arterial oxygen content
CvO ₂	ml/L	venous oxygen content
VO ₂	ml/min	oxygen consumption
RQ	None	respiratory quotient
ATMP	mmHg	atmospheric pressure
Height	cm	height
Weight	kg	weight

20.3.3 Calculated Parameters

Abbreviation	Unit	Full spelling
BSA	m²	body surface area
VO₂ calc	ml/min	oxygen consumption
C(a-v)O ₂	ml/L	arteriovenous oxygen content difference
O ₂ ER	%	oxygen extraction ratio
DO ₂	ml/min	oxygen transport
PAO ₂	mmHg	partial pressure of oxygen in the alveoli
AaDO ₂	mmHg	alveolar-arterial oxygen difference
CcO ₂	ml/L	capillary oxygen content
Qs/Qt	%	venous admixture
C.O. calc	L/min	calculated cardiac output

20.4 Ventilation Calculations

20.4.1 Performing Calculations

To perform a ventilation calculation:

- Select [Main Menu]→[Calculations >>]→[Ventilation >>], or select [Calculations]
 QuickKey→[Ventilation >>].
- 2. Enter values for calculation. If the patient monitor is connected to an anesthesia machine or a ventilator, the system automatically loads the supported parameter values to the [**Ventilation Calculation**] window.
- 3. Select the [**Calculate**] button. The system performs a calculation per the current settings and displays the calculated values.

- If a calculated value is outside the range, its background will highlight in yellow. You can select [**Range**] to view its normal range in the unit field.
- ♦ Invalid values are displayed as [---].

In the [Ventilation Calculation] window, you can:

- Change the pressure unit by selecting [Press. Unit] and then selecting the appropriate setting.
 Corresponding pressure values shall convert and update automatically.
- Trigger a recording by selecting the [**Record**] button. The currently displayed ventilation calculations are printed out by the recorder.
- Review the previously performed calculations by selecting [**Review**].

20.4.2 Entered Parameters

Abbreviation	Unit	Full spelling	
FiO ₂	%	percentage fraction of inspired oxygen	
RR	rpm	respiration rate	
PeCO ₂	mmHg	partial pressure of mixed expiratory CO2	
PaCO ₂	mmHg	partial pressure of carbon dioxide in the arteries	
PaO ₂	mmHg	partial pressure of oxygen in the arteries	
TV	ml	tidal volume	
RQ	None	respiratory quotient	
ATMP	mmHg	atmospheric pressure	

20.4.3 Calculated Parameters

Abbreviation	Unit	Full spelling	
PAO ₂	mmHg	partial pressure of oxygen in the alveoli	
AaDO ₂	mmHg	alveolar-arterial oxygen difference	
Pa/FiO ₂	mmHg	oxygenation ratio	
a/AO ₂	%	arterial to alveolar oxygen ratio	
MV	L/min	minute volume	
Vd	ml	volume of physiological dead space	
Vd/Vt	%	physiologic dead space in percent of tidal volume	
VA	L/min	alveolar volume	

20.5 Hemodynamic Calculations

20.5.1 Performing Calculations

To perform a hemodynamic calculation:

- Select [Main Menu]→[Calculations >>]→[Hemodynamic >>], or select [Calculations]
 QuickKey→[Hemodynamic >>].
- 2. Enter values for calculation.
 - ◆ For a patient who is being monitored, [HR], [Art mean], [PA mean] and [CVP] are automatically taken from the currently measured values. If you just have performed C.O. measurements, [C.O.] is the average of multiple thermodilution measurements. [Height] and [Weight] are the patient's height and weight you have entered. If the monitor does not provide these values, their fields appear blank.
 - For a patient who is not being monitored, confirm the values you have entered.
- 3. Select the [**Calculate**] button. The system performs a calculation per the current settings and displays the calculated values.
 - If a calculated value is outside the range, its background will highlight in yellow. You can select [**Range**] to view its normal range in the unit field.
 - ◆ Invalid values are displayed as [---].

In the [Hemodynamic Calculation] window, you can:

- Trigger a recording by selecting the [**Record**] button. The currently displayed renal calculations are printed out by the recorder.
- Review the previously performed calculations by selecting [**Review**].

20.5.2 Entered Parameters

Abbreviation	Unit	Full spelling
C.O.	L/min	cardiac output
HR	bpm	heart rate
PAWP	mmHg	pulmonary artery wedge pressure
Art Mean	mmHg	artery mean pressure
PA Mean	mmHg	pulmonary artery mean pressure
CVP	mmHg	central venous pressure
EDV	ml	end-diastolic volume
Height	cm	height
Weight	kg	weight

20.5.3 Calculated Parameters

Abbreviation	Unit	Full spelling
BSA	m ²	body surface area
C.I.	L/min/m ²	cardiac index
SV	ml	stroke volume

Abbreviation	Unit	Full spelling
SI	ml/m²	stroke index
SVR	DS/cm ⁵	systemic vascular resistance
SVRI	DS·m²/cm⁵	systemic vascular resistance index
PVR	DS/cm ⁵	pulmonary vascular resistance
PVRI	DS·m²/cm⁵	pulmonary vascular resistance index
LCW	kg⋅m	left cardiac work
LCWI	kg·m/m²	left cardiac work index
LVSW	g⋅m	left ventricular stroke work
LVSWI	g·m/m²	left ventricular stroke work index
RCW	kg⋅m	right cardiac work
RCWI	kg·m/m²	right cardiac work index
RVSW	g⋅m	right ventricular stroke work
RVSWI	g·m/m²	right ventricular stroke work index
EF	%	ejection fraction

20.6 Renal Calculations

20.6.1 Performing Calculations

To perform a renal calculation:

- Selecting [Main Menu]→[Calculations >>]→[Renal >>], or select [Calculations] QuickKey→[Renal >>].
- 2. Enter values for calculation.
- 3. Select the [**Calculate**] button. The system performs a calculation per the current settings and displays the calculated values.
 - If a calculated value is outside the range, its background will highlight in yellow. You can select [**Range**] to view its normal range in the unit field.
 - ◆ Invalid values are displayed as [---].

In the [**Renal Calculation**] window, you can:

- Trigger a recording by selecting the [**Record**] button. The currently displayed renal calculations are printed out by the recorder.
- Review the previously performed calculations by selecting [**Review**].

20.6.2 Entered Parameters

Abbreviation	Unit	Full spelling
URK	mmol/L	urine pstassium
URNa	mmol/L	urinary sodium
Urine	ml/24h	urine
Posm	mOsm/ kgH₂O	plasm osmolality
Uosm	mOsm/ kgH₂O	urine osmolality
SerNa	mmol/L	serum sodium
Cr	μmol/L	creatinine
UCr	μmol/L	urine creatinine
BUN	mmol/L	blood urea nitrogen
Height	cm	height
Weight	kg	weight

20.6.3 Calculated Parameters

Abbreviation	Unit	Full spelling
URNaEx	mmol/24h	urine sodium excretion
URKEx	mmol/24h	urine potassium excretion
Na/K	%	sodium potassium ratio
CNa	ml/24h	clearance of sodium
Clcr	ml/min	creatinine clearance rate
FENa	%	fractional excretion of sodium
Cosm	ml/min	osmolar clearance
CH₂O	ml/h	free water clearance
U/P osm	None	urine to plasma osmolality ratio
BUN/Cr	None*	blood urea nitrogen creatinine ratio
U/Cr	None	urine-serum creatinine ratio

^{*:} BUN/Cr is a ratio under the unit of mol.

20.7 Understanding the Review Window


With the review feature, you can review oxygenation, ventilation, hemodynamic and renal calculations. The review window for each calculation is similar. Take the hemodynamic calculations review window for example, you can access it by selecting [**Review**] in the [**Hemodynamic Calculation**] window.

In this review window:

- You can select ◀, ▶ ◀◀ or ▶▶ to view more values.
- The values that exceed the range are displayed in yellow background. The [**Unit**] field displays parameter units. If some parameter values are outside of their normal ranges, you can view their normal range in the [**Unit**] field by selecting [**Range**].
- You can review an individual calculation by selecting its corresponding column and then selecting [Original Calc]. You can record the currently displayed calculations or perform another calculation is this window.

21.1 Using a Recorder

The thermal recorder records patient information, measurement numerics, up to three waveforms, etc.



- 1. Start/Stop key: press to start a recording or stop the current recording.
- 2. Indicator
 - On: when the recorder works correctly.
 - Off: when the monitor is switched off.
 - Flashes: if an error occurred to the recorder, e.g., the recorder runs out of paper.
- 3. Paper outlet
- 4. Latch
- 5. Recorder door

21.2 Overview of Recording Types

By the way recordings are triggered, the recordings can be classified into the following categories:

- Manually-triggered realtime recordings.
- Timed recordings.
- Alarm recordings triggered by an alarm limit violation or an arrhythmia event.
- Manually-triggered, task-related recordings.

The task-related recordings include:

- Frozen wave recording
- Graphic trends recording
- Tabular trends recording
- Events recording:parameter alarm recording, arrh. alarm recording, manual events recording
- Wave review recording
- Titration table recording
- Hemodynamic calculations recording
- Oxygenation calculations recording
- Ventilation calculations recording
- Renal calculations recording
- oxyCRG recording
- C.O. curve recording
- PAWP recording
- Respiratory loops recording
- Monitor information recording

NOTE

- For details about alarm recording, refer to the chapter 7 Alarms.
- For details about task-related recordings, refer to respective sections of this manual.

21.3 Starting and Stopping Recordings

To manually start a recording, you can either:

- Select the 🖸 hardkey on the front of either the patient monitor or the recorder module, or
- Select the [**Record**] button from the current menu or window.

Automatic recordings will be triggered in the following conditions:

- Timed recordings will start automatically at preset intervals.
- If both [**Alarm**] and [**Alm Rec**] for a measurement are set on, an alarm recording will be triggered automatically as alarms occur.

To manually stop a recording, you can either:

- Select [Clear All Tasks] in the [Record Setup] menu.

Recordings stop automatically when:

- The runtime is over.
- The recorder runs out of paper.
- When the recorder has an alarm condition.

When a recording is stopped, the following markers will be added:

- Automatically stopped recording: print two columns of '*' at the end of the report.
- Manually or abnormally stopped recording: print one column of '*' at the end of the report.

21.4 Setting up the Recorder

21.4.1 Accessing the Record Setup Menu

By selecting [Main Menu]→[Record Setup >>], you can access the [Record Setup] menu.

21.4.2 Selecting Waveforms for Recording

The recorder can record up to 3 waveforms at a time. You can select, in turn, [**Waveform 1**], [**Waveform 2**] and [**Waveform 3**] in the [**Record Setup**] menu, and then select the waveforms you want. You can also turn off a waveform recording by selecting [**Off**]. These settings are intended for realtime and scheduled recordings.

21.4.3 Setting the Realtime Recording Length

After starting a realtime recording, the recording time depends on your monitor's settings. In the [**Record Setup**] menu, select [**Length**] and toggle between [**8 s**] and [**Continuous**].

- [8 s]: record 4-second waveforms respectively before and after current moment.
- [Continuous]: record the waveforms from the current moment until stopped manually.

21.4.4 Setting the Interval between Timed Recordings

Timed recordings start automatically at preset intervals. Each recording lasts 8 seconds. To set the interval between timed recordings: in the [**Record Setup**] menu, select [**Interval**] and then select the appropriate setting.

21.4.5 Changing the Recording Speed


In the [Record Setup] menu, select [Paper Speed] and toggle between [25 mm/s] and [50 mm/s]. This setting is for all recordings containing waveforms.

21.4.6 Clearing Recording Tasks

In the [**Record Setup**] menu, select [**Clear All Tasks**]. All queued recording tasks are cleared and the current recording is stopped.

21.5 Loading Paper

- 1. Use the latch at the upper right of the recorder door to pull the door open.
- 2. Insert a new roll into the compartment as shown below.
- 3. Close the recorder door.
- 4. Check if paper is loaded correctly and the paper end is feeding from the top.



ACAUTION

- Use only specified thermal paper. Otherwise, it may cause damage to the recorder's printhead, the recorder may be unable to print, or poor print quality may result.
- Never pull the recorder paper with force when a recording is in process. Otherwise, it may cause damage to the recorder.
- Do not leave the recorder door open unless you reload paper or remove troubles.

21.6 Removing Paper Jam

If the recorder works incorrectly or produces unusual sounds, check if there is a paper jam first. If a paper jam is detected, follow this procedure to remove it:

- 1. Open the recorder door.
- 2. Take out the paper and tear off the draped part.
- 3. Reload the paper and close the recorder door.

21.7 Cleaning the Recorder Printhead

If the recorder has been used for a long time, deposits of paper debris may collect on the printhead compromising the print quality and shortening the lifetime of the roller. Follow this procedure to clean the printhead:

- 1. Take measures against the static electricity such as Disposable Wrist Strap for the work.
- 2. Open the recorder door and take out the paper.
- ${\it 3.} \quad {\it Gently wipe around the printhead using cotton swabs dampened with alcohol.}$
- 4. After the alcohol has completely been dried, reload the paper and close the recorder door.

ACAUTION

- Do not use anything that may destroy the thermal element.
- Do not add unnecessary force to the thermal head.

22 Printing

22.1 Printer

The monitor can output patient reports via a connected printer. So far, the monitor supports the following printer:

- HP LaserJet 1505n
- HP LaserJet P2035n
- HP LaserJet P4015n
- HP LaserJet 1606dn

The specifications of the reports the monitor prints are:

- Paper: A4, Letter
- Resolution: 300 dpi
- Print on One/Both Sides: printing on one and both sides are supported if the printer supports

For more details about the printer, see the document accompanying the printer. With the upgrading of products, the monitor will support more printers and no prior notice will be given. If you have any doubt about the printer you have purchased, contact our company.

22.2 Connecting a printer

To print the reports or the trend data of a patient, you can choose either:

- the local printer
 - Connect the printer and the patient monitor directly with a network cable, and then start printing what you want, or
- the Central Monitoring System

If your monitor is connected to a central monitoring system, it is recommended to use the central monitoring system for printing.

22.3 Setting Up the Printer

To set the printer's properties, select [Main Menu]→[Print Setup >>]→[Printer Setup >>]. In the [Printer Setup] menu, you can:

Select a connected printer

Select [**Printer**] and then select a connected printer as the monitor's printer.

Search for a printer

If your selected printer is not in the list or a new printer is added into the network, you can select the [**Search Printer**] to re-search for all printers in the network.

Set up the paper

Select [Paper Size] and toggle between [A4] and [Letter].

Print on both sides

By default, the monitor prints out patient reports on one side. If you set [**Print On Both Sides**] to [**On**], the monitor will print out patient reports on both sides.

22.4 Starting Report Printouts

Reports	Contents	Procedures
ECG reports	ECG waveforms and relevant	Select [Main Menu]→[Print Setup >>]→[ECG Reports
	parameter values	>>]→[Print]
Tabular trends	Depend on the selected parameter	Select [Main Menu]→[Print Setup >>]→[Tabular
	group, resolution and time period	Trends Reports >>]→[Print], or select [Main
		Menu]→[Review >>]→[Tabular
		Trends]→[Print]→[Print]
Graphic trends	Depend on the selected parameter	Select [Main Menu]→[Print Setup >>]→[Graphic
	group, resolution and time period	Trends Reports >>]→[Print], or select [Main
		Menu]→[Review >>]→[Graphic
		Trends]→[Print]→[Print]
Arrh. alarm	ECG waveforms and relevant	Select [Print] in [Arrh. Events]
review	parameter values	
Parameter alarm	Depend on the selected alarms	Select [Main Menu]→[Review >>]→[Alarms]→[Print]
review		
Realtime waves	Depend on the selected waveforms	Select [Main Menu]→[Print Setup >>]→[Realtime
		Reports >>]→[Print]

22.5 Stopping Reports Printouts

To stop report printouts, select [Main Menu]→[Print Setup >>]→[Stop All Reports].

22.6 Setting Up Reports

22.6.1 Setting Up ECG Reports

You can print out ECG reports only under full-screen, half-screen or 12-lead monitoring screen. To set up ECG reports, select [Main Menu]→[Print Setup >>]→[ECG Reports >>].

- [Amplitude]: set the amplitude of the ECG waveforms.
- [Sweep]: set the wave print speed.
- [Auto Interval]: If [Auto Interval] is set to [On], the system will automatically adjust the space between waveforms to avoid overlapping.
- [**Gridlines**]: choose whether to show gridlines.

22.6.2 Setting Up Tabular Trends Reports

To set up tabular trends reports, select [Main Menu]→[Print Setup >>]→[Tabular Trends Reports >>].

- Start time: You can set a time period whose trend data will be printed out by setting [**From**] and [**Back**]. For example, if you set [**From**] as 2007-4-2 10:00:00 and [**Back**] as [**2 h**], the outputted data will be from 2007-4-2 08:00:00 to 2007-4-2 10:00:00. In addition, the [**Back**] can be set to either:
 - ◆ [Auto]: If [Report Layout] is set to [Time Oriented], the report will be printed by time. If [Report Layout] is set to [Parameter Oriented], the report will be printed by parameters.
 - ♦ [All]: If you select [All], all trend data will be printed out. In this case, it is no need to set [From].
- [Interval]: choose the resolution of the tabular trends printed on the report.
- [Report Layout]: If you select [Time Oriented], the report will be printed by time. If you select [Parameter Oriented], the report will be printed by parameters.
- [Select Parameter >>]: from the popup menu, you can:
 - ♦ [Currently Displayed Trended Parameters]: print the parameter trend data selected from the [Tabular Trends].
 - ◆ [**Standard Parameter Group**]: select the standard parameter group for printing.
 - [**Custom**]: You can define a parameter group for printing from the parameters displayed in the low part of the menu.

22.6.3 Setting Up Graphic Trends Reports

To set up graphic trends reports, select [Main Menu]→[Print Setup >>]→[Graphic Trends Reports >>]. As setting up graphic trends reports is similar with tabular trends reports, you can refer to the Setting Up Tabular Trend Reports section for details.

22.6.4 Setting Up Realtime Reports

To set up realtime reports, select [Main Menu]→[Print Setup >>]→[Realtime Reports >>].

- [Sweep]: set the wave print speed.
- [**Select Wave** >>]: from the popup menu, you can:
 - ◆ [Current]: select the currently displayed waves for printing.
 - ◆ [**Select Wave**]: select the desired waves for printing.

22.7 End Case Reports

ECG reports, tabular trends reports, graphic trends reports, NIBP review reports and realtime reports can be set as end case reports. When you discharge a patient, the system will automatically print out all contents that are set as end case reports.

For example, to set ECG report as end case report:

- 1. select [Main Menu]→[Print Setup >>]→[ECG Report >>].
- 2. select [End Case Report]→[Set as End Case Report] and then select [Ok] from the popup dialog box.
- 3. set as described in the **22.6.1 Setting Up ECG Reports**.

22.8 Printer Statuses

22.8.1 Printer Out of Paper

When the printer runs out of paper, the print request will not be responded. If there are too many print jobs that are not responded, a printer error may occur. In these cases, you need to install paper and then re-send the print request. Restart the printer if necessary.

Therefore, you'd better ensure that there is enough paper in the printer before sending a print request.

22.8.2 Printer Status Messages

If the monitor prompts that selected printer is not available, check that the printer is switched on, correctly connected, and installed with paper.

23 Other Functions

23.1 Analog Output

The patient monitor provides analog output signals to accessory equipment via the multifunctional connector on the rear of the monitor. To obtain analog output signals, connect the accessory equipment such as an oscillograph, etc. to the monitor.

NOTE

 The analog output feature is seldom applied in clinical applications. You can contact your service personnel for more details.

23.2 Transferring Data

You can transfer the patient data saved in the monitor to a PC via a crossover network cable or within a LAN for data management, review or print.

23.2.1 Data Export System

You must install the data export system on the intended PC before performing the data transfer operation. Refer to the document accompanying the installation CD-ROM for installation instructions.

The data transfer feature supports patient management, data review, data format conversion, print, etc. in addition to data transfer. Refer to the help file of the system software for more details.

23.2.2 Transferring Data by Different Means

NOTE

 Never enter the data transfer mode when the patient monitor is in normal operation or performs monitoring. You must re-start the patient monitor to exit the data transfer mode.

Transfer data via a crossover network cable

Before transferring data using a crossover network cable, do as follows:

- 1. Connect one end of the crossover network cable to the patient monitor and the other end to the PC.
- 2. Set the IP address of the PC. This IP address must be in the same network segment with that of the patient monitor.
- 3. Make sure that the data export system is active on the PC.

Then, follow this procedure to transfer data:

- Select [Main Menu]→[Patient Data >>]→[Transfer Data].
- 2. Select [**Yes**] from the popup message box.
- 3. Input the IP address already set on the PC.
- 4. Select [**Start**] to start transferring data.

Transfer data within a LAN

Before transferring data within a LAN, do as follows:

- 1. Connect the patient monitor and the intended PC into the same LAN and acquire the PC's IP address.
- 2. Make sure that the data export system is active on the PC.

Follow the same procedure as via a crossover network cable to transfer data.

23.3 Nurse Call

The patient monitor also provides nurse call signals to a nurse call system connected to the monitor via the multifunctional connector. To obtain nurse call signals, connect a nurse call system to the monitor and then follow this procedure:

- 1. Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password.
- 2. Select [Others >>] to access the [Others] menu.
- 3. Select [**Nurse Call Setup >>**] to change the nurse call settings as follows:
- Select [Signal Type] and toggle between [Pulse] and [Continuous].
 - [**Pulse**]: the nurse call signals are pulse signals and each pulse lasts 1 second. When multiple alarms occur simultaneously, only one pulse signal is outputted. If an alarm occurs but the previous one is not cleared yet, a new pulse signal will also be outputted.
 - [Continuous]: the nurse call signal lasts until the alarm ends, i.e. the duration of a nurse call signal equals to that of the alarm condition.
- Select [Contact Type] and toggle between [Normally Open] and [Normally Closed].
 - [Normally Open]: select if your hospital's nurse call relay contact is normally open.
 - [Normally Closed]: select if your hospital's nurse call relay contact is normally closed.
- Select [**Alm Lev**] and set the alarm level for nurse call-triggering alarms.
- Select [Alarm Cat.] and then select the category to which the nurse call-triggering alarms belong.

Alarm conditions are indicated to nurses only when:

- The nurse call system is enabled,
- An alarm that meets your preset requirements occurs, and
- The monitor is not in the alarm paused or silence status.

WARNING

• Do not rely exclusively on the nurse call system for alarm notification. Remember that the most reliable alarm notification combines audible and visual alarm indications with the patient's clinical condition.

NOTE

 If no setting is selected from [Alm Lev] or [Alarm Cat.], no nurse call signal will be triggered whatever alarms occur.

23.4 Network Connection

23.4.1 Setting the Network Type

The patient monitor supports both wired and wireless network. To set the network type, you can select [Main Menu]
→[Maintenance>>]→[User Maintenance>>]→enter the required password→[Network Setup >>].

23.4.2 Wireless Network

The patient monitors can be connected to a wireless network via a built-in Wi-Fi module. To set the wireless network:

- Select [Main Menu]→[Maintenance>>]→[User Maintenance>>]→enter the required password→
 [Network Setup >>].
- 2. Set the [Network Type] to [WLAN].
- 3. Select [WLAN Setup >>] to access the [WLAN Setup] menu.
- 4. Configure the [Network Name (SSID)], and [Password].
- 5. Click [**OK**] to confirm the setting.

To test the availability of the wireless network, follow this procedure:

- 1. Select [**WLAN Test >>**] in the [**WLAN Setup**] menu.
- 2. Enter the [**IP Address**] of wireless AP in the [**WLAN Test >>**] menu.
- 3. Click [Connection Test].

The Wi-Fi device used in the monitor is in compliance with IEEE 802.11b/g/n.

23.4.3 Setting the IP Address, Subnet Mask and Gateway

In the [**Network Setup**] menu, you can set IP address, subnet mask and gateway. You should not change the patient monitor's IP address randomly. If you want to know details about IP address setup, contact the technical personnel in charge of the CMS.

NOTE


- The design, installation, restruction and maintenance of the wireless network's distribution shall be performed by authorized service personnel of our company.
- The existence of obstacles (such as wall) will exert impact on data transferring or even cause network interruption.
- The Central Monitoring System is capable of connecting up to 32 bedside monitors via the wireless network.

24 Batteries

24.1 Overview

This monitor is designed to operate on rechargeable Lithium-ion battery power during intra-hospital patient transfer or whenever the power supply is interrupted. The battery is charged automatically when the monitor is connected to AC/DC power, no matter the monitor is powered on or not. Whenever the AC/DC power is interrupted during patient monitoring, the patient monitor automatically runs power from the internal battery. iPM 12 patieng monitor is available for up to two batteries. iPM 10 or iPM 8 patient monitor is available for only one battery.

On-screen battery symbols indicate the battery status as follows:



The capacity of the internal battery is limited. If the battery charge is too low, a technical alarm will be triggered and the message [Low Battery] or [Battery Depleted] displayed. At this moment, apply AC/DC power to the patient monitor. Otherwise, the patient monitor will power off automatically before the battery is completely depleted.


24.2 Replacing a Battery


When the iPM 12 patient monitor uses two battery packs, one battery pack can be easily exchanged while the patient monitor operates from the other. If the iPM 12 patient monitor uses one battery pack, you should insert a new battery pack before the old one depletes.

When the iPM 10 or iPM 8 patient monitor operates on battery power, make sure the patient monitor is power off before replacing a battery.

To replace a battery, follow this procedure:

1. Open the battery door.





iPM 12

iPM 10/iPM 8

- 2. Push aside the latch fixing the battery to be replaced and remove the battery.
- 3. Insert a battery into the slot with its contact point inward.
- 4. Close the battery door.

24.3 Battery Guidelines

Life expectancy of a battery depends on how frequent and how long it is used. For a properly maintained and stored lithium ion battery, its life expectancy is about 3 years. For more aggressive use models, life expectancy can be less. We recommend replacing lithium ion batteries every 3 years.

To get the most out of the battery, observe the following guidelines:

- The battery performance test must be performed every two years, before monitor repairs, or whenever the battery is suspected as being the source of the problems.
- Condition a battery once when it is used or stored for 3 months, or when its operating time becomes noticeably shorter.
- Take out the battery before the monitor is transported or will not be used for more than 3 months.
- Remove the battery from the monitor if it is not being used regularly. (Leaving the battery in a monitor that is not in regular use will shorten the life of the battery).
- The shelf life of a Lithium Ion battery is about 6 months when the battery is stored with the battery power being 50% of the total power. In 6 months the battery power must be depleted before the Lithium Ion battery is fully charged. Then run the monitor on this fully charged battery .When its battery power becomes 50% of the total power, take out the battery from the monitor and store it.

MWARNING

- Keep the battery out of the reach of children.
- Use only the battery specified by the manufacturer.
- If the battery shows signs of damage or signs of leakage, replace it immediately. Do not use a faulty battery in the monitor.

24.4 Battery Maintenance

Conditioning a Battery

A battery should be conditioned before it is used for the first time. A battery conditioning cycle is one uninterrupted charge of the battery, followed by an uninterrupted battery discharge and charge. Batteries should be conditioned regularly to maintain their useful life.

NOTE

• The actual battery capacity will decrease over time with use of batteries. When a monitor operates on batteries that have been used before, the full capacity battery symbol does not indicate the capacity and operating time of this battery can still fulfill battery specifications in the operator's manual. When conditioning a battery, please replace the battery if its operating time is significantly lower than the specified time.

To condition a battery, follow this procedure:

- 1. Disconnect the monitor from the patient and stop all monitoring or measuring.
- 2. Insert the battery in need of conditioning in the battery slot of the monitor.
- 3. Apply AC/DC power to the monitor and allow the battery to charge uninterrupted for 10 hours.
- 4. Remove AC/DC power and allow the monitor to run from the battery until it shuts off.
- 5. Apply AC/DC power again to the monitor and allow the battery to charge uninterrupted for 10 hours.
- 6. This battery is now conditioned and the monitor can be returned to service.

Checking a Battery

The battery performance test must be performed every two years, before monitor repairs, or whenever the battery is suspected as being the source of the problems. The performance of a rechargeable battery may deteriorate over time. To check the performance of a battery, follow this procedure:

- 1. Disconnect the monitor from the patient and stop all monitoring or measuring.
- 2. Apply AC/DC power to the monitor and allow the battery to charge uninterrupted for 10 hours.
- 3. Remove AC/DC power and allow the monitor to run from the battery until it shuts off.
- 4. The operating time of battery reflects its performance directly.

Please replace the battery or contact with the maintenance personnel if its operating time is significantly lower than the specified time.

NOTE

- The battery might be damaged or malfunctioned if its operating time is too short after being fully charged.
 The operating time depends on the configuration and operation. For example, measuring NIBP more frequently will also shorten the operating time.
- When a battery has visual signs of damage, or no longer holds a charge, it should be replaced. Remove the old battery from the monitor and recycle it properly.

24.5 Battery Recycling

When a battery has visual signs of damage, or no longer holds a charge, it should be replaced. Remove the old battery from the monitor and recycle it properly. To dispose of the batteries, follow local laws for proper disposal.



Do not disassemble batteries, or dispose of them in fire, or cause them to short circuit. They may ignite,
 explode, leak or heat up, causing personal injury.

25 Care and Cleaning

Use only the substances approved by us and methods listed in this chapter to clean or disinfect your equipment.

Warranty does not cover damage caused by unapproved substances or methods.

We make no claims regarding the efficacy of the listed chemicals or methods as a means for controlling infection. For the method to control infection, consult your hospital's Infection Control Officer or Epidemiologist.

25.1 General Points

Keep you equipment and accessories free of dust and dirt. To avoid damage to the equipment, follow these rules:

- Always dilute according the manufacturer's instructions or use lowest possible concentration.
- Do not immerse part of the equipment into liquid.
- Do not pour liquid onto the equipment or accessories.
- Do not allow liquid to enter the case.
- Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).



WARNING

• Be sure to shut down the system and disconnect all power cables from the outlets before cleaning the equipment.



∕!∆ CAUTION

• If you spill liquid on the equipment or accessories, contact us or your service personnel.

NOTE

• To clean or disinfect reusable accessories, refer to the instructions delivered with the accessories.

25.2 Cleaning

Your equipment should be cleaned on a regular basis. If there is heavy pollution or lots of dust and sand in your place, the equipment should be cleaned more frequently. Before cleaning the equipment, consult your hospital's regulations for cleaning the equipment.

Recommended cleaning agents are:

- Sodium hypochlorite bleach (diluted)
- Hydrogen peroxide (3%)
- Ethanol (70%)
- Isopropanol (70%)

To clean your equipment, follow these rules:

- 1. Shut down the patient monitor and disconnect it from the power line.
- 2. Clean the display screen using a soft, clean cloth dampened with a glass cleaner.
- 3. Clean the exterior surface of the equipment using a soft cloth dampened with the cleaner.
- 4. Wipe off all the cleaning solution with a dry cloth after cleaning if necessary.
- 5. Dry your equipment in a ventilated, cool place.

25.3 Disinfecting

Disinfection may cause damage to the equipment and is therefore not recommended for this patient monitor unless otherwise indicated in your hospital's servicing schedule. Cleaning equipment before disinfecting is recommended.

The recommended disinfectants include: ethanol 70%, isopropanol 70%, Perform® classic concentrate OXY (KHSO₄ solution).



Never use EtO or formaldehyde for disinfection.

26 Maintenance

⚠ WARNING

- Failure on the part of the responsible individual hospital or institution employing the use of this equipment to implement a satisfactory maintenance schedule may cause undue equipment failure and possible health hazards.
- The safety checks or maintenance involving any disassembly of the equipment should be performed by professional servicing personnel. Otherwise, undue equipment failure and possible health hazards could
- If you discover a problem with any of the equipment, contact your service personnel or us.

26.1 Safety Checks

Before every use, after your patient monitor has been used for 6 to 12 months, or whenever your patient monitor is repaired or upgraded, a thorough inspection should be performed by qualified service personnel to ensure the reliability.

Follow these guidelines when inspecting the equipment:

- Make sure that the environment and power supply meet the requirements.
- Inspect the equipment and its accessories for mechanical damage.
- Inspect all power cords for damage, and make sure that their insulation is in good condition.
- Make sure that only specified accessories are applied.
- Inspect if the alarm system functions correctly.
- Make sure that the recorder functions correctly and the recorder paper meets the requirements.
- Make sure that the batteries meet the performance requirements.
- Make sure that the patient monitor is in good working condition.
- Make sure that the grounding resistance and leakage current meet the requirement.

In case of any damage or abnormity, do not use the patient monitor. Contact the hospital's biomedical engineers or your service personnel immediately.

26.2 Checking Monitor and Module Information

To view the information about system start time, selftest, etc., select [Main Menu]→[Maintenance >>]→[Monitor Information >>]. You can print out the information for the convenience of troubleshooting. The information will not be saved during shut down.

You can also view the information about the monitor configuration and system software version by selecting [Main Menu]→[Maintenance >>]→[Software Version >>].

26.3 Calibrating ECG

The ECG signal may be inaccurate due to hardware or software problems. As a result, the ECG wave amplitude becomes greater or smaller. In that case, you need to calibrate the ECG module.

- Select the ECG parameter window or waveform area→[Filter]→[Diagnostic].
- Select [Main Menu]→[Maintenance >>]→[Calibrate ECG]. A square wave appears on the screen and the message [ECG Calibrating] is displayed.
- 3. Compare the amplitude of the square wave with the wave scale. The difference should be within 5%.
- 4. After the calibration is completed, select [**Stop Calibrating ECG**]

You can print out the square wave and wave scale and then measure the difference between them if necessary. If the difference exceeds 5%, contact your service personnel.

26.4 Calibrating the Touchscreen

- Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password→[Cal. Touchscreen].
- 2. ## will, in turn, appear at different positions of the screen.
- 3. Select each + as it appears on the screen.
- 4. After the calibration is completed, the message [**Screen Calibration Completed!**] is displayed. Select [**Ok**] to confirm the completion of the calibration.

26.5 Calibrating CO₂

For sidestream and microstream CO_2 modules, a calibration is needed every year or when the measured values have a great deviation. For maintream CO_2 module, no calibration is needed. Calibration for sidestream CO_2 module can be performed only when the sidestream module enters the full accuracy mode.


Tools required:

■ A steel gas cylinder with 6±0.05% CO₂ and balance gas N₂

- T-shape connector
- Tubing

Follow this procedure to perform a calibration:

- 1. Make sure that the sidestream or microstream CO₂ module has been warmed up or started up.
- 2. Check the airway for leakage and perform a leakage test as well to make sure the airway has no leakage.
- Select [Main Menu] → [Maintenance >>] → [User Maintenance >>] → enter the required password → [Maintain CO₂ >>] → [Calibrate CO₂ >>].
- 4. In the [Calibrate CO₂] menu, select [Zero].
- 5. After the zero calibration is finished successfully, connect the equipment as follows:



- 6. Turn on and adjust the relief valve to make the flowmeter reads within 10-50mL/min and keeps stable as well.
- 7. In the [Calibrate CO_2] menu, enter the vented CO_2 concentration in the [CO_2] field.
- 8. In the [Calibrate CO₂] menu, the measured CO₂ concentration is displayed. After the measured CO₂ concentration becomes stable, select [Calibrate CO₂] to calibrate the CO₂ module.
- 9. If the calibration is finished successfully, the message [Calibration Completed!] is displayed in the [Calibrate CO₂] menu. If the calibration failed, the message [Calibration Failed!] is displayed. In this case, perform another calibration.

26.6 Calibrating AG


Calibrate the AG module every year or when the measured value has a great deviation.

Tools required:

- Gas bottle, with a certain standard gas or mixture gas. Gas concentration should meet the following requirements: AA>1.5%, CO₂>1.5%, N₂O>40%, O₂>40%, of which AA represents an anesthetic agent. $a/c \le 0.01$ (a is the gas absolute concentration accuracy; c is the gas concentration)
- T-shape connector
- Tubing

Follow this procedure to perform a calibration:

- Select [Main Menu]→[Maintenance >>]→[User Maintenance >>]→enter the required password→[Calibrate AG >>].
- 2. Check the airway and make sure that there are no occlusions or leaks.
 - ◆ Vent the tubing to the air and check if the [**Current FlowRate**] and [**Set FlowRate**] are approximately the same. If the deviation is great, it indicates that there is an occlusion in the tubing. Check the tubing for an occlusion.
 - Perform a leakage test to make sure that the airway has no leakage.
- 3. Connect the test system as follows:



- 4. Open the relief valve and vent a certain standard gas or gas mixture. Adjust the relief valve to make the flowmeter reads within 10-50mL/min and keeps stable as well.
- 5. In the [Calibrate AG] menu, the concentration and flowrate of each measured gas are displayed
 - If the difference between the measured gas concentration and the actual one is very small, a calibration is not needed.
 - ♦ If the difference is great, you should perform a calibration. Select [Calibrate >>] to enter the calibrate menu.
- 6. Enter the vented gas concentration. If you use only one gas for calibration, set other gases' concentration to 0.
- 7. Select [**Calibrate** >>] to start calibration.
- 8. If the calibration is finished successfully, the message [Calibration Completed!] is displayed. If the calibration failed, the message [Calibration Failed!] is displayed. Perform another calibration.



If the O₂ module has been transported for long distance, calibrate it when installing the monitor.

27 Accessories

The accessories listed in this chapter comply with the requirements of IEC 60601-1-2 when in use with the patient monitor.



⚠ WARNING

- Use accessories specified in this chapter. Using other accessories may cause damage to the patient monitor or not meet the claimed specifications.
- Single-use accessories are not designed to be reused. Reuse may cause a risk of contamination and affect the measurement accuracy.
- Check the accessories and their packages for any sign of damage. Do not use them if any damage is detected.

27.1 ECG Accessories

ECG Electrodes

Model	Quantity	Patient Category	Part No.
31499224	10 pieces	Adult	0010-10-12304
2245	50 pieces	Pediatric	9000-10-07469
2258-3	3 pieces	Neonate	900E-10-04880

12-Pin Integrative Trunk Cables

Leadwire	Compatible	Туре	Patient Category	Model	Part No.
supported	with				
5-leadwire	AHA			EA6251B	040-000961-00
5-leadwire	IEC	- Snap, Defibrillator-proof		EA6252B	040-000963-00
3-leadwire	AHA	- Shap, Delibrillator-proof		EA6231B	040-000965-00
3-leadwire	IEC		Adult padiatric	EA6232B	040-000967-00
5-leadwire	AHA		Adult, pediatric	EA6251A	040-000960-00
5-leadwire	IEC	Clip, Defibrillator-proof		EA6252A	040-000962-00
3-leadwire	АНА	Clip, Delibrillator-proof		EA6231A	040-000964-00
3-leadwire	IEC			EA6232A	040-000966-00

12-Pin Separable Trunk Cables

Leadwire	Commetible with	Toma	Detient Catemany	Part No.	
supported	conted Compatible with Type		Patient Category	Part No.	
3-leadwire	AHA, IEC	Defibrillator-proof		0010-30-42720	
3-leadwire	AHA, IEC	ESU-proof	Infant, neonate	0010-30-42724	
3-leadwire	/	Defibrillator-proof		040-000754-00	
3/5-leadwire	AHA, IEC	Defibrillator-proof	Adult padiatric	0010-30-42719	
3/5-leadwire	AHA, IEC	ESU-proof	- Adult, pediatric	0010-30-42723	
12-leadwire	АНА	Defibrillator-proof	- Adult	0010-30-42721	
12-leadwire	IEC	Defibrillator-proof	Adult	0010-30-42722	

Cable Sets

3-Electr	3-Electrode Cable Sets							
Туре	Compatible with	Model	Patient Category	Part No.	Length	Remark		
		EL6304A	· Adult, pediatric	0010-30-42732	1m	Long		
		EL6302A	Addit, pediatric	0010-30-42725	0.6m	/		
	IEC	EL6308A	Pediatric	0010-30-42899	0.6m	/		
		EL6306A	Infant, neonate	0010-30-42897	1m	Long		
Clip		EL6312A	illiant, neonate	040-000149-00	1m	Long		
Clip	EL6303	EL6303A	Adult, pediatric	0010-30-42731	1m	Long		
		EL6301A		0010-30-42726	0.6m	/		
	АНА	EL6307A	Pediatric	0010-30-42898	0.6m	/		
		EL6305A	Infant, neonate	0010-30-42896	1m	Long		
		EL6311A	illiant, neonate	040-000148-00	1m	Long		
		EL6302B	Adult, pediatric	0010-30-42733	1m	Long		
	IEC	EL6308B	Pediatric	0010-30-42901	0.6m	/		
Snap		EL6312B	Infant, neonate	040-000147-00	1m	Long		
Зпар		EL6301B	Adult, pediatric	0010-30-42734	1m	Long		
	АНА	EL6307B	Pediatric	0010-30-42900	0.6m	/		
		EL6311B	Infant, neonate	040-000146-00	1m	Long		

5-Electr	5-Electrode Cable Sets							
Туре	Compatible with	Model	Patient Category	Part No.	Length	Remark		
	IEC	EL6502A		0010-30-42728	0.6m	/		
Clin	iec	EL6504A		0010-30-42730	1m to 1.4m	Long		
Спр	AHA	EL6501A	Adult, pediatric	0010-30-42727	0.6m	/		
		EL6503A		0010-30-42729	1m to 1.4m	Long		
		EL6502B			1.4m for F			
	IEC			0010-30-42736	and N; 1m for	Long		
Snap					others			
					1.4m for RL			
	АНА	EL6501B		0010-30-42735	and LL; 1m	Long		
					for others			

12-Electrode Cable Sets						
_ Compatible	Model	Patient	Part No.	Lawath		
Type	with	Model	Category	Part No.	Length	Remark
	Clip IEC E	EL6802A		0010-30-42903	0.8m	Limb
Clin		EL6804A	Adult	0010-30-42905	0.6m	Chest
Спр		EL6801A	Addit	0010-30-42902	0.8m	Limb
	АПА	EL6803A		0010-30-42904	0.6m	Chest
	IEC	EL6802B		0010-30-42907	0.8m	Limb
Cnan	IEC	EL6804B	Adult	0010-30-42909	0.6m	Chest
Snap	EL6801B	Addit	0010-30-42906	0.8m	Limb	
	AHA	EL6803B		0010-30-42908	0.6m	Chest

27.2 SpO₂ Accessories

Extension Cable

Module type	Remarks	Part No.
Mindray	/	0010-20-42710
Masimo	8 pins, purple connector	040-000332-00
Nellcor	8 pins	0010-20-42712

SpO₂ Sensors

The SpO_2 sensor material that patients or other staff will come into contact with have undertaken the bio-compatibility test and is verified to be in compliance with ISO 10993-1.

Mindray SpO₂ Module					
Туре	Model	Patient Category	Part No.		
	MAX-A	Adult (>30Kg)	0010-10-12202		
Disposable	MAX-P	Pediatric (10 to 50Kg)	0010-10-12203		
Disposable	MAX-I	Infant (3 to 20Kg)	0010-10-12204		
	MAX-N	Neonate (<3Kg), Adult (>40Kg)	0010-10-12205		
	520A	Adult	520A-30-64101		
Cingle patient use	520P	Pediatric	520P-30-64201		
Single patient use	5201	Infant	5201-30-64301		
	520N	Neonate	520N-30-64401		
	DS-100A	Adult	9000-10-05161		
	OXI-P/I	Pediatric, infant	9000-10-07308		
	OXI-A/N	Adult, neonate	9000-10-07336		
	ES-3212-9	Adult	0010-10-12392		
Reusable	518B	Neonate (Multi-sites)	518B-30-72107		
Reusable	518C	Neonate (Multi-sites)	040-000330-00		
	512E	Adult (Finger type)	512E-30-90390		
	512F	Adult (Filiger type)	512F-30-28263		
	512G	Pediatric (Finger type)	512G-30-90607		
	512H	rediatric (Filiger type)	512H-30-79061		

Masimo SpO₂ Module					
Туре	Model	Patient Category	Part No.		
	LNCS NeoPt-L	Pediatric, neonate	0010-10-42626		
	LNCS Neo-L	Neonate	0010-10-42627		
Disposable	LNCS Inf-L	Infant	0010-10-42628		
	LNCS Pdtx	Pediatric	0010-10-42629		
	LNCS Adtx	Adult	0010-10-42630		
	LNCS DCI	Adult	0010-10-42600		
Reusable	LNCS DCIP	Pediatric	0010-10-42634		
	LNCS YI	Adult, pediatric, neonate	0010-10-43016		
Nellcor SpO ₂ Module					
Туре	Model	Patient Category	Part No.		
	MAX-A	Adult (>30Kg)	0010-10-12202		
Disposable	MAX-P	Pediatric (10 to 50Kg)	0010-10-12203		
Disposable	MAX-I	Infant (3 to 20Kg)	0010-10-12204		
	MAX-N	Neonate (<3Kg), Adult (>40Kg)	0010-10-12205		
	DS-100A	Adult	9000-10-05161		
Damakla	OXI-P/I	Pediatric, infant	9000-10-07308		
Reusable	OXI-A/N	Adult, neonate	9000-10-07336		
	D-YS	Adult, pediatric, infant, neonate	0010-10-12476		

Wavelength emitted by the sensors:

- 520A, 520P, 520I, 520N, 518B, 518C, 512E, 512F, 512G, 512H, and ES-3212-9: red light: 660 nm, infrared light: 905 nm.
- LNCS NeoPt-L, LNCS Neo-L, LNCS Inf-L, LNCS Pdtx, LNCS Adtx, LNCS DCI, LNCS DCIP, and LNCS YI: red light: 660 nm, infrared light: 940 nm.
- MAX-A, MAX-P, MAX-I, MAX-N, DS-100A, OXI-P/I, OXI-A/N, and D-YS: red light: 660 nm, infrared light: 890 nm.

The maximum photic output consumption of the sensor is less than 18 mW. The information about the wavelength range and maximum photic output consumption can be especially useful to clinicians, for example, clinicians performing photodynamic therapy.

27.3 NIBP Accessories

Tubing

Туре	Patient Category	Part No.
Reusable	Adult, pediatric, infant	6200-30-09688
neusable	Neonate	6200-30-11560

Reusable Cuff

Model	Patient	Measurement	Limb Circumference (cm)	Bladder Width	Part No.
oue.	Category	Site		(cm)	- urerror
CM1200	Small infant		7 to 13	5.8	115-002480-00
CM1201	Infant		10 to 19	9.2	0010-30-12157
CM1202	Pediatric	Arm	18 to 26	12.2	0010-30-12158
CM1203	Adult		24 to 35	15.1	0010-30-12159
CM1204	Large adult	1	33 to 47	18.3	0010-30-12160
CM1205	Thigh	Thigh	46 to 66	22.5	0010-30-12161
CM1300	Small infant		7 to 13	5.8	040-000968-00
CM1301	Infant		10 to 19	9.2	040-000973-00
CM1302	Pediatric	Arm	18 to 26	12.2	040-000978-00
CM1303	Adult		24 to 35	15.1	040-000983-00
CM1304	Large adult	7	33 to 47	18.3	040-000988-00
CM1305	Adult thigh	Thigh	46 to 66	22.5	040-000993-00

Single-Patient Cuff

Model	Patient	Measurement	Limb Circumference (cm)	Bladder Width	Part No.
Model	Category	Site	Limb circumference (cm)	(cm)	Tarrito.
CM1500A			3.1 to 5.7	2.2	001B-30-70677
CM1500B			4.3 to 8.0	2.9	001B-30-70678
CM1500C	Neonate		5.8 to 10.9	3.8	001B-30-70679
CM1500D			7.1 to 13.1	4.8	001B-30-70680
CM1500E		Arm	8 to 15	/	001B-30-70681
CM1501	Infant		10 to 19	7.2	001B-30-70682
CM1502	Pediatric		18 to 26	9.8	001B-30-70683
CM1503	Adult		25 to 35	13.1	001B-30-70684
CM1504	Large adult		33 to 47	16.5	001B-30-70685
CM1505	Adult	Thigh	46 to 66	20.5	001B-30-70686

27.4 Temp Accessories

Temp Cable

Туре	Model	Remark	Part No.
Extension cable (reusable)	MR420B	Applicable to sensor MR411 and MR412	0011-30-37391
TEMP adapter cable (2-pin to	MR421	,	0010 20 42056
audio)	WR421	7	0010-30-43056

Temp Probes

Туре	Model	Patient Category	Measurement Site	Part No.
	MR401B	Adult	Esophageal/Rectal	0011-30-37392
eusable	MR403B	Adult	Skin	0011-30-37393
eusable	MR402B	D. P. C. C.	Esophageal/Rectal	0011-30-37394
	MR404B	Pediatric, infant	Skin	0011-30-37395
Disposable	MR411	Adult, pediatric, infant	Esophageal/Rectal	0011-30-37398
Disposable	MR412	Addit, pediatric, illiant	Skin	0011-30-37397

27.5 IBP/ICP Accessories

Material	Part No.	
IBP adapter cable		0010-30-43055
IBP extended cable with d	IBP extended cable with dual-receptacle	
Accessories Kit No.	Components	Part No.
	IM2201 12Pin IBP Cable	001C-30-70759
6800-30-50876	Disposable Transducer	0010-10-42638
(Hospira)	Steady Rest for IBP Transducer and Clamp	M90-000133
	Steady Rest for IBP Transducer and Clamp	M90-000134
6800-30-50877	IM2202 12Pin IBP Cable	001C-30-70757
	Disposable Pressure Transducer	6000-10-02107
(BD)	Transducer/Manifold Mount	0010-10-12156
ICP		
Model	Material	Part No.
Gaeltec TYPE.S13	12Pin ICP cable	0010-30-42742
Gaeltec ICT/B	Intracranial Pressure Transducer	0010-10-12151

It is proved through tests that the following accessories are compatible with the patient monitor. Only the accessories proceeded by "*" are available from our company. If you want to purchase other accessories, contact respective manufacturers and make sure if these accessories are approved for sale in local.

Manufacturer	Accessories		
	MX961Z14 Logical Cable, to be used in connection with the Adapter Cable (0010-20-42795)		
	MX960 Reusable Transducer Kit		
Smith Medical	MX261 Logical Clamp For Transducer Bracket		
(Medex)	MX262 Logical Clamp For 2 Transducer Mount Plates		
	MX960E6441 Logical Transducer Mounting Plate		
	(More Logical Clamps are available from Medex. For detailed information, contact Medex.)		
	IBP Reusable Cable (REF: 5203511), to be used in connection with the Adapter Cable (0010-20-42795)		
Braun	Combitrans Monitoring Set (contact Braun for detailed information)		
Diduii	Combitrans Attachment Plate Holder (REF:5215800)		
	Combitrans Attachment Plate (contact Braun for detailed information)		
	*Truck cable (0010-21-43082)		
Memscap	SP844 Physiological Pressure Transducer		
	844-26 Monitoring Line Set		

	84X-49 Mounting Bracket	
	Reusable Blood Pressure Monitor Interface Cable (REF: 650-206)	
	Deltran Disposable Pressure Transducer System	
Hank	(More Deltran sensors are available from Utah. For detailed information, contact Utah.)	
Utah	Pole Mount Unit (ERF: 650-150)	
	Deltran Three Slot Organizer, Attaches to I.V. Pole Mount (REF: 650-100)	
	Deltran Four Slot Organizer, Attaches to I.V. Pole Mount (REF: 650-105)	
	* IBP Truwave Reusable Cable (0010-21-12179)	
	Pressure Monitoring Kit With Truwave Disposable Pressure Transducer.	
Edwards	(More Truwave sensors are available from Edwards. For detailed information, contact Edwards.)	
	DTSC IV Pole Clamp for Model DTH4 Backplate Holder	
	DTH4 Disposable Holder for DPT	

27.6 C.O. Accessories

Model	Material	Part No.
CO7702	12Pin C.O. cable.	0010-30-42743
SP4042	IT Sensor	6000-10-02079
SP5045	IT Sensor Housing	6000-10-02080
MX387	12CC Control Syringe W/1CC Stop W/Rotator	6000-10-02081

27.7 CO₂ Accessories

Sidestream CO₂ module

Material	Patient Category	Remark	Part No.
DRYLINE Watertrap	Adult, pediatric	Reusable	9200-10-10530
DRYLINE Watertrap	Neonate	Neusable	9200-10-10574
Sampling Line Adult 2.5m	Adult, pediatric		9200-10-10533
Sampling Line, Neonate, 2.5m	Neonate	Disposable	9200-10-10555
Adult Nasal CO₂ Sample Cannula	Adult		M02A-10-25937
Pediatric Nasal CO ₂ Sample Cannula	Pediatric		M02A-10-25938
Infant Nasal CO2 Sample Cannula	Neonate		M02B-10-64509
DRYLINE Airway Adapter	/	Straight,	9000-10-07486
DNTLINE All Way Adapter	1	disposable	9000-10-07460

Microstream CO₂ Module

Disposable Airway Sampling Line			
Model	Patient Category	Remark	Part No.
XS-04620		/	0010-10-42560
XS-04624	A dult in a diatria	Humidified	0010-10-42561
007768	Adult, pediatric	Long	0010-10-42563
007737		Long, humidified	0010-10-42564
006324	Infant, Neonate	Humidified	0010-10-42562
007738	miant, Neonate	Long, humidified	0010-10-42565

Disposable Nasal Sampling Line			
Model	Patient Category	Remark	Part No.
009818		/	0010-10-42566
009822	Adult, intermediate	Plus O ₂	0010-10-42568
009826		Long, plus O ₂	0010-10-42570
008174		/	0010-10-42577
008177	Adult	Humidified	0010-10-42572
008180		Humidified, plus O ₂	0010-10-42575
007266		/	0010-10-42567
008175		/	0010-10-42578
008178	- Pediatric	Humidified	0010-10-42573
008181	Pediatric	Humidified, plus O ₂	0010-10-42576
007269		Plus O ₂	0010-10-42569
007743		Long, plus O ₂	0010-10-42571
008179	Infant, Neonate	Humidified	0010-10-42574

Mainstream CO₂ Module

Material	Model	Patient Category	Remark	Part No.
	6063		Disposable	0010-10-42662
	6421	Adult, pediatric	Disposable, with	0010-10-42663
	0421	Addit, pediatric	mouthpiece	0010-10-42003
Airway adapter	7007		Reusable	0010-10-42665
	6312	Neonate, pediatric	Disposable	0010-10-42664
	7053		Reusable	0010-10-42666
	9960STD	Adult	/	0010-10-42670
Mask	9960LGE	Adult	Adult large	0010-10-42669
	9960PED	Pediatric	/	0010-10-42671
Cable management straps	6934-00	/	/	0010-10-42667
Sensor holding clips	8751	/	/	0010-10-42668
Sensor	1022386	Adult, pediatric,	Reusable	6800-30-50760
	1022300	neonate	Reusable	

27.8 AG Accessories

Material	Patient Category	Remark	Part No.
Watertrap	Adult, pediatric Reusable		9200-10-10530
watertrap	Neonate	neusable	9200-10-10574
Compiling line	Adult, pediatric	- Disposable	9200-10-10533
Sampling line	Neonate	Disposable	9200-10-10555
Airway adapter	Adult, pediatric, neonate	Disposable, straight	9000-10-07486
All way adapter	Adult, pediatric, neonate	Disposable, elbow	9000-10-07487

27.9 Others

Material	Part No.
Lithium battery	022-00008-00
Power cord	509B-10-05996
U.K. power cord	DA8K-10-14453
European power cord	DA8K-10-14454
U.S. power cord	DA8K-10-14452
Brazil power cord (250V, 10A, 3M)	009-001075-00
South Africa Power Cable (250V, 16A, 3M)	009-001791-00
India power cord	0000-10-10903
Grounding cable	1000-21-00122
LCD display, 17"	0000-10-11284
USB drive, 4G	023-000217-00
03b unve, 4d	023-000218-00
Recorder	TR6F-30-67306
Thermal paper	A30-000001
Wall mount bracket for external display	0010-30-42956
Rolling bracket	045-000670-00
Wall mount	045-000672-00
Bedrail Hook subassembly (iPM 10/iPM 12)	115-012698-00
Bedrail Hook subassembly (iPM 8)	115-012697-00
Transition Plate Kit	115-012695-00
Beneview data output package (CD, Cable, User's Guide)	6800-30-51213

FOR YOUR NOTES		

A Product Specifications

A.1 Monitor Safety Specifications

A.1.1 Classifications

The patient monitor is classified, according to IEC60601-1:

Type of protection against electrical shock	Class I, equipment energized from an external and internal electrical	
Type of protection against electrical shock	power source.	
Doggoo of protection against alectrical sheek	Type BF defibrillation proof for CO₂ and AG monitoring.	
Degree of protection against electrical shock	Type CF defibrillation proof for ECG, RESP, TEMP, SpO ₂ , NIBP, IBP and C.O	
Mode of operation	Continuous	
Degree of protection against harmful ingress	IPX1	
of water	IFAI	

A.1.2 Environmental Specifications

Main unit			
Item	Operating conditions	Storage conditions	
Temperature (°C) 0 to 40	iPM 8: -30 to 70		
	0 10 40	iPM 12/ iPM 10: -20 to 60	
Relative humidity (noncondensing)	15% to 95%	10% to 95%	
Barometric (kPa)	57.0 to 107.4	16.0 to 107.4	

Microstream CO ₂ module		
Item	Operating conditions	Storage conditions
Temperature ($^{\circ}\mathbb{C}$)	0 to 40	-20 to 60
Relative humidity (noncondensing)	15% to 95%	10% to 95%
Barometric (kPa)	57.3 to 105.3	57.3 to 105.3

Sidestream CO ₂ module		
Item	Operating conditions	Storage conditions
Temperature ($^{\circ}$ C)	5 to 40	-20 to 60
Relative humidity (noncondensing)	15% to 95%	10% to 95%
Barometric (kPa)	57.3 to 105.3	57.3 to 105.3

Mainstream CO₂ module		
Item	Operating conditions	Storage conditions
Temperature ($^{\circ}\mathbb{C}$)	0 to 40	-20 to 60
Relative humidity (noncondensing)	15% to 90%	10% to 90%
Barometric (kPa)	57.0 to 107.4	53.3 to 107.4

AG module			
Item	Operating conditions	Storage conditions	
Temperature (°C)	10 to 40	-20 to 60	
Relative humidity (noncondensing)	15% to 95%	10% to 95%	
Barometric (kPa)	70 to 107.4	70 to 107.4	

NOTE

The environmental specifications of unspecified parameters are the same as those of the main unit of iPM
 12/iPM 10.

A.2 Power Supply Specifications

AC power		
Line voltage		100 to 240 VAC
Current		1.3 to 0.5 A
Frequency		50/60 Hz
DC power (availa	ıble for iPM 8 oı	nly)
Voltage		12 VDC
Current		3.5 A
Battery		
Battery Type		Chargeable Lithium-lon, 11.1DVC, 2.6 Ah
		4 hours
		when powered by a new fully-charged battery (25 °C, SpO₂, ECG, disconnected from
	iPM 12	Temp cable, Auto NIBP measurements at intervals of 15 minutes)
	IFIVI 12	8 hours
		when powered by two new fully-charged batteries (25 °C , SpO ₂ , ECG, disconnected from
Run time		Temp cable, Auto NIBP measurements at intervals of 15 minutes)
nuii tiille		4 hours
	iPM 10	when powered by a new fully-charged battery (25 °C, SpO ₂ , ECG, disconnected from
		Temp cable, Auto NIBP measurements at intervals of 15 minutes)
		6 hours
	iPM 8	when powered by a new fully-charged battery (25 °C, SpO₂, ECG, disconnected from
		Temp cable, Auto NIBP measurements at intervals of 15 minutes)
Charma time		Less than 3 hours to 90%, and less than 4 hours to 100% when the monitor is off.
Charge time		Less than 8 hours to 90%, and less than 12 hours to 100% when the monitor is on.
Shutdown delay		at least 20 min (after a low battery alarm first occurs)

A.3 Physical Specifications

Model	Size (Width×Hight×Thickness)	Weight	Remark
iPM 12	318mm×274mm×128mm	≤4.2 kg	Standard parameters,
iPM 10	282mm×252mm×128mm	≤3.6 kg	including touchscreen and
iPM 8	238mm×225mm×128mm	≤3.2 kg	recorder.

A.4 Hardware Specifications

A.4.1 Display

Host display				
	Screen Size (diagonal)		Screen type	Resolution
iPM 12	12.1"			
iPM 10	10.4" 8.4"		color LED display	800×600 pixels
iPM 8				
External display				
Screen type	een type Med		ical-grade TFT LCD	
Screen Size	15", 1		, 17" or above	
Resolution	800×		×600 pixels	
EMC		MPR II, CISPR 11B		
Third certificate		UL, C-UL, TUV, CE, FCC		

A.4.2 Recorder

Method	Thermal dot array
Horizontal resolution	16 dots/mm (25 mm/s paper speed)
Vertical resolution	8 dots/mm
Paper width	50 mm
Paper length	20 m
Paper speed	25 mm/s or 50 mm/s with accuracy within ±5%
Number of waveform channels	1, 2, or 3 (optional)

A.4.3 LEDs

Alarm lamp	1 (two color coded: yellow and red)
Power on LED	1 (green)
AC power LED (iPM 12/ iPM 10)	1 (green)
AC/DC power LED (iPM 8)	1 (green)
Battery LED	1 (green)

A.4.4 Audio Indicator

Speaker	Give alarm tones (45 to 85 dB), key tones, QRS tones; support PITCH TONE and	
Speaker	multi-level tone modulation; alarm tones comply with IEC60601-1-8.	

A.4.5 Monitor Interface Specifications

Power	1 AC power input connector	
	1 DC power input connector (for iPM 8 only)	
Wired network	1 RJ45 connector, 100 Base-TX, IEEE 802.3	
USB	2 connectors, USB 2.0	
Equipotential Grounding Terminal	1	
Multifunctional connector	1	
VGA connector	1	

A.4.6 Outputs

Standard Meets the requirements of IEC60601-1 for short-circuit protection and leakage current ECG Analog Output Bandwidth	Analog Output				
ECG Analog Output Bandwidth (-3dB; reference frequency: 10Hz) Surgical mode: 0.05 to 40 Hz Surgical mode: 1 to 20 Hz ST mode: 0.05 to 40 Hz ST mode: 0.05 to 40 Hz ST mode: 0.05 to 40 Hz Max transmission delay SE ms (in diagnostic mode, and with Notch off) Sensitivity IV/mV ±5% Pace enhancement Signal amplitude: Voh≥2.5V Pulse width: 10ms±5% Signal rising and falling time: ≤100µs BBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) Max transmission delay 30 ms Sensitivity IV/100 mmHg ±5% Nurse Call Signal Output mode Electrical requirements Se0W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage Sensitivity Normally open or normally contact (optional) Defib Sync Pulse Output impedance Max time delay S ms (R-wave peak to leading edge of pulse) Amplitude Pulse width 100 ms ±10% Rising and falling time ≤1 ms Alarm output (Network connector) Alarm output (Network connector) Alarm output (Network connector) Alarm output (Network connector) The alarm delay time from patient monitor) The alarm delay time from patient monitor	Analog Output	Martatha manipusasi CIEC	COCOL 1 for the state size it was to still a stall a		
ECG Analog Output Bandwidth (-3dB; reference frequency: 10Hz) Diagnostic mode: Diagnostic mode; Diagnostic mode, and with Notch off) Max transmission delay (PACE rejection/enhancement Signal amplitude: Voh≥2.5V Pulse width: 10ms±5% Signal rising and falling time: ≤100µs Bandwidth (-3dB; reference frequency: 1Hz) (Prevented of the patient mode) Diagnostic mode and provided of the patient mode of the patien	Standard				
Bandwidth (Annitor mode: 0.05 to 150 Hz Monitor mode: 0.5 to 40 Hz Monitor mode: 0.05 to 40 Hz Monitor mode, and with Notch off) Sensitivity 10 Monitor mode, and with Notch off) Pace enhancement Signal amplitude: Voh≥2.5V Pulse width: 10ms±5% Signal rising and falling time: ≤100μs IBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 11/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defits Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: <0.5 V, receiving a maximum of 5 mA input current; Low level: <0.5 V, receiving a maximum of 5 mA input current. Pulse width 100 ms ±10% Rising and falling time ≤1 ms Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time from the patient monitor to remote equipment is ≤2		current			
Bandwidth (-3dB; reference frequency: 10Hz) Monitor mode: 1 to 20 Hz (-3 to 40 Hz) Max transmission delay 25 ms (in diagnostic mode, and with Notch off) Sensitivity 1V/mV±5% PACE rejection/enhancement Signal amplitude: Voh≥2.5V Pulse width: 10ms±5% (signal rising and falling time: ≤100μs IBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse 100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	ECG Analog Output				
Surgical mode: 1 to 20 Hz		Diagnostic mode:	0.05 to 150 Hz		
Max transmission delay 25 ms (in diagnostic mode, and with Notch off) Sensitivity 1V/mV ±5% Pace enhancement Signal amplitude: Voh≥2.5V Pulse width: 10ms±5% Signal arising and falling time: ≤100µs IBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Bandwidth		0.5 to 40 Hz		
Max transmission delay 25 ms (in diagnostic mode, and with Notch off) Sensitivity 1V/mV ±5% PACE rejection/enhancement Signal amplitude: Voh≥2.5V Pulse width: 10ms±5% Signal rising and falling time: ≤100µs IBP Analog Output (For iPM 12 only) DC to 40 Hz Bandwidth (-3dB; reference frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Vulue mode Belay Electrical requirements ≤ 60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Some (R-wave peak to leading edge of pulse) Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	(-3dB; reference frequency: 10Hz)	Surgical mode:	1 to 20 Hz		
Sensitivity 1V/mV±5% PACE rejection/enhancement Pace enhancement Signal amplitude: Voh≥2.5V Pulse width: 10ms±5% Signal rising and falling time: ≤100µs IBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 1 V/100 mmHg±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width 100 ms ±10% Rising and falling time ≤1 ms Alarm output (Network connector) The alarm delay time form the patient monitor to remote equipment is ≤2		ST mode:	0.05 to 40 Hz		
PACE rejection/enhancement Signal amplitude: Voh≥2.5V Pulse width: 10ms±5% Signal rising and falling time: ≤100μs IBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) Max transmission delay 30 ms Sensitivity 1 1/100 mmHg±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width Rising and falling time Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤2	Max transmission delay	25 ms (in diagnostic mode, and	d with Notch off)		
Signal amplitude: Voh≥2.5V PLISE width: 10ms±5% Signal rising and falling time: ≤100µs IBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width 100 ms ±10% Rising and falling time ≤1 ms Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor	Sensitivity	1V/mV ±5%			
Pulse width: 10ms±5% Signal rising and falling time: ≤100μs IBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) Max transmission delay Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance Max time delay 35 ms (R-wave peak to leading edge of pulse) High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width 100 ms ±10% Rising and falling time Alarm output (Network connector) The alarm delay time form the patient monitor to remote equipment is ≤2		Pace enhancement			
Pulse width: 10ms±5% Signal rising and falling time: ≤100μs Bandwidth (-3dB; reference frequency:1Hz) Max transmission delay Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance Max time delay 35 ms (R-wave peak to leading edge of pulse) High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width Rising and falling time Alarm output (Network connector) The alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤2	DACE valenties /only a server	Signal amplitude: Voh≥2.5V			
IBP Analog Output (For iPM 12 only) Bandwidth (-3dB; reference frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	PACE rejection/enhancement	Pulse width: 10ms±5%			
Bandwidth (-3dB; reference frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width 100 ms ±10% Rising and falling time ≤1 ms Alarm output (Network connector) The alarm delay time form the patient monitor to remote equipment is ≤2		Signal rising and falling time: s	≤100μs		
frequency:1Hz) DC to 40 Hz Max transmission delay 30 ms Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	IBP Analog Output (For iPM 12 only)				
frequency:1Hz) 30 ms Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Bandwidth (-3dB; reference	DC +- 40 H-			
Sensitivity 1 V/100 mmHg ±5% Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width 100 ms ±10% Rising and falling time ≤1 ms Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤2	frequency:1Hz)	DC to 40 HZ			
Nurse Call Signal Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Max transmission delay	30 ms			
Output mode Relay Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Sensitivity	1 V/100 mmHg ±5%			
Electrical requirements ≤60W, ≤2A, ≤36VDC, ≤25VAC Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Nurse Call Signal				
Isolation voltage 1500 VAC Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Output mode	Relay			
Contact type Normally open or normally contact (optional) Defib Sync Pulse Output impedance $\leq 100\Omega$ Max time delay 35 ms (R-wave peak to leading edge of pulse) High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width $100 \text{ ms} \pm 10\%$ Rising and falling time $\leq 1 \text{ ms}$ Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤ 2	Electrical requirements	≤60W, ≤2A, ≤36VDC, ≤25VAC			
Defib Sync Pulse Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Isolation voltage	1500 VAC			
Output impedance ≤100Ω Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Contact type	Normally open or normally co	ntact (optional)		
Max time delay 35 ms (R-wave peak to leading edge of pulse) Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current.	Defib Sync Pulse				
Amplitude High level: 3.5 to 5 V, providing a maximum of 10 mA output current; Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width 100 ms ±10% Rising and falling time ≤1 ms Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤2	Output impedance	≤100Ω			
Amplitude Low level: $< 0.5 \text{ V}$, receiving a maximum of 5 mA input current. Pulse width $100 \text{ ms} \pm 10\%$ Rising and falling time $\le 1 \text{ ms}$ Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤ 2	Max time delay	35 ms (R-wave peak to leading	gedge of pulse)		
Low level: < 0.5 V, receiving a maximum of 5 mA input current. Pulse width $100 \text{ ms} \pm 10\%$ Rising and falling time $\le 1 \text{ ms}$ Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤ 2		High level: 3.5 to 5 V, providing a maximum of 10 mA output current;			
Rising and falling time ≤1 ms Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤2	Amplitude				
Alarm output (Network connector) Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤2	Pulse width				
Alarm delay time from patient monitor The alarm delay time form the patient monitor to remote equipment is ≤2	Rising and falling time	≤1 ms			
	Alarm output (Network connector)				
to remote equipment seconds, measured at the patient monitor's signal output connector.	Alarm delay time from patient monitor	The alarm delay time form the	patient monitor to remote equipment is ≤2		
	to remote equipment	seconds, measured at the patient monitor's signal output connector.			

A.5 Data Storage

	Trends: 120 hours, at 1 min resolution	
Trends	Mid-length trends: 4 hours, at 5 s resolution	
	Minitrends: 1 hour, at 1 s resolution	
Parameter alarms	100 alarms and manual events and related parameter waveforms. The waveform	
Parameter alarms	recording length can be 8s, 16s, or 32s.	
A will a consider	100 arrhythmia events and relate waveforms and parameters. The waveform	
Arrh. events	recording length can be 8s, 16s, or 32s.	
NIBP measurements	1000 sets	
5 11 15 1	48 hours at maximum. The specific storage time depends on the waveforms stored	
Full-disclosure waveforms	and the number of stored waveforms.	

A.6 Wireless Network

Standards	IEEE 802.11b/g, support Wi-Fi	
Frequency range	2.400 to 2.4835 GHz	
Transmission range	≥50 m (visible range)	

A.7 Measurement Specifications

The adjustable range of alarm limits is the same with the measurement range of signals unless otherwise specified.

A.7.1 ECG

ECG			
Standards	Meet standards of EC11, EC13, EN60601-2-27/IEC60601-2-27 and IEC60601-2-25		
	3-lead: I, II, III 5-lead: I, II, III, aVR, aVL, aVF, V		
Lead set			
	12-lead: I, II, III, aVR, aVL, aVF, V1	to V6 (for iPM 12 and iPM 10 only)	
ECG standard	AHA, IEC		
Display consitivity	1.25 mm/mV (X0.125), 2.5 mm/	mV (X0.25), 5 mm/mV (X0.5), 10 mm/mV (X1), 20	
Display sensitivity	mm/mV (X2), 40 mm/mV (X4),	Auto	
Sweep speed	6.25 mm/s, 12.5 mm/s, 25 mm/	s, 50 mm/s	
	Diagnostic mode:	0.05 to 150 Hz	
Bandwidth (-3dB)	Monitor mode:	0.5 to 40 Hz	
bandwidth (-5db)	Surgical mode:	1 to 20 Hz	
	ST mode:	0.05 to 40 Hz	
	Diagnostic mode:	>90 dB	
Common mode rejection ratio	Monitor mode:	>105 dB	
Common mode rejection ratio	Surgical mode:	>105 dB	
	ST mode:	>105 dB	
	50/60 Hz		
Notch	Monitor, ST and surgical mode: Notch turns on automatically. Diagnostic mode:		

	Notch is turned on/off manually		
Differential input impedance	≥5 MΩ		
Input signal range	±8 mV (peak-to-peak value)		
Accuracy of reappearing input signal	Based on EC11 to determine	system total error and frequency response.	
Electrode offset potential tolerance	±500 mV		
Lead-off detection current	Measuring electrode: <0.1 μ	A	
Lead-on detection current	Drive electrode: <1 μA		
Input offset current	Measuring electrode: ≤0.1 μ	A	
input onset current	Drive electrode: ≤1 μA		
Baseline recovery time	<5 s (after defibrillation)		
Patient leakage current	<10 uA		
Calibration signal	1mV (peak-to-peak value)		
	Cut mode: 300 W		
ESU protection	Coagulate mode: 100 W		
L30 protection	Recovery time: ≤10 s		
	In compliance with the requirements in clause 4.2.9.14 of ANSI/AAMI EC 13:2002		
	Based on the test method in clause 5.2.9.14 of EC 13, use ECG lead wires which a in compliance with AAMI. Compared with ECG baseline, the noise of peak to pea		
ESU noise suppression			
	value ≤2 mV.		
Pace Pulse			
	Pace pulses meeting the following conditions are labelled with a PACE marker:		
Pace pulse markers	Amplitude:	±2 to ±700 mV	
race paise markers	Width:	0.1 to 2 ms	
	Rise time:	10 to 100 μs	
	When tested in accordance v	with the ANSI/AAMI EC13-2002: Sections 4.1.4.1 and	
	4.1.4.3, the heart rate meter	rejects all pulses meeting the following conditions.	
Pace pulse rejection	Amplitude:	±2 to ±700 mV	
	Width:	0.1 to 2 ms	
	Rise time:	10 to 100 μs	
Pacer pulse detector rejection of fast	10V/s RTI when measured in accordance with ANSI/AAMI EC13-2002 Section		
ECG signals	4.1.4.3.		

Mindray algorithm

HR			
	Neonate: 15 to 350 bpm		
Measurement range	Pediatric:	15 to 350 bpm	
	Adult:	15 to 300 bpm	
Resolution	1 bpm		
Accuracy	±1 bpm or ±1%, whichever is greater.		
Sensitivity	200μV (lead II)		
	In compliance with the requirements in Clause 4.1.2.1 d)of ANSI/AAMI EC13-2002,		
	the following method is used:		
HR averaging method	If the last 3 consecutive RR intervals are greater than 1200 ms, the 4 most recent RR		
	intervals are averaged to compute the HR. Otherwise, heart rate is computed by subtracting the maximum and minimum ones from the most recent 12 RR interv		

	and then averaging them.			
	The HR value displayed on the monitor screen is updated every second.			
		In compliance with the requirements in Clause 4.1.2.1 e)of ANSI/AAMI EC13-2002,		
		the heart rate after 20 seconds of stabilization is displayed as follows:		
Response to irregular rhythm	Ventricular bigeminy (3a):	80±1 bpm		
	Slow alternating ventricular			
		Rapid alternating ventricular bigeminy (3c): 120±1 bpm		
	Bidirectional systoles (3d):	90±2 bpn	1	
	Meets the requirements of A	ANSI/AAMI EC13-2002: Se	ection 4.1.2.1 f).	
Response time to heart rate change	From 80 to 120 bpm: less th	an 11 s		
	From 80 to 40 bpm: less tha	n 11 s		
	Meets the requirements of	ANSI/AAMI EC13-2002: se	ction 4.1.2.1 g).	
	Waveform			
	4ah - range:	11 s		
Time to alarm for tachycardia	4a - range:	11 s		
(not available in USA)	4ad - range:	4ad - range: 11 s		
	Waveform 4bh - range:	Waveform 4bh - range: 11 s		
	4b - range:	4b - range: 11 s		
	4bd - range:	4bd - range: 11 s		
	When the test is performed	based on part 4.1.2.1 c)o	f ANSI/AAMI EC 13-2002, the	
	heart rate meter will reject a	all 100 ms QRS complexes	s with less than 1.2 mV of	
Tall T-wave rejection capability	amplitude, and T waves with T-wave interval of 180 ms and those with Q-T interval			
	of 350 ms.			
ST Segment Analysis (Not available	in USA)			
Measurement range	-2.0 to 2.0 mV			
Accuracy	-0.8 to 0.8 mV:	± 0.02 mV or $\pm 10\%$, whichever is greater.		
Accuracy	Beyond this range:	Beyond this range: Not specified.		
Refreshing rate	10 s	10 s		
Resolution	10 μV			
Alarm limit	Range	Range Step		
HR High	(low limit + 2) to 300 bpm		· 1bpm	
HR Low	15 to (high limit – 2) bpm			
ST High	(low limit +0.2) to 2.0 mV		0.1mV	
ST Low	-2.0 to (high limit – 0.2) mV		0.11111	

Mortara algorithm

Only the differences from the Mindray algorithm are listed.

HR	
	In compliance with the requirements in Clause 4.1.2.1 d)of ANSI/AAMI EC13-2002,
	the following method is used:
HR averaging method	Heart rate is computed by averaging the most recent 16 RR intervals, unless the HR
	by averaging the most recent 4 heart beats is less than or equals to 48.
	The HR value displayed on the monitor screen is updated every second.

	Meets the requirements of ANSI/AAMI EC13-2002: section 4.1.2.1 g).		
	Waveform		
	4ah – range:	11 s	
Time to alarm for tachycardia	4a – range:	11 s	
	4ad – range:	11 s	
	4bh – range:	11 s	
	4b – range:	11 s	
	4bd – range:	11 s	
Arrhythmia Analysis Classifications	Asystole, Vfib, Vtac, Vent. Rhythm, Couplet, VT>2, Bigeminy, Trigeminy, R on T,		
Armytiinia Analysis Classifications	Multif. PVC, Irr. Rhythm, Tachy, Brady, Missed Beats, PNP, PNC		
ST Segment Analysis			
Refreshing rate	per 16 heartbeats		

A.7.2 Resp

7117 12 11C3P				
Technique	Trans-thoracic impedance			
Lead	Options are lead I and II. The default is lead II.			
Respiration excitation waveform	<300 μA RMS, 62.8 kHz	z (±10%)		
Respiration impedance range	0.3 to 5Ω			
Baseline impedance range	200 to 2500Ω (using ar	n ECG cable with 1kΩ resista	ance)	
Bandwidth	0.2 to 2.5 Hz (-3 dB)			
Sweep speed	6.25 mm/s, 12.5 mm/s	or 25 mm/s		
Respiration Rate	Respiration Rate			
Management	Adult: 0 to 120 rpm			
Measurement range	Pediatric, neonate: 0 to 150 rpm			
Resolution	1 rpm			
Accuracy	7 to 150 rpm: ±2 rpm or ±2%, whichever is greater			
Accuracy	0 to 6 rpm: Not specified.			
Apnea alarm time	10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s			
Alarm limit	Range (rpm) Step (rpm)		Step (rpm)	
DD Himb	Adult, pediatric: (I	ow limit + 2) to 100		
RR High	Neonate: (I	ow limit + 2) to 150	1	
RR Low	0 to (high limit – 2)			

A.7.3 SpO₂

Alarm limit	Range (%)	Step (%)	
SpO₂ High	(low limit + 2) to 100		
	Mindray, Masimo: Desat to (high limit – 2)		
SpO₂ Low	Nellcor: Desat or20 (whichever is greater) to (high limit – 2)		1
Desat	Mindray, Masimo:	0 to (high limit – 2)	
Desat	Nellcor:	0 to (high limit – 2)	

Mindray SpO₂ Module

Standards	Meet standards of ISO9919				
*Measurement accuracy verification: The SpO ₂ accuracy has been verified in human experiments by comparing with arterial					
blood sample reference measured with a	CO-oximeter. Pulse oximeter measurement are star	tistically distributed	d and about		
two-thirds of the measurements are expe	ected to come within the specified accuracy range o	compared to CO-ox	rimeter		
measurements.					
Measurement range	0 to 100%				
Resolution	1%				
70 to 100%: ±2% (measured without motion in adult/pediatric mode)					
A 551,175 51.	70 to 100%: ±3% (measured without motion in neonate mode)				
Accuracy	70 to 100%: ±3% (measured with motion)				
	0% to 69%: Not specified.				
*Studies were performed to validate the accuracy of Pulse Oximeter with neonatal SpO ₂ sensors by contrast with a					
CO-Oximeter. Some neonates aged from 1 day to 30 days with a gestation age of 22 weeks to full term were involved in this					
study. The statistical analysis of data of this study shows the accuracy (Arms) is within the stated accuracy specification. Please					
see the following table.					
Sensor type	Totally neonates Data Arms				
518B	97 (51 male & 46 female) 200 pairs 2.38%				

Sensor type	Totally neonates	Data	Arms
518B	97 (51 male & 46 female)	200 pairs	2.38%
520N	122 (65 male & 57 female)	200 pairs	2.88%
The Pulse Oximeter with neonatal SpO₂ sensors was also validated on adult subjects.			
Refreshing rate	1 s		
	7 s (When the sensitivity is set to High)		
SpO₂ averaging time	9 s (When the sensitivity is set to Medium)		
	11 s (When the sensitivity is set to Low)		

Masimo SpO₂ Module

Standards	Meet standards of ISO9919
Measurement range	1 to 100%
Resolution	1%
	70 to 100%: ±2% (measured without motion in adult/pediatric mode)
Accuracy	70 to 100%: ±3% (measured without motion in neonate mode)
Accuracy	70 to 100%: ±3% (measured with motion)
	1% to 69%: Not specified.
Refreshing rate	1 s
SpO₂ averaging time	2-4 s, 4-6 s, 8 s, 10 s, 12 s, 14 s, 16 s
Low particion conditions	Pulse amplitude: >0.02%
Low perfusion conditions	Light penetration: >5%
Low perfusion SpO₂ accuracy	±2%

Nellcor SpO₂ Module

Standards	Meet standards of ISO9919
Measurement range	0 to 100%
Resolution	1%
	70 to 100%: ±2% (adult/pediatric)
Accuracy	70 to 100%: ±3% (neonate)
	0% to 69%: Not specified.
*: When the SpO ₂ sensor is applied for neonatal patients as indicated, the specified accuracy range is increased by $\pm 1\%$, to	

^{*:} When the SpO₂ sensor is applied for neonatal patients as indicated, the specified accuracy range is increased by $\pm 1\%$, to compensate for the theoretical effect on oximeter measurements of fetal hemoglobin in neonatal blood.

A.7.4 PR

Alarm limit	Range (bpm)	Step (bpm)
PR High	(low limit +2) to 300	1
PR Low	15 to (high limit-2)	

PR from Mindray SpO₂ Module

Measurement range	20 to 254 bpm
Resolution	1 bpm
Accuracy	±3 bpm (measured without motion)
Accuracy	±5 bpm (measured with motion)
Refreshing rate	1 s
	7 s (when sensitivity is set to High)
SpO₂ averaging time	9 s (when sensitivity is set to Medium)
	11 s (when sensitivity is set to Low)

PR from Masimo SpO₂ Module

Measurement range	25 to 240 bpm
Resolution	1 bpm
Accuracy	±3 bpm (measured without motion)
Accuracy	±5 bpm (measured with motion)
Refreshing rate	1 s
SPO ₂ averaging time	2-4 s, 4-6 s, 8 s, 10 s, 12 s, 14 s, 16 s
Laurantinian annditiona	Pulse amplitude: >0.02%
Low perfusion conditions	Light penetration: >5%
Low perfusion PR accuracy	±3 bpm

PR from Nellcor SpO₂ Module

Measurement range	20 to 300 bpm
Resolution	1 bpm
A	20 to 250 bpm: ±3 bpm
Accuracy	251 to 300 bpm, not specified
Refreshing rate	1 s

PR from NIBP Module

Measurement range	40 to 240 bpm
Resolution	1 bpm
Accuracy	±3bpm or ±3%, whichever is greater

PR from IBP Module

Measurement range	25 to 350 bpm
Resolution	1 bpm
Accuracy	±1 bpm or ±1%, whichever is greater
Refreshing rate	1 s

A.7.5 NIBP

	Meet standards of EN60601-2-30/IEC60601-2-30, EN1060-1, EN1060-3, SP10 and			
Standards	EN1060-4			
Technique	Oscillometry			
Mode of operation	Manual, Auto and ST	AT		
Auto mode repetition intervals	1, 2, 2.5, 3, 5, 10, 15, 2	20, 30, 60, 90, 120, 18	0, 240 or 480 min	
STAT mode cycle time	5 min			
	Adult, pediatric:	180 s		
Max measurement time	Neonate:	90 s		
		Adult	Pediatric	Neonate
Measurement ranges	Systolic:	40 to 270	40 to 200	40 to 135
(mmHg)	Diastolic:	10 to 210	10 to 150	10 to 100
	Mean:	20 to 230	20 to 165	20 to 110
A a a	Max mean error: ±5 mmHg			
Accuracy	Max standard deviation: 8 mmHg			
Resolution	1mmHg			
Initial cuff inflation procesure range	Adult:	80 to 280		
Initial cuff inflation pressure range (mmHg)	Pediatric:	Pediatric: 80 to 210		
(IIIIII Ig)	Neonate:	ate: 60 to 140		
Default initial cuff inflation pressure	Adult:	160		
(mmHg)	Pediatric:	Pediatric: 140		
(IIIIIIII)	Neonate: 90			
	Adult:	297±3 mmHg		
Software overpressure protection	Pediatric:	240±3 mmHg		
	Neonate:	147±3 mmHg		

Alarm limit	Range (mmHg)	Step (mmHg)
	Adult: (low limit+5) to 270	
Sys High	Pediatric: (low limit+5) to 200	
	Neonate: (low limit+5) to 135	
Sys Low	40 to (high limit-5)	
	Adult: (low limit+5) to 230	
Mean High	Pediatric: (low limit+5) to 165	5
	Neonate: (low limit+5) to 110	5
Mean Low	20 to (high limit-5)	
	Adult: (low limit+5) to 210	
Dia High	Pediatric: (low limit+5) to 150	
	Neonate: (low limit+5) to 100	
Dia Low	10 to (high limit-5)	

^{**}Measurement accuracy verification: In adult and pediatric modes, the blood pressure measurements measured with this device are in compliance with the American National Standard for Electronic or Automated Sphymomanometers (ANSI/AAMI SP10-1992) in terms of mean error and stardard deviation by comparing with intra-arterial or auscultatory measurements (depending on the configuration) in a typical patient population. For auscultatory reference, the 5th Korotkoff sound was used to determine the diastolic pressure.

In neonatal mode, the blood pressure measurements measured with this device are in compliance with the American National Standard for Electronic or Automated Sphymomanometers (ANSI/AAMI SP10-1992 and AAMI/ANSI SP10A-1996) in

terms of mean error and stardard deviation by comparing with intra-arterial measurements (depending on the configuration) in a typical patient population.

A.7.6 Temp

<u> </u>	
Standards	Meet standard of EN12470-4
Technique	Thermal resistance
Measurement range	0 to 50 ℃ (32 to 122 °F)
Resolution	0.1 ℃
Accuracy	$\pm 0.1~^{\circ}\mathrm{C}$ or $\pm 0.2~^{\circ}\mathrm{F}$ (without probe)
Refreshing rate	1 s
Minimum time for accurate	Body surface: <100 s
measurement	Body cavity: <80 s

Alarm limit	Range	Step
T1/T2 High	(low limit +1) to 50 $^{\circ}\mathrm{C}$	
T1/T2 High	(low limit +1.8) to 122 °F	
T1/T2 Low	0 to (high limit -1) $^{\circ}\mathrm{C}$	0.1 ℃
11/12 LOW	32 to (high limit -1.8) °F	0.1 °F
TD High	0 to 50 ℃	
TD High	0 to 90 °F	

A.7.7 IBP

Standards	Meet standard of EN60601-2-34/IEC60601-2-34.
Technique	Direct invasive measurement
IBP	
Measurement range	-50 to 300 mmHg
Resolution	1 mmHg
Accuracy	±2% or ±1 mmHg, whichever is greater (without sensor)
Refreshing rate	1 s
Pressure transducer	
Excitement voltage	5 VDC, ±2%
Sensitivity	5 μV/V/mmHg
Impedance range	300 to 3000Ω
Volume displacement (ABBOTT)	<0.04 mm³ /100 mmHg

Alarm limit		Range (mmHg)	Step (mmHg)
ART	Sys High		1
Ao	Mean High	(low limit + 2) to 300	
FAP BAP	Dia High		
UAP	Sys Low	0 to (high limit – 2)	
LV	Mean Low		

Alarm limit		Range (mmHg)	Step (mmHg)
	Dia Low		
PA	Sys High		
	Mean High	(low limit + 2) to 120	
	Dia High		1
	Sys Low		
	Mean Low	-6 to (high limit – 2)	
	Dia Low		
CVP, LAP	Mean High	(low limit + 2) to 40	1
RAP, ICP	Mean Low	-10 to (high limit – 2)	
P1 to P4	Sys High		
	Mean High	(low limit + 2) to 300	
	Dia High		1
	Sys Low		
	Mean Low	-50 to (high limit – 2)	
	Dia Low		

A.7.8 C.O.

Measurement method	Thermodilution method	
	C.O.:	0.1 to 20 L/min
Measurement range	TB:	23 to 43 ℃
	TI:	0 to 27 ℃
Resolution	C.O.:	0.1 L/min
	TB, TI:	0.1 ℃
Accuracy	C.O.:	±5% or ±0.1 L /min, whichever is greater
Accuracy	TB, TI:	$\pm 0.1~^{\circ}\mathrm{C}$ (without sensor)
Repeatability	C.O.:	±2% or ±0.1 L/min, whichever is greater
Alarm range	TB:	23 to 43 °C

Alarm limit	Range	Step
TB High	(low limit + 1) to 43 °C	0.1 ℃
16 rigii	(low limit + 1.8) to 109.4 $^{\circ}\mathrm{F}$	0.1 °F
TB Low	23 to (high limit - 1) $^{\circ}\mathrm{C}$	0.1 1
16 Low	73.4 to (high limit - 1.8) $^{\circ}\mathrm{F}$	

A.7.9 CO₂

Measurement mode	Sidestream, microstream, mainstream
Technique	Infrared absorption

Sidestream CO₂ Module

Standard	Meet standard of ISO 21647		
CO ₂ Measurement range	0 to 99 mmHg		
	0 to 40 mmHg: ±2 mmHg		
Accuracy*	41 to 76 mmHg: ±5% of the reading		
	77 to 99 mmHg: ±10% of the reading		
Accuracy drift	Meet the requirement for measurement accuracy within 6 hours		
Resolution	1 mmHg		
Sample flowrate	70 ml/min, 100 ml/min		
Sample flowrate tolerance	15% or 15 ml/min, whichever is greater.		
Warm-up time	45s, enter the iso accuracy mode		
warm-up time	After 10 min, enters the full accuracy mode,		
	Measured with a neonatal watertrap and a 2.5-meter neonatal sampling line:		
	<4.5 s @ 100 ml/min		
Response time	<5 s @ 70 ml/min		
nesponse time	Measured with a neonatal watertrap and a 2.5-meter adult sampling line:		
	<6 s @ 100 ml/min		
	<7 s @ 70 ml/min		
	Measured with a neonatal watertrap and a 2.5-meter neonatal sampling line:		
	<4 s @ 100 ml/min		
Gas sampling delay time	<4.5 s @ 70 ml/min		
Gus sampling delay time	Measured with a neonatal watertrap and a 2.5-meter adult sampling line:		
	<5.5 s @ 100 ml/min		
	<6.5 s @ 70 ml/min		
awRR measurement range	0 to 120 rpm		
awRR measurement precision	±2 rpm		
Apnea time	10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s		

Effect of interference gases on CO₂ measurements		
Gas	Concentration (%)	Quantitive effect*
N ₂ O	≤60	
Hal	≤4	
Sev	≤5	±1 mmHg
Iso	≤5	
Enf	≤5	
Des	≤15	±2 mmHg
* means an extra error should be added in case of gas interference when CO measurements are performed between		

^{*:} means an extra error should be added in case of gas interference when CO₂ measurements are performed between 0-40mmHg.

Alarm limit	Range		Step
EtCO₂ High	(low limit + 2) to 99 mmHg		
EtCO ₂ Low	1 to (high limit - 2)mmHg		1 mmHg
FiCO₂ High	1 to 99 mmHg		
au DD Lliab	Adult, pediatric: (lov	w limit + 2) to 100 rpm	
awRR High	Neonate: (low limit + 2) to 150 rpm		1 rpm
awRR Low	0 to (high limit - 2) rpm		

^{*} Accuracy applies for the following conditions:

- 1. Measurements begin after the CO₂ module warms up;
- 2. Ambient pressure is from 750 to 760 mmHg, and ambient temperature from 22 to 28°C;
- 3. The measured gas is a dry gas and the balance gas N_2 ;
- 4. Gas sample flow rate is 100 ml/min, respiration rate is no greater than 50 rpm with a fluctuation between ± 3 rpm, and I:E is 1:2.

When the operating temperature (near the module detector) is between 5° C and 25° C or between 50° C and 65° C, or the respiration rate is between 50° rpm and 60° rpm, the measurement accuracy is: $\pm 4^{\circ}$ mmHg (0 to 40° mmHg) or 12° 0 of the reading (41 to 99° mmHg).

Microstream CO₂ Module

Standard	Meet standard of ISO 21647		
CO ₂ Measurement range	0 to 99 mmHg		
A savus sv.*	0 to 38 mmHg:	±2 mmHg	
Accuracy* 39 to 99 mmHg: ±5% of the reading+0.08% of ($\pm 5\%$ of the reading+0.08% of (the reading-38)	
Accuracy drift	Meet the requirement for measurement accuracy within 6 hours		

^{*} This accuracy is applied to respiration rate no greater than 80 rpm. For respiration rate greater than 80 rpm and $EtCO_2$ value greater than 18 mmHg, the accuracy is 4 mmHg or $\pm 12\%$ of the reading, whichever is greater. For respiration rate greater than 60 rpm, the above accuracy can be achieved by using the CapnoLine H Set for Infant/Neonatal. In the presence of interfering gases, the accuracy specification deteriorates by 4% of the above accuracy.

Resolution	1 mmHg		
Sample flow rate	$50_{+15}^{-7.5}$ ml/min		
Initialization time	30 s (typical)		
	2.9 s (typical)		
	(The response time is the sum of	the rise time and the delay t	time when using a
Response time	FilterLine of standard length)		
	Rise time: <190 ms (10% to 90%)		
	Delay time: 2.7 s (typical)		
awRR measurement range	0 to 150 rpm		
	0 to 70 rpm:	±1 rpm	
awRR measurement accuracy	71 to 120 rpm:	±2 rpm	
	121 to 150 rpm:	±3 rpm	
Apnea alarm time	10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s		
Alarm limit	Range		Step

EtCO₂ High	(low limit + 2) to 99 mmHg	
EtCO ₂ Low	1 to (high limit – 2)mmHg	1 mmHg
FiCO ₂ High	1 to 99 mmHg	
awRR High	Adult, pediatric: (low limit + 2) to 100 rpm	
awnn riigii	Neonate: (low limit + 2) to 150 rpm	1 rpm
awRR Low	0 to (high limit – 2) rpm	

Mainstream CO₂ Module

Standard	Meet standard of ISO 21647		
CO ₂ Measurement range	0 to 150 mmHg		
	0 to 40 mmHg:	±2 mmHg	
Accuracy	41 to 70 mmHg:	±5% of the	reading
Accuracy	71 to 100 mmHg: $\pm 8\%$ of the reading		reading
	101 to 150 mmHg:	±10% of th	e reading
Accuracy drift	Meet the requirement for	or measurement accuracy wit	hin 6 hours
Resolution	1 mmHg		
Response time	<60 ms		
awRR measurement range	0 to 150 rpm		
awRR measurement accuracy	$\pm 1 \text{rpm}$		
Apnea alarm time	10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s		
Alarm limit	Range Step		Step
EtCO ₂ High	(low limit + 2) to 99 mmHg		
EtCO ₂ Low	1 to (high limit - 2)mmHg 1 mmHg		1 mmHg
FiCO ₂ High	1 to 99 mmHg		
awDD High	Adult, pediatric:	(low limit + 2) to 100 rpm	
awRR High	Neonate:	(low limit + 2) to 150 rpm	1 rpm
awRR Low	0 to (high limit - 2) rpm		

A.7.10 AG (For iPM 12/ iPM 10 only)

Standards	Meet standard of ISO 21647		
Technique	Infrared absorption		
	Iso accuracy mode:	45 s	
Warm-up time	Full accuracy mode:	10 min	
	Adult, pediatric:	120, 150, 200 ml/min	
Sample flow rate	Neonate:	70, 90, 120 ml/min	
	Accuracy:	± 10 ml/min or $\pm 10\%$, whichever is greater	
	CO ₂ :	0 to 30%	
	O ₂ :	0 to 100%	
	N₂O:	0 to 100%	
Measurement range	Des:	0 to 30%	
	Sev:	0 to 30%	
	Enf:	0 to 30%	
	Iso:	0 to 30%	

	Hal:	0 to 30%		
	awRR:	2 to 100 rpm		
	CO ₂ :	1 mmHg		
Resolution	awRR:	1 rpm		
	CO ₂ :	CO ₂ full accuracy ± 0.3% _{ABS}		
Iso accuracy	N ₂ O:	± (8% _{REL} +2% _{ABS})		
	Other anesthetic gases:	Full accuracy of other anesthe	etic gases ± 8% _{REL}	
	Gases	Range (% _{REL})	Accuracy (% _{ABS})	
		0 to 1	±0.1	
		1 to 5	±0.2	
	CO ₂	5 to 7	±0.3	
		7 to 10	±0.5	
		>10	Not specified	
	N O	0 to 20	±2	
	N ₂ O	20 to 100	±3	
		0 to 25	±1	
	O ₂	25 to 80	±2	
Full accuracy		80 to 100	±3	
		0 to 1	±0.15	
		1 to 5	±0.2	
	Des	5 to 10	±0.4	
	l Des	10 to 15	±0.6	
		15 to 18	±1	
		>18	Not specified	
		0 to 1	±0.15	
	Sev	1 to 5	±0.2	
		5 to 8	±0.4	
		>8	Not specified	
		0 to 1	±0.15	
	Enf, Iso, Hal	1 to 5	±0.2	
		>5	Not specified	
	awRR	2 to 60 rpm	±1 rpm	
		>60 rpm	Not specified	
Accuracy drift	-	easurement accuracy within 6 h	ours	
Apnea alarm time	10 s, 15 s, 20 s, 25 s, 30 s, 35 s	, 40 s		
Refreshing rate	1 s		<u> </u>	
Rise time		min, using the DRYLINE™ water	trap and neonatal	
(10 % ~ 90%)	DRYLINE™ sampling line (2.5			
	CO ₂	≤250 ms (fall time: ≤200ms)		
	N ₂ O	≤250 ms		
	O ₂	≤600 ms		
	Hal, Iso, Sev, Des	≤300 ms		
	Enf	≤350 ms		
		min, using the DRYLINE™ water	trap and adult DRYLINE™	
	sampling line (2.5m):			

	CO ₂	≤250 ms (fall time: ≤200 ms)	
	N ₂ O	≤250 ms	
	O ₂	≤500 ms	
	Hal, Iso, Sev, Des	≤300 ms	
	Enf	≤350 ms	
Delay time	<4 s		
	Primary anesthetic agent		
	In full accuracy mode: 0.15%,		
	In ISO accuracy mode: 0.4%		
Anesthetic agent limit	Second anesthetic agent:		
	In full accuracy mode: 0.3% or 5% REL (10% in ISO accuracy mode) of primary agent if		
	primary agent is greater than 10%		
	In ISO accuracy mode: 0.5%		

Effect of interference gases on AG measurements					
Gas	C(0/)	Quantitive effect(%ABS) ³⁾			
Gas	Concentration(%)	CO ₂	N ₂ O	Agent 1)	O ₂
CO ₂	/	/	0.1	0	0.2
N ₂ O	/	0.1	/	0.1	0.2
Agent 1) 2)	/	0.1	0.15)	0.14)	1
Xenon	<100%	0.1	0	0	0.5
Helium	<50%	0.1	0	0	0.5
Ethanol	<0.1%	0	0	0	0.5
Acetone	<1%	0.1	0.1	0	0.5
Methane	<1%	0.1	0.1	0	0.5
Saturated Isopropanol vapour	/	0.1	0	0	0.5
Metered dose inhaler propellants,	/	Unspecified	Unspecified	Unspecified	0.5

- 1) Agent represents one of Des, Iso, Enf, Sev, and Hal.
- 2) Multiple agent interference on CO_2 , N_2O and O_2 is typically the same as single agent interference.
- 3) For CO_2 , N_2O and Agents, maximum interference from each gas at concentrations within specified accuracy ranges for each gas. The total interference of all gases is never larger than $5\%_{REL}$.
- 4) Applicable to type A AG module, representing the interference effect of secondary anesthetic agents on primary anesthetic agent.
- 5) Measurement interference to type M AG module originates from the applied anesthetic agent that is configured manually.

Alarm limit	Range	Range	
EtCO ₂ High	(low limit + 2) to 99 mr	(low limit + 2) to 99 mmHg	
EtCO ₂ Low	1to (high limit - 2)mmHg		1 mmHg
FiCO ₂ High	1 to 99 mmHg		
awRR High	Adult, pediatric: Neonate:	(low limit + 2) to 100 rpm (low limit + 2) to 150 rpm	1 rpm
awRR Low	0 to (high limit - 2)rpm		
EtO ₂ High	(low limit + 2) to 100 %		
EtO ₂ Low	18 to (high limit - 2)%		0.1%
FiO ₂ High	(low limit + 2) to 100 %		0.1%
FiO ₂ Low	18 to (high limit - 2)%		
EtN₂O High	(low limit + 2) to 100 %		
EtN ₂ O Low	0 to (high limit - 2)%		1%
FiN₂O High	(low limit + 2) to 100 %		170
FiN₂O Low	0 to (high limit - 2)%		
EtHal/Enf/Iso High	(low limit + 0.2) to 5.0 %		
EtHal/Enf/Iso Low	0 to (high limit - 0.2)%		0.1%
FiHal/Enf/Iso High	(low limit + 0.2) to 5.0 %		0.170
FiHal/Enf/Iso Low	0 to (high limit - 0.2)%		
EtSev High	(low limit + 0.2) to 8.0 9	%	
EtSev Low	0 to (high limit - 0.2)%		0.1%
FiSev High	(low limit + 0.2) to 8.0 9	%	0.1%
FiSev Low	0 to (high limit - 0.2)%		
EtDes High	(low limit + 0.2) to 18.0	%	
EtDes Low	0 to (high limit - 0.2)%		0.1%
FiDes High	(low limit + 0.2) to 18.0	%	0.170
FiDes Low	0 to (high limit - 0.2)%		

B EMC And Radio Regulatory Compliance

B.1 EMC

The device meets the requirements of IEC 60601-1-2.

Note

- Using accessories, transducers and cables other than those specified may result in increased electromagnetic emission or decreased electromagnetic immunity of the patient monitoring equipment.
- The device or its components should not be used adjacent to or stacked with other equipment. If adjacent
 or stacked use is necessary, the device or its components should be observed to verify normal operation in
 the configuration in which it will be used.
- The device needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.
- Other devices may interfere with this monitor even though they meet the requirements of CISPR.
- When the inputted signal is below the minimum amplitude provided in technical specifications, erroneous measurements could result.
- Portable and mobile communication equipment may affect the performance of this monitor.

Guidance and Declaration - Electromagnetic Emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emission tests	Compliance	Electromagnetic environment - guidance
Radio frequency (RF) emissions	Group 1	The device uses RF energy only for its internal function. Therefore, its
CISPR 11		RF emissions are very low and are not likely to cause any interference
		in nearby electronic equipment.
RF emissions CISPR 11	Class A	The device is suitable for use in all establishments other than
Harmonic emissions	Class A	domestic and those directly connected to the public low-voltage
IEC61000-3-2		power supply network that supplies buildings used for domestic
Voltage Fluctuations/Flicker	Complies	purposes, provided the following warning is heeled.
Emissions IEC 61000-3-3		



! WARNING

This equipment/system is intended for use by healthcare professionals only. This equipment/ system may
cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take
mitigation measures, such as re-orienting or relocating the [ME EQUIPMENT or ME SYSTEM] or shielding
the location.

Guidance and Declaration - Electromagnetic Immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment -
minumity test	iEC00001 test level	Compnance level	guidance
Electrostatic discharge	±6 kV contact	±6 kV contact	Floors should be wood, concrete or
(ESD) IEC 61000-4-2	±8 kV air	±8 kV air	ceramic tile. If floors are covered
			with synthetic material, the relative
			humidity should be at least 30%.
Electrical fast	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that
transient/burst IEC	±1 kV for input/output lines	±1 kV for input/output lines	of a typical commercial or hospital
61000-4-4	(>3 m)	(>3 m)	environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s)	±1 kV line(s) to line(s)	
	±2 kV line(s) to earth	±2 kV line(s) to earth	
Voltage dips, short	<5 % U _T (>95 % dip in U _T) for	<5 % U _T (>95 % dip in U _T) for	Mains power quality should be that
interruptions and	0.5 cycle	0.5 cycle	of a typical commercial or hospital
voltage variations on			environment. If the user of our
power supply input	40 % U _T (60 % dip in U _T) for 5	40 % U _T (60 % dip in U _T) for 5	product requires continued
lines IEC 61000-4-11	cycles	cycles	operation during power mains
			interruptions, it is recommended
	70 % U _T (30 % dip in U _T) for	70 % U _T (30 % dip in U _T) for	that our product be powered from
	25 cycles	25 cycles	an uninterruptible power supply or
			a battery.
	<5 % U _T (>95 % dip in U _T) for	<5 % U _T (>95 % dip in U _T) for	
	5 s	5 s	
Power frequency	3 A/m	3 A/m	Power frequency magnetic fields
(50/60 HZ) magnetic			should be at levels characteristic of
field IEC 61000-4-8			a typical location in a typical
			commercial or hospital
			environment.
Note: U _T is the AC mains voltage prior to application of the test level.			

Guidance and Declaration - Electromagnetic Immunity

The device is intended for use in the specified electromagnetic environment. The customer or the user of the device should assure that it is used in such an environment as described below.

Immunity test	IEC60601 test	Compliance	Electromagnetic environment - guidance
	level	level	
Conduced RF	3 Vrms	3Vrms	Portable and mobile RF communications equipment should be
IEC61000-4-6	150 kHz to 80 MHz		used no closer to any part of the system, including cables, than
			the recommended separation distance calculated from the
			equation appropriate for the frequency of the transmitter.
			Recommended separation distances:
			$d = 1.2\sqrt{P}$
Radiated RF	3V/m	3V/m	Recommended separation distances:
IEC61000-4-3	80MHz to 2.5GHz		80 MHz∼800 MHz
			$d = 1.2\sqrt{P}$
			800MHz-2.5GHz
			$d = 2.3\sqrt{P}$
			Where, P is the maximum output power rating of the
			transmitter in watts (W) according to the transmitter
			manufacturer and d is the recommended separation distance
			in meters (m).b
			Field strengths from fixed RF transmitters, as determined by an
			electromagnetic site survey c, should be less than the
			compliance level in each frequency range d.
			Interference may occur in the vicinity of
			equipment marked with the following
			symbol: ♣

Note 1: At 80 MHz to 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the [ME EQUIPMENT or ME SYSTEM] is used exceeds the applicable RF compliance level above, the [ME EQUIPMENT or ME SYSTEM] should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the [ME EQUIPMENT or ME SYSTEM].

 $_{\mbox{\scriptsize b}}$ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.



∆ WARNING

• The monitor is configured with a wireless network connector to receive wireless signal. Other devices may interfere with this monitor even though they meet the requirements of CISPR.

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum	Separation distance in meters (m) according to frequency of the transmitter		
output power of	150 kHz \sim 80 MHz	80 MHz \sim 800 MHz	800 MHz \sim 2.5 GHz
transmitter (W)	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.00

- ◆ For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
- ◆ Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- ◆ **Note 2**: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

B.2 Radio Regulatory Compliance

This Wi-Fi device complies with Canadian ICES-001. Cet appareil ISM est conforme a la norme NMB-001 du Canada.

The patient monitor including Wi-Fi module (US only) complies with part 15 of the FCC Rules. Operation is subject to the condition that this device does not cause harmful interference.

Operation of this equipment requires the prior coordination with a frequency coordinator designated by the FCC for the Wireless Medical Telemetry Service.

The patient monitor including Wi-Fi module - FCC and Industry Canada Radio Compliance: This device complies with Part 15 of the FCC Rules and RSS-210 of Industry Canada. Operation is subject to the following two conditions: (1) this device may not cause harmful

interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications to this equipment not expressly approved by Mindray may cause harmful radio frequency interference and void your authority to operate this equipment.

The maximum antenna gain permitted complies with the e.i.r.p. limits as stated in RSS-210.

The maximum antenna gain permitted complies with the e.i.r.p. limits specified for point-to-point operation, as stated in RSS-210.



The radio device used in this product is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC (Radio Equipment and Telecommunications Terminal Equipment Directive). The product is in compliance with ETSI EN 300 328 and ETSI EN 301 489.



WARNING

Keep a distance of at least 20cm away from the monitor when Wi-Fi function is in use.

FOR YOUR NOTES		

C Default Configurations

This chapter lists some of the most important factory default settings for each department in configuration management. You cannot change the factory default configuration itself. However, you can make changes to the settings from the factory default configuration and then save the changed configuration as a user configuration. The last column of the following tables is for your notes and review.

Note: In this chapter, O.M means the monitor's operating mode. Column C refers to the settings that can be changed in configuration management. Column M refers to the settings that can be changed in monitoring mode.

C.1 Parameters Configuration

C.1.1 ECG

ECG Setup

L. N		0.1	1		0.0	ICII	NICH	CCII	
Item Name		C	М	General	OR	ICU	NICU	CCU	User Defaults
Lead Set			*	Auto					
Alm Source		*	*	HR					
Alarm		*	*	On					
Alm Lev		*	*	Med					
	Adu			120					
HR/PR High	Ped	*	*	160					
	Neo			200					
	Adu			50					
HR/PR Low	Ped	*	*	75					
	Neo			100					
Sweep		*	*	25 mm/s					
Beat Vol		*	*	2		1			
Paced			*	No					
Notch Filter		*	*	On					
Gain		*	*	X1					
Filter		*	*	Monitor	Surgery	Monitor		Diagnostic	
ECG Display		*	*	Normal	•				
Pacemaker R	ate		*	60					

ST Analysis

Item Name	0.1	N	General	OR	ICU	NICU	CCU	User Defaults
item Name	С	M	General	OK	ico	NICO	cco	Oser Delauits
ST Analysis	*	*	Off				On	
Alarm	*	*	Off					
Alm Lev	*	*	Med					
ST-X High	*	*	when ST Un	it is mV:		0.20		
31-X High			when ST Un	it is mm:		2.0		
ST-X Low	*	*	when ST Un	it is mV:		-0.20		
31-X LOW			when ST Un	it is mm:		-2.0		
ISO			-80 ms					
J	*	*	48 ms					
ST			J + 60 ms					

X represents I, II, III, aVR, aVL, aVF, V, V1, V2, V3, V4, V5 or V6.

Arrh. Analysis

		0.1	M						
Item Name	Algorithm	С	М	Genral	OR	ICU	NICU	CCU	User Defaults
Arrhythmia Thresh	old Settings								
PVCs High	Mindray	*	*	Adu, Ped:	10				
T VC3Triigit	lviiiiaiay			Neo:	N/A	1			
				Adu:	120				
Tachy		*	*	Ped:	160				
				Neo:	N/A	ı			
				Adu:	50				
Brady		*	*	Ped:	75				
				Neo:	N/A				
Asys. Delay		*	*	Adu, Ped:	5				
Asys. Delay				Neo:	N/A				
Vtac Rate		*	*	Adu, Ped:	130				
viac Rate				Neo:	N/A	<u>.</u>			
Vtac PVCs		*	*	Adu, Ped:	6				
VIACPVCS		"	,	Neo:	N/A	ı			
Multif. PVC's		*	*	Adu, Ped:	15				
Window		Î	_	Neo:	N/A				
				Adu:	160				
Extreme Tachy		*	*	Ped:	180				
				Neo:	N/A	ı.			
				Adu:	35				
Extreme Brady		*	*	Ped:	50				
				Neo:	N/A				

Vbrd Rate			0.1	И						
Neo: N/A	Item Name	Algorithm	С	м	Genral	OR	ICU	NICU	CCU	User Defaults
Neo: N/A	Vhrd Rate		*	*	Adu, Ped:	40				
Pause Time	Void Nate				Neo:	N/A				
Neo: N/A	Vbrd PVCs		*	*	Adu, Ped:	5				
Pause Time		1				N/A				
PVCs High	Pause Time		*	*						
Neo: N/A										
Asys. Delay	PVCs High		*	*						
* * Neo: N/A			_							
Vtac Rate * * Adu, Ped: N/A 130 Neo: N/A Vtac PVCs * * Adu, Ped: 6 Neo: N/A 6 Neo: N/A Multif. PVC's Window * * Adu, Ped: 15 Neo: N/A 15 Neo: N/A Tachy * * Ped: 160 Neo: N/A 160 Neo: N/A Brady * * Ped: 75 Neo: N/A 75 Neo: N/A Arrhythmia Alarm Settings * * Off On PVCs/min Alarm Mindray Alarm * * Off On Nonsus. Vtac Alarm Vent. Rhythm Alarm Bigeminy Alarm * * Off On Trigeminy Alarm * * Off On	Asys. Delay		*	*						
Vtac PVCs Mortara * * Neo: N/A Multif. PVC's Window Mortara * * Adu, Ped: 6 Neo: N/A Tachy * * Neo: N/A Brady * * Ped: 15 Neo: N/A Adu: 120 Neo: N/A Neo: N/A Adu: 50 Neo: N/A Neo: N/A Arrhythmia Alarm Settings * * Off PVCs/min Alarm Mindray Alarm * * Off R on T Alarm Mindray Nonsus. Vtac Alarm * * Off Vent. Rhythm Alarm Bigeminy Alarm * * Off On Trigeminy Alarm * * Off On			_							
Vtac PVCs	Vtac Rate		*	*						
Window * * Neo: N/A Multif. PVC's Window * Adu, Ped: 15 Neo: N/A Tachy Adu: 120 Neo: N/A Brady * * Ped: 160 Neo: N/A Adu: 50 Ped: 75 Neo: N/A Arrhythmia Alarm Settings PVCs/min Alarm Mindray R on T Alarm * * Off Nonsus. Vtac Alarm * * Off Vent. Rhythm Alarm * * Off Bigeminy Alarm * * Off Trigeminy Alarm On										
Multif. PVC's Window * * Adu, Ped: N/A 15 Neo: N/A Tachy Adu: 120 * * Ped: 160 Neo: N/A Brady * * Ped: 75 Neo: N/A * * Ped: 75 Neo: N/A Arrhythmia Alarm Settings PVCs/min Alarm Mindray * * Off On R on T Alarm * * Off On Nonsus. Vtac Alarm * * Off On Vent. Rhythm Alarm * * Off On Bigeminy Alarm * * Off On Trigeminy Alarm * * Off On	Vtac PVCs		*	*						
Neo: N/A Adu: 120	Multif. PVC's	- Mortara			Adu, Ped:	15				
* * Ped: 160 Neo: N/A	Window		*	*	Neo:	N/A				
Neo: N/A					Adu:	120				
Adu: 50	Tachy		*	*	Ped:	160				
R on T Alarm Mindray * * Off On					Neo:	N/A				
Neo: N/A										
Arrhythmia Alarm Settings PVCs/min Alarm Mindray * * Off On R on T Alarm * * Off On Nonsus. Vtac Alarm * * Off On Vent. Rhythm Alarm * * Off On Bigeminy Alarm * * Off On Trigeminy Alarm On On	Brady		*	*						
PVCs/min Alarm Mindray * * Off On R on T Alarm * * Off On Nonsus. Vtac Alarm * * Off On Vent. Rhythm Alarm * * Off On Bigeminy Alarm * * Off On Trigeminy Alarm On On					Neo:	N/A				
X		Settings			I					
R on T Alarm * * Off On Nonsus. Vtac Alarm * * Off On Vent. Rhythm Alarm * * Off On Bigeminy Alarm * Off On Trigeminy Alarm * Off On		Mindray	*	*	Off				On	
* * Off			_	¥	0#				0	
Vent. Rhythm Alarm * * Off On Bigeminy Alarm * * Off On Trigeminy Alarm Off On On On On			_							
Bigeminy Alarm * * Off On Trigeminy Alarm Off On On		_	_						+	
Trigeminy Alarm * * Off On			<u> </u>						-	
		_	_							
Asystole Alarm * * On									OII	
VFib/VTac					JII					
Alarm * * On			*	*	On					
Vtac Alarm * * On			*	*	On					
Vent. Brady Alarm * * On			*	*						
Extreme Tachy	·									
Alarm * * On			*	*	On					

		0.1	M						
Item Name	Algorithm	С	М	Genral	OR	ICU	NICU	CCU	User Defaults
Extreme Brady		*	*	On					
Alarm			"	On					
X Alarm		*	*	Off					
Asystole Alm Lev		*	*	High					
VFib/VTac		*	*	High					
Alm Lev				riigii					
Vtac Alm Lev		*	*	High					
Vent. Brady Alm Lev		*	*	High					
Extreme Tachy Alm Lev		*	*	High					
Extreme Brady Alm Lev		*	*	High					
VT>2 Alm Lev		*	*	Low					
Pause Alm Lev		*	*	Low					
Couplet Alm Lev		*	*	Prompt					
PVC Alm Lev		*	*	Prompt					
Irr. Rhythm Alm Lev		*	*	Prompt					
PNP Alm Lev		*	*	Prompt					
PNC Alm Lev		*	*	Prompt					
Missed Beats Alm Lev		*	*	Prompt					
X Alm Lev		*	*	Med					
X Alm Rec		*	*	Off					
PVCs/min									
Alarm	Mortara	*	*	Off				On	
R on T Alarm		*	*	Off				On	
Vent. Rhythm Alarm		*	*	Off				On	
Bigeminy Alarm		*	*	Off				On	
Trigeminy Alarm		*	*	Off				On	
Asystole Alarm				On					
VFib Alarm				On					
VTac Alarm	*			On					
X Alarm		*	*	On					
Asystole Alm Lev		*	*	High					
VFib Alm Lev		*	*	High					
VTac Alm Lev		*	*	High					
VT>2 Alm Lev		*	*	Low					

Item Name	Algorithm	0.1	И	Genral	OR	ICU	NICU	ccu	User Defaults
item Name	Aigoritiiii	c	M	Geniai			Nico		Osei Delauits
Couplet Alm Lev		*	*	Prompt					
PVC Alm Lev		*	*	Prompt					
Irr. Rhythm Alm Lev		*	*	Prompt					
PNP Alm Lev		*	*	Prompt					
PNC Alm Lev		*	*	Prompt					
Missed Beats Alm		*	*	Prompt					
Lev				Trompt					
X Alm Lev		*	*	Med					
X Alm Rec		*	*	Off		·	·		

X represents a certain arrhythmia event. Refer to the Specifications chapter for details. The X in "X Alm Lev" refers to all arrhythmia events except for those specially marked ones.

C.1.2 RESP

Item Name	0.1	Λ	General	OR	ICU	NICU	ccu	User Defaults
item Name	С	М	General	OR		NICO	cco	Oser Delauits
Alarm	*	*	On					
Alm Lev	*	*	Med					
Sweep	*	*	6.25 mm/s					
Lead	*	*	II					
Gain	*	*	X2					
RR High	*	*	Adu, Ped:	3	30			
nn riigii			Neo:	1	00			
RR Low	*	*	Adu, Ped:	8	3			
NN LOW			Neo:	3	30			
Apnea Delay	*	*	Adu, Ped:	2	20			
Aprilea Delay			Neo:	1	5			
Detection Mode	*	*	Auto					
RR Source		*	Auto					

C.1.3 PR

Itom Namo	Item Name			General	OR	ICU	NICU	ccu	User Defaults		
item Name	СМ				OK	ico	NICO	cco	Oser Delauits		
Alarm		*	*	On							
Alm Lev		*	*	Med							
	Adu			120							
HR/PR High	Ped	*	*	160	60						
	Neo			200							
	Adu			50							
HR/PR Low	Ped	*	*	75							
	Neo			100	100						
PR Source		*	*	SpO ₂							
Beat Vol		*	*	2	·	1					

C.1.4 SpO₂

Ham Nama	O.M		General	OB	ICII	NICH	CCII	User Defaults
Item Name	С	М	General	OR	ICU	NICU	CCU	User Detaults
Alarm	*	*	On					
Alm Lev	*	*	Med					
SpO ₂ High	*	*	Adu, Ped:			100		
SpO₂ nigii			Neo:			95		
SpO₂ Low	*	*	90					
Desat Limit	*	*	80					
Sweep	*	*	25 mm/s					
NIBP Simul		*	Off					
Sensivity (Mindray)	*	*	Med					
Sensivity (Masimo)	*	*	Normal					
Averaging (Masimo)	*	*	8 s					
Sat-Seconds (Nellcor)	*	*	0 s					

C.1.5 Temp

Item Name	O.M		General	OR	ICU	NICU	ccu	User Defaults
Item Name	C	М	General		ico	IIICO	CCO	Osei Delauits
Alarm	*	*	On					
Alm Lev	*	*	Med					
T1/T2 High (℃)	*	*	38.0					
T1/T2 Low(°C)	*	*	35.0					
TD High (°C)	*	*	2.0					

C.1.6 NIBP

Item Name		0.	М	General	OR	ICU	NICU	ccu	User Defaults
item Name		С	М	General	JOK J		Nico		Osei Deiaults
Alarm			*	On					
Alm Lev		*	*	Med					
Interval		*	*	15 min	5 min	15 min	30 min	15 min	
Adu				80					
Cuff Press. (mmHg)	Ped	*	*	60					
(mining)	Neo			40					
Initial Drossura	Adu			160					
Initial Pressure Ped			*	140					
(mmHg) Neo			90						

Alarm Limits					
NIDD C I I: «h	Adu			160	
NIBP-S High	Ped	*	*	120	
(mmHg)	Neo			90	
NIBP-S Low	Adu			90	
(mmHg)	Ped	*	*	70	
(IIIIIIIIIII)	Neo			40	
NIBP-M High	Adu			110	
(mmHg)	Ped	*	*	90	
(IIIIIIIIIII)	Neo			70	
NIBP-M Low	Adu			60	
(mmHg)	Ped	*	*	50	
(mining)	Neo			25	
NIBP-D High	Adu			90	
(mmHg)	Ped	*	*	70	
(IIIIIIIIII)	Neo			60	
NIBP-D Low	Adu			50	
	Ped	*	*	40	
(mmHg)	Neo			20	

C.1.7 IBP

C.1./ IBP			_											
Item Name		O.N	и М	General	OR	ICU	NICU	сси	User Defaults					
Alarm		*	*	On	On									
Alm Lev		*	*	Med										
P1 Measure		*	*	All										
P2 Measure		*	*	All										
P3 Measure		*	*	Mean										
P4 Measure		*	*	Mean	ean									
Sensitivity		*	*	Med										
Sweep		*	*	25 mm/s										
Filter		*		12.5 Hz										
Art, Ao, UAP, BAP	, FAP, LV,	P1-	P2 A	rterial Press	ure Alarm	Limits								
IBP-S High	Adu			160										
(mmHg)	Ped	*	*	120										
(пшну)	Neo			90										
IBP-S Low	Adu			90										
(mmHg)	Ped	*	*	70										
(IIIIIIII)	Neo			55										
IBP-M High	Adu			110										
(mmHg)	Ped	*	*	90										
(illilling)	Neo			70										
IBP-M Low	Adu	*	*	70										

Item Name		0.1		General	OR	ICU	NICU	ccu	User Defaults
		C	M		ON	ico	NICO	CCO	Osei Delaults
(mmHg)	Ped			50					
	Neo			35					
IBP-D High	Adu			90					
(mmHg)	Ped	*	*	70					
	Neo			60					
IBP-D Low	Adu			50					
(mmHg)	Ped	*	*	40					
	Neo			20					
PA Alarm Limits	A also			25					
PA-S High	Adu	*	*	35					
(mmHg)	Ped		^	60					
	Neo			60					
DA CL (II)	Adu	*	*	10					
PA-S Low(mmHg)	Ped	*	*	24					
	Neo			24					
PA-M High	Adu	*	*	20					
(mmHg)	Ped	*	*	26					
	Neo			26					
PA-M Low	Adu	*	*	0					
(mmHg)	Ped	*	*	12					
	Neo			12					
PA-D High	Adu	*	*	16					
(mmHg)	Ped			4					
	Neo			4					
PA-D Low	Adu	*	*	0					
(mmHg)	Ped	*	*	-4					
	Neo			-4		•.			
CVP, LAP, RAP, ICP		3-P4 	Ven	I	Alarm Lin	nits			
IBP-M High	Adu	*	*	10					
(mmHg)	Ped	*	*	4					
	Neo			4					
IBP-M Low	Adu	*	*	0					
(mmHg)	Ped		^	0					
A A - DAD FAD	Neo	2 0	•• -	0	-1-				
Art, Ao, BAP, FAP,	LV, P1-P	2 Ar	teria 	T	.ale				
Scale (mmHg)		L^	<u> </u>	0-160					
PA Scale			*	0-30					
Scale (mmHg)	11\/D.C-	*	1	0-30					
CVP, LAP, RAP, ICP	, 007 50	ale *	*	0.20					
Scale (mmHg)	e Drocess			0-20					
UAP, P3-P4 Venou	s Pressul	re So	:aie *	0.90					
Scale (mmHg)		, ·		0-80					

C.1.8 C.O.

Item Name	O.M		General	OR	ICU	NICU	ccu	User Defaults		
item Name	С	М	General	OK	ico	NICO	CCO	Oser Delauits		
Alarm	*	*	On	On						
Alm Lev	*	*	Med							
TB High (°C)	*	*	39.0	39.0						
TB Low (°C)	*	*	36.0	36.0						
Comp. Const	*	*	0.542	D.542						
Auto TI	*	*	Auto	Auto						
Manual TI(°C)	*	*	2.0	2.0						
Measuring mode	*	*	Manual							

C.1.9 CO₂

C.1.9 CO ₂	0.1	Λ									
Item Name	C	М	Ge	nera	ı	OR	ICU	NICU	ccu	User Defaults	
Alarm	*	*	On								
Alm Lev	*	*	Me	d							
Operating Mode	*	*	Me	asure	2						
Sweep	*	*	6.2	5 mm	n/s						
Scale (mmHg)	*	*	50								
	*	*	Adı	u, Pe	d:						
Apnea Delay	*	*	Ne	o:			15				
RR Source		*	Aut	to							
Sidestream CO ₂ Setup											
Flow Rate	*	Adı	Adu, Ped: 100 ml/min								
FIOW Rate	*	"	Ne	o:			70	ml/min			
BTPS Compen	*	*	Off								
N₂O Compen	*	*	0								
O ₂ Compen	*	*	21			50	21				
Des Compen	*	*	0								
Microstream CO ₂ Setu	цр										
BTPS Compen				*	Off						
Max Hold			*	*	20 s						
Auto Standby (min)			*	*	0						
Mainstream CO₂ Setup											
Max Hold			*	*	10 s						
O ₂ Compen			*	*	Off						
Balance Gas			*	*	Room Air	•					
AG Compen			*	*	0						

Alarm Limits					
EtCO High (mmHg)	*	*	Adu, Ped:	50	
EtCO ₂ High (mmHg)			Neo:	45	
EtCO ₂ Low (mmHg)	*	*	Adu, Ped:	25	
EtCO2 LOW (IIIIIIIII)			Neo:	30	
FiCO ₂ High (mmHg)	*	*	Adu, Ped, Neo:	4	
RR High	*	*	Adu, Ped:	30	
Mittigii			Neo:	100	
RR Low	*	*	Adu, Ped:	8	
INI LOW			Neo:	30	

C.1.10 AG

Itom Nama	0.1	И	Comoval	OB	ICU	NICU	ccu	Heav Deferrite			
Item Name	С	М	General	OR	ICU	NICU	CCU	User Defaults			
Alarm	*	*	On			•	•				
Alm Lev	*	*	Med	Med							
Sweep	*	*	6.25 mm/s								
O ₂ Compen	*	*	Off	On	Off						
Operating Mode	*	*	Measure		•						
Flow Rate	*	*	Adu, Ped:	120 ml/ı	min						
riow kate			Neo:	70 ml/m	iin						
Auto Standby	*	*	Off								
Apnea Time	*	*	20 s								
RR Source		*	Auto								
CO₂ Setup											
Wave Type	*	*	Draw								
Caala	*	*	when unit is mmHg:		50						
Scale	Î		when unit is % or k	(Pa:	7.0						
TACO Himb (manallar)	*	*	Adu, Ped:	50							
EtCO ₂ High (mmHg)			Neo:	45							
T#CO Lavy (mamalla)	*	*	Adu, Ped:	25							
EtCO ₂ Low (mmHg)			Neo:	30							
FiCO ₂ High(mmHg)	*	*	4								
RR High	*	*	Adu, Ped:	30							
Mittigii			Neo:	100							
RR Low	*	*	Adu, Ped:	8							
NN LOW			Neo:	30							
Gas Setup											
Agent	*	*	AA								
N₂O Scale	*	*	50								
O ₂ Scale	*	*	when unit is mmHg: 400								
O ₂ Scale			when unit is % or k	(Pa:	50						
AA Scale	*	*	9.0								
Hal/Enf/Iso Scale	*	*	2.5								

Itam Nama	٥.٨	1	CI	OD	ICII	NICH	CCII	Han Dafaulta
Item Name	С	М	General	OR	ICU	NICU	CCU	User Defaults
Des Scale	*	*	9.0	•	•	•	•	
Sev Scale	*	*	4.0					
EtO₂ High	*	*	88					
EtO ₂ Low	*	*	18					
FiO₂ High	*	*	Adu, Ped: Neo:	100 90				
FiO ₂ Low	*	*	18					
EtN₂O High	*	*	55					
EtN ₂ O Low	*	*	0					
FiN₂O High	*	*	53	53				
FiN ₂ O Low	*	*	0	0				
EtHal/Enf/Iso High	*	*	3.0					
EtHal/Enf/Iso Low	*	*	0.0					
FiHal/Enf/Iso High	*	*	2.0					
FiHal/Enf/Iso Low	*	*	0.0					
EtSev High	*	*	6.0					
EtSev Low	*	*	0.0					
FiSev High	*	*	5.0					
FiSev Low	*	*	0.0					
EtDes High	*	*	8.0					
EtDes Low	*	*	0.0					
FiDes High	*	*	6.0					
FiDes Low	*	*	0.0					

C.2 Routine Configuration

C.2.1 Alarm

Item Name	O.M		General	OR	ICU	NICU	ccu	User Defaults
item Name	C	M	General	OK	ico	NICO	CCO	Oser Delauits
Alm Volume	*	*	2	1	2			
Reminder Vol	*	*	Low	_ow				
Recording Length	*	*	16 s					
Annes Delay	*	*	Adu, Ped: 20 s					
Apnea Delay			Neo: 15 s					
Alarm Delay	*	*	6 s					
ST Alarm Delay	*	*	30 s					

C.2.2 Screens

Item Name		0.1	И	General	OR	ICU	NICU	CCU	User	
item Name		C	M	General	OK	lco	NICO		Defaults	
Choose Screen		*	*	Normal Screen						
Display the ST segments on ECG screen			*	Unselected						
Select QuickKeys		*		NIBP Measure-	→Stop A	ll→Revi	ew→Standb	у		
	1			ECG1						
	2			ECG2						
Calast Move Converse	3	*		SpO ₂ +PR	SpO₂+PR					
Select Wave Sequence for Normal Screen	4		*	Any IBP						
Tor Normal Screen	5			Any IBP						
	6			CO ₂						
	7			Resp						
	Parameter 1			ECG						
Select Parameters for	Parameter 2	*	*	SpO ₂ +PR						
Big Numerics Screen	Parameter 3	Î		Resp						
	Parameter 4			NIBP						

C.2.3 Waveform

Itama Nama		O.N	1	Camaral	OD	ICII	NICH	CCII	Hear Defends												
Item Name		С	М	General	OR	ICU	NICU	ccu	User Defaults												
Parameter/	ECG		*	Green																	
Wave Colour	NIBP			White																	
	SpO ₂			Cyan																	
	PR			Cyan																	
	TEMP			White																	
	Art/Ao/UAP/FAP																				
	/BAP/LV/P1~P4			Red																	
	(arterial pressure)																				
	PA			Yellow																	
	CVP/ICP/P1~P4			Blue																	
	(venous pressure)			Blue																	
	LAP			Purple																	
	RAP															Orange					
	UVP			Cyan																	
	CO ₂			Yellow																	
	RESP			Yellow																	
	AA			Yellow																	
	N ₂ O			Blue																	
	O ₂			Green																	
	Hal			Red																	
	Enf		Orange																		
	Iso			Purple																	
	Des			Cyan																	

	Sev		Yellow	
	C.O.		White	

X represents a waveform label, such as ECG, RESP, CO₂ and so forth. The ECG waveform cannot be set off.

C.2.4 Review

Item Name		O.M	ı	General	OR	ICU	NICU	ccu	User Defaults
		С	М	General	ON		NICO	cco	Oser Delauits
Tabular Trends	Interval	*	*	30 min	5min	30 min			
Tabular Herius	Trend Group	*	*	Standard					
Graphic Trends	Trend Group	*	*	Standard					
Minitrend Length	1		*	2 h					
Full Disclosure	Save Waves	*	*	Save ECG1 by d	efault.				

C.2.5 Event

Item Name	O.M		General	OR	ICU	NICU	ccu	User Defaults		
item Name	c	М	General	OK	ico	NICO		Oser Delauits		
Waveform 1		*	II	l						
Waveform 2		*	I	l Pleth I						
Waveform 3		*	Pleth	Pleth Resp Pleth						

C.2.6 Record

Item Name		O.M		General	OR	ICU	NICU	ccu	User Defaults	
		c	М	General			luco	-	Oser Delauits	
Length			*	8 s	3 s					
Interval			*	Off	Off					
Paper Speed			*	25 mm/s						
Alm Rec	Х		*	Off						

X represents a parameter label.

C.2.7 Print

I4 NI		0.1	М	General	00	ICU	NICU	CCII	User Defaults	
Item Name		С	М	General	OR	ICO	NICO	CCU	Oser Defaults	
Paper Size			*	A4						
	Amplitude		*	10 mm/mV						
ECG Reports	Sweep		*	25 mm/s						
LCG Reports	Auto Interval		*	Off						
	12-Lead Format		*	12X1						
	Set as End Case Report		*	Unselected						
	Back		*	Auto						
	Spacing		*	Auto						
Tabular Trends	Report Layout		*	Parameter 0	Oriente	d				
Reports	Currently Displayed Trended Parameters	d		Selected						
	Standard Parameter Group		*	Unselected						
	Custom		*	Unselected						
Cuambia Tuamda	Set as End Case Report		*	Unselected						
Graphic Trends Reports	Back		*	Auto	Auto					
vehous	Zoom		*	Auto						
	Set as End Case Report		*	Unselected						
Realtime Report	Sweep		*	Auto						
	Select Wave		*	Current						

C.2.8 Others

Item Name		O.M		General	OR	ICU	NICU	ccu	User Defaults
		C	М	General	OK		NICO		Oser Delauits
Brightness			*	5					
Key Volume			*	2					
View Other	Auto Alarm		*	On					
Patient	Auto Alami			On					

C.3 User Maintenance Items

Item Name		1	General	OR	ICU	NICU	CCU	User Defaults		
item vanie	C	М	General	OK	ico	illeo	cco	Oser Delauits		
Changing Bed No.		*	Protected	Protected						
Atmospheric Pressure		*	760 mmHg	760 mmHg						
Height Unit		*	cm							
Weight Unit		*	kg							
ST Unit		*	mV	mV						
Press. Unit		*	mmHg							
CVP Unit		*	cmH₂O	cmH₂O						

Item Name		0.1	1	General	OR	ICU	NICU	ccu	User Defaults		
item Name		С	М	General	OK		NICO		Oser Delauits		
CO ₂ Unit			*	mmHg	mmHg						
O ₂ Unit			*	%							
Temp Unit			*	$^{\circ}$ C	$^{\circ}$						
Network Type			*	LAN	LAN						
Latching Alarms		*	*	No							
Alarm Pause Time		*	*	2 min							
Minimum Alarm Vo	lume	*	*	2	1	2					
Reminder Tone			*	Off							
Reminder Interval			*	1 min							
ECGLeadOff Lev.			*	Low							
SpO₂SensorOff Lev.			*	Low							
Alaysa Taya Justawa	•		*	High Level Alar	m:	•	10 s				
Alarm Tone Interval				Med/Low Level	l Alarm:	2	20 s				
Lethal Arrh. OFF			*	Disable							
Silence Other Bed			*	On							
E xtended Arrh.			*	Enable							
Wave Line			*	Mediate							
ECG Standard			*	АНА							
Notch Freq.			*	50 Hz							
Data Transfer Metho	od		*	Module							
Transferred Data Le	ngth		*	4 h							
Data Transfer Strate	egy		*	Always Ask							
SpO ₂ Tone			*	Mode 1							
Signal Type		*	Continuous								
Contact Type			*	Normally Close	d						
	Signal Type		**	Continuous							
Nivers Call	Contact Type		*	Normally Close	d						
Nurse Call	Alm Lev	*	*	High, Med, Low							
	Alarm Cat.	*	*	Phys., Tech.							

D Alarm Messages

This chapter lists only the most important physiological and technical alarm messages. Some messages appearing on your monitor may not be included.

In this chapter:

The "I" field indicates how alarm indications are cleared: "A" means all alarm indications are cleared after the A hardkey is pressed, "B" indicates alarm lamp flashing and alarm tones are cleared and the alarm messages change to prompt messages after the A hardkey is pressed, and "C" indicates alarm lamp flashing and alarm tones are cleared and \checkmark appears before the alarm message after the A hardkey is pressed.

The "L" field indicates the alarm level: H means high, M means medium and L means low. "*" means the alarm level is user-adjustable.

XX represents a measurement or parameter label, such as ECG, NIBP, HR, ST-I, PVCs, RR, SpO₂, PR, etc.

In the "Cause and Solution" column, corresponding solutions are given instructing you to troubleshoot problems. If the problem persists, contact your service personnel.

D.1 Physiological Alarm Messages

Measurement	Alarm messages	L	Cause and solution
	XX Too High	M*	XX value has risen above the high alarm limit or fallen below the low
XX	XX Too Low	M*	alarm limit. Check the patient's condition and check if the patient
	AX 100 LOW	IVI	category and alarm limit settings are correct.
ECG	ECG Weak Signal	Н	The ECG signal is so weak that the monitor can't perform ECG
LCG	ECG Weak Signal	"	analysis. Check the patient's condition and the ECG connections.
	Asystole	Н	Arrhythmia has occurred to the patient. Check the patient's condition
	VFib/VTac	Н	and the ECG connections.
	Vtac	Н	
	Vent. Brady	Н	
	Extreme Tachy	Н	
	Extreme Brady	Н	
	RonT	M*	
	VT>2	M*	
	Couplet	M*	_
	PVCs/min	M*]
	Bigeminy	M*]
	Trigeminy	M*]

Measurement	Alarm messages	L	Cause and solution
	Tachy	M*	
	Brady	M*	
	Missed Beats	M*	
	Irr. Rhythm	M*	
	Vent. Rhythm	M*	
	Multif. PVC	M*	
	Nonsus. Vtac	M*	
	Pause	M*	1
	PNP	M*	The massy are according to the second
	PNC	M*	The pacer appears abnormal. Check the pacer.
			The respiration signal was so weak that the monitor cannot perform
	Resp Apnea	Н	respiration analysis. Check the patient's condition and the Resp
Resp			connections.
	Door Autifort	1,,	The patient's heartbeat has interfered with his respiration. Check the
	Resp Artifact	Н	patient's condition and the Resp connections.
			The SpO2 value has fallen below the desaturation alarm limit. Check
	SpO₂ Desat	Н	the patient's condition and check if the alarm limit settings are
SpO ₂			correct.
3pO ₂			The pulse signal was so weak that the monitor cannot perform pulse
	No Pulse	Н	analysis. Check the patient's condition, SpO₂ sensor and
			measurement site.
CO ₂	CO ₂ Apnea	Н	The patient stops breathing, or the respiration signal was so weak
	AG Apnea	Н	that the monitor cannot perform respiration analysis. Check the
AG	ла прпеа	''	patient's condition and the RM connections.
Au	FiO₂ Too Low	Н	Check the patient's condition, the ventilated O₂ content and the AG
	1102 100 LOW		connections.

D.2 Technical Alarm Messages

Measurement	Alarm message	L	I	Cause and solution
	XX SelfTest Err		С	
	XX Init Err		Α	An error occurred to the XX module, or there is a problem
	XX Init Err N	Н	Α	with the communications between the module and the
	N is within 1 to 8			monitor. Re-plug the module and restart the monitor, or
XX	XX Comm Err		Α	plug the module into another monitor.
**	XX Comm Stop	Н	С	
	XX Limit Err		С	XX parameter limit is accidentally changed. Contact your
	AA LIIIII LII	L		service personnel.
	XX Overrange	L	С	The measured XX value is not within the specified range
	AX Overlange	_	C	for XX measurement. Contact your service personnel.
ECG	ECG Lead Off	L*	В	The electrode has become detached from the patient or
	ECG YY Lead Off	L*	В	the lead wire has become disconnected from the adapter

Measurement	Alarm message	L	I	Cause and solution
	Note: YY represents the leadwires,	, V (V1,	V2,	cable. Check the connections of the electrodes and
	V3, V4, V5, V6,), LL, LA, RA, as per A	·ΗΑ		leadwires.
	standard, or C (C1, C2, C3, C4, C5,	C6), F, L	-	
	and R as per IEC standard.			
				The ECG signal is noisy. Check for any possible sources of
	ECG Noisy	L	Α	signal noise around the cable and electrode, and check the
				patient for great motion.
				Artifacts are detected on the ECG analysis lead and as a
				result heart rate cannot be calculated and Asystole, Vfib
				and Vtac cannot be analyzed. Check the connections of
	ECG Artifact	L	Α	the electrodes and leadwires and check for any possible
				source of interference around the cable and electrode.
				Check the patient's condition and check the patient for
				great motion.
				High frequency signals are detected on the ECG analysis
	ECG High Freq. Noise	L	Α	lead. Check for any possible source of interference around
				the cable and electrode.
				Low frequency signals are detected on the ECG analysis
	ECG Low Freq. Noise	L	Α	lead. Check for any possible source of interference around
				the cable and electrode.
				The ECG amplitude didn't reach the detected threshold.
	ECG Amplitude Too Small	L	С	Check for any possible source of interference around the
				cable and electrode.
				ECG configuration is wrongly downloaded. Check the
	ECG Config. Err	L	C	downloaded configuration and re-download the correct
				configuration.
Resp	Resp Disturbed	L	Α	The respiration circuit is disturbed. Restart the monitor.
	Temp Cal. Err	Н	С	A calibration failed. Restart the monitor.
Temp	T1 Sensor Off	L	Α	The Temp sensor has become detached from the patient or
	T2 Sensor Off	L	Α	the module. Check the sensor connections.
SpO ₂	SpO ₂ Sensor Off	L*	В	The SpO ₂ sensor has become detached from the patient or
	SpO₂ Sensor Fault	L	С	the module, or there is a fault with the SpO ₂ sensor, or an
	SpO₂ No Sensor	L	В	unspecified SpO ₂ sensor has been used. Check the sensor
	SpO ₂ Unknown Sensor	L	C	application site and the sensor type, and make sure if the
	SpO₂ Sensor Incompatible	L	С	sensor is damaged. Reconnect the sensor or use a new
	традости по разоти			sensor.
				There is too much light on the SpO₂ sensor. Move the
	SpO₂ Too Much Light	L	С	sensor to a place with lower level of ambient light or cover
				the sensor to minimize the ambient light.
	SpO₂ Low Signal	L	С	The SpO₂ signal is too low or too weak. Check the patient's
	SpO₂ Weak Signal	L	С	condition and change the sensor application site. If the
	SpO₂ Weak Pulse	L	С	error persists, replace the sensor.
				error persists, replace the sensor.
		<u> </u>	l	

Measurement	Alarm message	L	1	Cause and solution
				The SpO2 signal has been interfered. Check for any
	SpO₂ Interference	L	C	possible sources of signal noise around the sensor and
				check the patient for great motion.
	C.O. Daniel Family		_	There is a problem with the SpO ₂ measurement board. Do
	SpO₂ Board Fault	L	С	not use the module and contact your service personnel.
	NIBP Loose Cuff	L	Α	The NIBP cuff is not properly connected, or there is a leak
	NIBP Air Leak	L	Α	in the airway.
	NIBP Pneumatic Leak	L	Α	Check the NIBP cuff and pump for leakages.
	AUDD Coff Torre Washing			The cuff type applied mismatches the patient category.
	NIBP Cuff Type Wrong	L	Α	Verify the patient category and replace the cuff.
				An error occurred to the air pressure. Verify that the
	NIDD Air Droccure Err	١,	_	monitor application site meets the environmental
	NIBP Air Pressure Err	L	Α	requirements and check if there is any source that affects
				the air pressure.
				The patient's pulse is weak or the cuff is loose. Check the
	NIBP Weak Signal	L	Α	patient's condition and change the cuff application site. If
				the error persists, replace the cuff.
NIBP	NIBP Signal Saturated		Α	The NIBP signal is saturated due to excess motion or other
	Tribi Signal Saturated	L	^	sources.
	NIBP Overrange	L	Α	The measured NIBP value is not within the specified range.
	NIBP Excessive Motion	L	Α	Check the patient's condition and reduce the patient
	THE EXCESSIVE MOTION	_		motion.
	NIBP Cuff Overpress.	L	Α	The NIBP airway may be occluded. Check the airway and
			'`	measure again.
	NIBP Equip Err	Н	Α	An error occurred during NIBP measurement and therefore
	NIBP Timeout	L	Α	the monitor cannot perform analysis correctly. Check the
	NIBP Measure Failed	L	Α	patient's condition and NIBP connections, or replace the
				cuff.
	NIBP Illegally Reset	L	Α	An illegal reset occurred during NIBP measurement. Check
				if the airway is occluded.
	YY Sensor Off			Check the sensor connection and reconnect the sensor.
IBP	YY Non-Pulsatile	L	Α	The catheter may be occluded. Please flush the catheter.
	YY represents an IBP label.	1	1	·
C.O.	TB Sensor Off	L	Α	Check the sensor connection and reconnect the sensor.
CO ₂	CO₂ Sensor High Temp	L	С	Check, stop using or replace the sensor.
	CO₂ Sensor Low Temp	L	С	Check, stop using or replace the sensor.
				The operating temperature of the CO2 module goes
	CO₂ Temp Overrange	L	С	beyond the specified range. After it restores within the
		<u> </u>		specified range, the module will restart automatically.
	CO ₂ Airway High Press.	L	С	An error occurred in the airway pressure. Check the patient
	CO₂ Airway Low Press.	L	С	connection and patient circuit, and then restart the
	•	ļ	_	monitor.
	CO ₂ High Barometric Press.	L	С	Check the CO2 connections, make sure that the monitor

Measurement	Alarm message	L	1	Cause and solution		
				application site meets the requirements, and check for		
	CO ₂ Low Barometric Press.	L	С	special sources that affect the ambient pressure. Restart		
				the monitor.		
	CO ₂ FilterLine Occluded	L	С	The airway or watertrap was occluded. Check the airway		
				and remove the occlusion.		
	CO ₂ No Watertrap	L	В	Check the watertrap connections.		
	CO ₂ Check Adapter	L	Α	There is a problem with the airway adapter. Check, clean or		
		L		replace the adapter.		
	CO ₂ FilterLine Err	L	С	Check if there is a leak in the CO ₂ sample line or the CO ₂		
		_		sample line has been occluded.		
	CO₂ Zero Failed	L	Α	Check the CO₂ connections. After the sensor's temperature		
				becomes stabilized, perform a zero calibration again.		
	CO ₂ System Err	L	Α	Re-plug the module or restart the monitor.		
	CO₂ Check Cal.	L	С	Perform a calibration.		
	CO ₂ Check Airway	L	С	An error occurred to the airway.		
	CO ₂ No Filterline	L	Α	Make sure that the filterline is connected.		
	CO ₂ No Sensor	L	Α	Make sure that the sensor is connected.		
	CO ₂ Main Board Err	Н	С			
	CO₂ Checking Sensor	L	С	There is a problem with the CO ₂ module. Re-plug the		
	CO ₂ Replace Scrubber&Pump	L	С	module or restart the monitor.		
	CO ₂ 15V Overrange	Н	С			
	CO₂ Hardware Err	Н	С			
	AG No Watertrap	L	В	Check the connections of the watertrap and re-connect it.		
	AG Change Watertrap	L	Α	Wait until the change is completed.		
	AG Watertrap Type Wrong	L	Α	Make sure that a correct watertrap has been used.		
	O ₂ Accuracy Unspecified	L	Α			
	N₂O Accuracy Unspecified	L	Α	The measured value has exceeded the specified accuracy range.		
	CO ₂ Accuracy Unspecified	L	Α			
	Enf Accuracy Unspecified	L	Α			
	Iso Accuracy Unspecified	L	Α			
AG	Sev Accuracy Unspecified	L	Α			
	Hal Accuracy Unspecified	L	Α			
	Des Accuracy Unspecified	L	Α			
	awRR Accuracy Unspecified	L	Α			
	AG Hardware Err	Н	Α	Remove the AG module. Stop using the module and		
	AC Aimmen C. I. I. I.	,	_	contact your service personnel.		
	AG Airway Occluded AG Zero Failed	L	Α	Check the airway and remove the occlusion.		
		L	Α	Re-plug the module or restart the monitor, and then		
Dawe	13//	,.	-	perform a zero calibration again.		
Power	12V Too High	Н	С	There is a problem with the system power supply. Restart		
	12V Too Low	Н	С	the monitor.		
	5V Too High	Н	С	-		
	5V Too Low	Н	С			

Measurement	Alarm message	L	I	Cause and solution		
	3.3V Too High	Н	С			
	3.3V Too Low	Н	С			
	Battery Too Low	Н	С	Connect the monitor to an AC power source and allow the batteries to charge.		
	Different Battery Voltages	М	С	The two batteries have different charge capacity, or the batteries unspecified have been used, or there is a problem with the batteries. Make sure that correct batteries are used and the batteries are not damaged, or replace the batteries.		
	Battery Power Overload	Н	С	The power consumption of the equipment is too high. Power the monitor using an AC power source.		
	RT Clock Not Exist	Н	С	Contact your service personnel.		
	Recorder Init Err N	L	Α	Restart the monitor.		
	N is within 1 to 8.			Restart the monitor.		
	Recorder SelfTest Err	L	Α			
	Recorder Comm Err	L	Α	Stop the recording and restart the monitor.		
	Recorder S. Comm Err	L	Α			
Recorder	Recorder Unavailable	L	Α			
necoraci	Recorder VIt High	L	С	An error occurred to the system power supply. Restart the		
	Recorder VIt Low	L	С	monitor.		
	Recorder Head Hot	L	С	The recorder has been working for too long time. Stop the recording and resume the recording till the recorder's		
				printhead cools down.		
	Rec Paper Wrong Pos.	L	Α	Re-load the recorder paper.		
System	System Watchdog Err	Н	С	An error occurred to the system. Restart the monitor.		
	System Software Err	Н	С			
	System CMOS Full	Н	С			
	System CMOS Err	Н	С			
	System FPGA Err	Н	С			
	System Err N	Н С				
	N is within 2 to 12.					

E Symbols and Abbreviations

E.1 Symbols

μΑ microampere μ۷ microvolt Microsecond μs ampere Α Ah ampere hour bpm beat per minute bps bit per second ٥C centigrade cubic centimeter cc centimeter cm dB decibel DS dyne second ٥F fahrenheit gram g GHz gigahertz GTT gutta hour h hertz Hz inch in kg kilogram kPa kilopascal L litre

mAh milliampere hour
Mb mega byte
mcg microgram
mEq milli-equivalents
mg milligram

pound

meter

min minute
ml milliliter
mm millimeter

lb

m

mmHg millimeters of mercury cmH2O centimeters of water

 $\begin{array}{ll} ms & millisecond \\ mV & millivolt \\ mW & milliwatt \\ M\Omega & megaohm \end{array}$

nm nanometer

rpm breath per minute

s second V volt

VA volt ampere

 $\begin{array}{ccc} \Omega & & \text{ohm} \\ W & & \text{watt} \end{array}$

- minus, negative

% percent

/ per; divide; or

+ plus
= equal to
< less than
> greater than

≤ less than or equal to≥ greater than or equal to

 $\begin{array}{ll} \pm & \text{plus or minus} \\ \times & \text{multiply} \end{array}$

E.2 Abbreviations

AaDO₂ alveolar-arterial oxygen gradient

AAMI Association for Advancement of Medical Instrumentation

AC alternating current ACI acceleration index

Adu adult

AG anaesthesia gas

AHA American Heart Association

ANSI American National Standard Institute

Ao aortic pressure

Art arterial

aVF left foot augmented lead
aVL left arm augmented lead
aVR right arm augmented lead
awRR airway respiratory rate
BAP brachial arterial pressure

BIS bispectral index BP blood pressure

BPSK binary phase shift keying

BSA body surface area
BT blood temperature

BTPS body temperature and pressure, saturated

C.I. cardiac index

CCI Continuous Cardiac Index

C.O. cardiac output

CCO Continuous Cardiac Output
CaO2 arterial oxygen content
CCU cardiac (coronary) care unit
CE Conformité Européenne
CFI cardiac function index
CIS Clinical Information System

CISPR International Special Committee on Radio Interference

CMOS complementary metal oxide semiconductor

CMS central monitoring system

CO2 carbon dioxide

COHb carboxyhemoglobin

CP cardiopulmonary

CPI cardiac power index

CPO Cardiac Power Output

CVP central venous pressure

DC direct current

Des desflurane

Dia diastolic

DPI dot per inch

dPmx left ventricular contractility
DVI digital video interface

DO₂ oxygen delivery

DO₂l oxygen delivery index ECG electrocardiograph EDV end-diastolic volume

EEC European Economic Community

EEG electroencephalogram

EMC electromagnetic compatibility

EMG electromyography

EMI electromagnetic interference

Enf enflurane

ESU electrosurgical unit

Et end-tidal

EtCO₂ end-tidal carbon dioxide EtN₂O end-tidal nitrous oxide

 $\begin{array}{ll} \text{EtO} & \text{ethylene oxide} \\ \text{EtO}_2 & \text{end-tidal oxygen} \end{array}$

EVLW extravascular lung water
ELWI extravascular lung water index

FAP femoral arterial pressure

FCC Federal Communication Commission

FDA Food and Drug Administration

FEV1.0% first second forced expiratory volume ratio

Fi fraction of inspired

FiCO₂ fraction of inspired carbon dioxide

FiN₂O fraction of inspired nitrous oxide

FiO₂ fraction of inspired oxygen FPGA field programmable gate array

FV flow-volume

GEDV global end diastolic volume

GEDI global end diastolic volume index

GEF global ejection fraction

Hal halothane
Hct haematocrit
Hb hemoglobin

Hb-CO carbon mono-oxide hemoglobin

 \mbox{HbO}_2 oxyhemoglobin \mbox{HR} heart rate

I:E inspiratory-expiratory ratio
IBP invasive brood pressure
ICG impedance cardiography
ICP intracranial pressure

ICT/B intracranial catheter tip pressure transducer

ICU intensive care unit

ID identification

IEC International Electrotechnical Commission
IEEE Institute of Electrical and Electronic Engineers

Ins inspired minimum
IP internet protocol

Iso isoflurane

IT injectate temperature

ITBI Intrathoracic Blood Volume Index

ITBV Intrathoracic Blood Volume

LA left arm

LAP left atrial pressure

Lat lateral

LCD liquid crystal display LCW left cardiac work

LCWI left cardiac work index LED light emitting diode

LL left leg

LVD low voltage directive

LVDS low voltage differential signal
LVET left ventricular ejection time
LVSW left ventricular stroke work
LVSWI left ventricular stroke work index
MAC minimum alveolar concentration

Art mean mean arterial pressure
MDD Medical Device Directive

MetHb methemoglobin

MRI magnetic resonance imaging

MVe expiratory minute volume

MVi inspiratory minute volume

 $\begin{array}{ll} N/A & \text{not applied} \\ N_2 & \text{nitrogen} \\ N_2O & \text{nitrous oxide} \\ Neo & \text{neonate} \end{array}$

NIBP noninvasive blood pressure

O₂ oxygen

 O_2CI oxygen consumption index O_2R oxygen extraction ratio

OR operating room

oxyCRG oxygen cardio-respirogram

PA pulmonary artery

pArt-D diastolic artery pressure
pArt-M mean artery pressure
pArt-S systolic artery pressure

Paw airway pressure

PAWP pulmonary artery wedge pressure

PD photodetector Ped pediatric

PEEP positive end expiratory pressure

PEF peak expiratory flow
PEP pre-ejection period
PIF peak inspiratory flow

PIP peak inspiratory pressure

Pleth plethysmogram
Pmean mean pressure
Pplat plateau pressure

PPV Pulse Pressure Variation

PR pulse rate

PVC premature ventricular contraction
PVR pulmonary vascular resistance

PVRI pulmonary vascular resistance index PVPI pulmonary vascular permeability index

PPV pulse pressure variation

pArt artery pressure

pCVP central venous pressure

R right RA right arm

RAM random access memory
RAP right atrial pressure
Rec record, recording

Resp respiration

RHb reduced hemoglobin

RL right leg

RM respiratory mechanics

RR respiration rate

RSBI rapid shallow breathing index SaO_2 arterial oxygen saturation SEF spectral edge frequency

Sev sevoflurane

SFM self-maintenance SI stroke index

SMR satellite module rack

SpO₂ arterial oxygen saturation from pulse oximetry

SQI signal quality index
SR suppression ratio
STR systolic time ratio
SV stroke volume

SVI Stroke Volume Index

SVR systemic vascular resistance

SVRI systemic vascular resistance index

SVV stroke volume variation

 SvO_2 mixed venous oxygen saturation $ScvO_2$ central venous oxygen saturation

Sync synchronization
Sys systolic pressure
Taxil axillary temperature
TB Blood Temperature
TD temperature difference

Temp temperature

TFC thoracic fluid content
TFI thoracic fluid index
TFT thin-film technology
Toral oral temperature

TP total power

Trect rectal temperature

TVe expiratory tidal volume

TVi inspiratory tidal volume

UAP umbilical arterial pressure

UPS uninterruptible power supply

USB universal serial bus

UVP umbilical venous pressure VAC volts alternating current

VEPT volume of electrically participating tissue

VI velocity index

VO₂ oxygen consumption

VO₂I oxygen consumption index WLAN wireless local area network